

101.418-12/04 – including *Taq* polymerase
 101.418-12u/04u – without *Taq* polymerase

Lot No.: **60E**

Olerup SSP[®] HLA-A*68

Product number: 101.418-12/04 – including *Taq* pol.
 101.418-12u/04u – without *Taq* pol.

Lot number: 60E

Expiry date: 2010-June-01

Number of tests: 12 tests – Product No. 101.418-12
 4 tests – Product No. 101.418-04

Number of wells per test: 32

CHANGES COMPARED TO THE PREVIOUS HLA-A*68 LOT (X65):

Well	5'-primer	3'-primer	rationale
13	-	Modified	Increased specificity of specific primer pair.
14	Added	Added	Primer pair added for the A*6841 allele.
16	-	Added	Primer added for the A*6840 allele.
29	Added	Added	Primer pair added for the A*6842 allele.
32	Exchanged	-	Primer exchange to decrease primer dimer formation tendency.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*6801 to A*6842, i.e. all the currently recognized HLA-A*68 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-A*68 HOMO- AND HETEROZYGOTES:

The A*6801,6801 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-A*68 SUBTYPINGS BY NON-HLA-A*68 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.