

101.435-06 – including *Taq* polymerase101.435-06u – without *Taq* polymeraseLot No.: **3E9****Olerup SSP® HLA-A\*69**

**Product number:** 101.435-06 – including *Taq* polymerase  
101.435-06u – without *Taq* polymerase

**Lot number:** 3E9

**Expiry date:** 2019-04-01

**Number of tests:** 6

**Number of wells per test:** 5+1

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
6	-	-	Updated Negative Control

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

All the HLA-A\*69 alleles, i.e. **A\*69:01 to A\*69:03**, recognized by the HLA Nomenclature Committee in July 2016<sup>1</sup> will be amplified by the primers in the HLA-A\*69 SSP kit.

The HLA-A\*69 kit enables separation of the confirmed HLA-A\*69 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A\*69 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HLA-A\*69 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*69 SUBTYPINGS BY NON-HLA-A\*69 ALLELES:**

None of importance