

## **Olerup SSP® HLA-A\*74**

<b>Product number:</b>	<b>101.433-06 – including <i>Taq</i> polymerase</b> <b>101.433-06u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>9D3</b>
<b>Expiry date:</b>	<b>2018-10-01</b>
<b>Number of tests:</b>	<b>6</b>
<b>Number of wells per test:</b>	<b>15+1</b>

### **CHANGES COMPARED TO THE PREVIOUS HLA-A\*74 LOT (13Y):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
16	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged

#### **ALLELE COVERAGE:**

All the HLA-A\*74 alleles, i.e. **A\*74:01 to A\*74:27 alleles**, recognized by the HLA Nomenclature Committee in January 2016<sup>1</sup> will be amplified by the primers in the HLA-A\*74 subtyping kit.

The HLA-A\*74 kit enables separation of the confirmed HLA-A\*74 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-A\*74 alleles is listed below.

The HLA-A\*74 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The HLA-A\*74 subtyping kit cannot distinguish the silent mutations in the A\*74:01:01-74:01:03, the A\*74:02:01:01-74:02:01:02 alleles or the A\*74:16:01-74:16:02 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

#### **RESOLUTION IN HLA-A\*74 HOMO- AND HETEROZYGOTES:**

The A\*74:01,74:01 genotype gives rise to a unique amplification pattern.

#### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*74 SUBTYPINGS BY NON-HLA-A\*74 ALLELES:**

None frequently occurring.

#### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.