

101.524-12 – including *Taq* polymerase

101.524-12u – without *Taq* polymerase

Lot No.: **4D1**

## **Olerup SSP® HLA-B\*14**

**Product number:** 101.524-12 – including *Taq* polymerase  
 101.524-12u – without *Taq* polymerase

**Lot number:** 4D1

**Expiry date:** 2018-07-01

**Number of tests:** 12

**Number of wells per test:** 21+1

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*14 LOT (61X):**

| Well | 5'-primer | 3'-primer | rationale   |
|------|-----------|-----------|---|
| 2    | Added     | -         | 5'-primer added for the B*14:01:06 allele.                                    |
| 8    | Added     | -         | 5'-primer added for the B*14:01:06 allele.                                    |
| 19   | Added     | Added     | Negative control moved to well 22, primer pair added for the B*14:41N allele. |
| 20   | Added     | Added     | New primer pair added for the B*14:33 allele.                                 |
| 21   | Added     | Added     | New primer pair added for the B*14:46 allele.                                 |
| 22   | -         | -         | Negative control added from well 19.  |

**THE NUMBER OF WELLS** is increased from 19 to 22 wells.

#### **ALLELE COVERAGE:**

B\*14:01 to B\*14:50 i.e. all the currently recognized HLA-B\*14 alleles, will be amplified by the primers in the HLA-B\*14 subtyping kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-October-10, release 3.22.0.

The HLA-B\*14 kit enables separation of the confirmed HLA-B\*14 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B\*14 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The following HLA-B\*14 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

| Alleles        | Primer mix |
|----------------|------------|
| B*14:16, 14:23 | 16         |

The HLA-B\*14 primer set cannot distinguish the following silent mutations: the B\*14:01:01-14:01:06, the B\*14:02:01-14:02:13, the B\*14:06:01-14:06:02 alleles or the B\*14:08:01-14:08:02 alleles.

101.524-12 – including *Taq* polymerase

101.524-12u – without *Taq* polymerase

Lot No.: **4D1**

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page

<http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The B\*14:08:01-14:08:02 and the B\*39:43 alleles will give rise to identical amplification patterns with the HLA-B\*14 subtyping kit. These alleles can be distinguished by the HLA-B low resolution kit and/or the HLA-B\*39 subtyping kit.

**RESOLUTION IN HLA-B\*14 HOMO- AND HETEROZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*14 SUBTYPINGS BY NON-HLA-B\*14 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.