

101.523-24/03 – including *Taq* polymerase

101.523-24/03u – without *Taq* polymerase

Lot No.: **5F8**

## Olerup SSP® HLA-B\*40

**Product number:** 101.523-24/03 – including *Taq* pol.  
101.523-24/03u – without *Taq* pol.

**Lot number:** 5F8

**Expiry date:** 2021-10-01

**Number of tests:** 24 tests – Product No. 101.523-24/24u  
3 tests – Product No. 101.523-03/03u

**Number of wells per test:** 94+1

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*40 LOT (9D1):

Well	5'-primer	3'-primer	rationale
2	Added	-	5'-primer added for the B*40:06:15 allele.
24	Added	-	5'-primer added for the B*40:334 allele.
31	Modified	-	5'-primer modified for improved HLA-specific amplification.
32	Added	-	5'-primer added for the B*40:06:15 allele.
36	Added	-	5'-primer added for the B*40:273 allele.
42	Added	-	5'-primer added for the B*40:273 allele.
47	Added	Modified	5'-primer added for the B*40:318 allele, 3'-primer modified for decreased tendency of primer oligomer formation.
51	Added	-	5'-primer added for the B*40:345N allele.
57	Exchanged	Exchanged	Primer pair exchanged for improved HLA-specific amplification.
59	Added	-	5'-primers added for the B*40:331 and B*40:346 alleles.
67	-	Added	3'-primer added for the B*40:320 allele.
68	-	Added	3'-primer added for the B*40:320 allele.
73	Added	-	5'-primer added for the B*40:329 allele.
80	Added	-	5'-primer added for the B*40:334 allele.
82	Added	-	5'-primers added for the B*40:311, B*40:331 and B*40:346 alleles.
83	-	Added	3'-primer added for the B*40:348 allele.
85	Added	Added	Primer pair added for the B*40:337N allele.
90	Added	-	5'-primer added for the B*40:329 allele.
91	-	Added	3'-primer added for the B*40:338N allele.
93	-	Added	3'-primer added for the B*40:335 allele.
94	Added	Added	Negative control moved to well 95, primer pair added for the B*40:353 allele.
95	-	-	Negative control added from well 94.

**THE NUMBER OF WELLS** is increased from 94 to 95 wells.

### ALLELE COVERAGE:

B\*40:01 to B\*40:355, i.e. all the currently recognized HLA-B\*40 alleles, will be amplified by the primers in the HLA-B\*40 subtyping kit<sup>1</sup>, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2017-April-13, release 3.28.0.

101.523-24/03 – including *Taq* polymerase

101.523-24/03u – without *Taq* polymerase

Lot No.: **5F8**

The HLA-B\*40 kit enables separation of the confirmed HLA-B\*40 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B\*40 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The following HLA-B\*40 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
B*40:69, 40:221	56
B*49:179, 40:253	83

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

#### RESOLUTION IN HLA-B\*40 HOMO- AND HETEROZYGOTES:

Good.

#### INFLUENCE ON THE INTERPRETATION OF HLA-B\*40 SUBTYPINGS BY NON-HLA-B\*40 ALLELES:

None frequently occurring.

#### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

In primer mix 57, the primer pair has been exchanged for improved HLA-specific amplification.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.