

## Olerup SSP<sup>®</sup> HLA-B\*41

**Product number:** 101.542-06 – including *Taq* polymerase  
101.542-06u – without *Taq* polymerase  
**Lot number:** 0R3  
**Expiry date:** 2026-03-01  
**Number of tests:** 6  
**Number of wells per test:** 23+1

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*41 LOT (9K0).

Well	5'-primer	3'-primer	rationale
2	-	Added	3'-primer added for the B*41:71 allele.

THE NUMBER OF WELLS is unchanged.

#### ALLELE COVERAGE:

**B\*41:01 to B\*41:75** i.e. all the currently recognized HLA-B\*41 alleles, will be amplified by the primers in the HLA-B\*41 SSP kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2022-January-13, release 3.47.0.

The HLA-B\*41 kit enables separation of the confirmed HLA-B\*41 alleles as listed in the IMGT/HLA database 3.27.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B\*41 kit also enables identification of many null and alternatively expressed alleles.

The following HLA-B\*41 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
B*41:04, 41:36	4
B*41:10, 41:45N	3

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page

<http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The following alleles give rise to identical amplification patterns with the HLA-B\*41 high resolution kit. These alleles can be distinguished by the HLA-B low resolution and/or the respective kits.

#### Alleles

B\*41:08, B\*40:136, 40:231, 40:279

B\*41:23, B\*42:27

B\*41:46-41:47, B\*44:166

B\*41:62, B\*40:02:20



0197

HLA-B\*41

101.542-06 – including *Taq* polymerase

101.542-06u – without *Taq* polymerase

Lot No.: **OR3**

**RESOLUTION IN HLA-B\*41 HOMO- AND HETEROZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*41 SUBTYPINGS BY NON-HLA-B\*41 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.



0197

For *In Vitro* Diagnostic Use

MA125 v01 SSP Release Note Template

Date: March 2022, Rev. No: 00