

101.563-24/03 – including *Taq* polymerase

101.563-24u/03u – without *Taq* polymerase

Lot No.: 11Y

Olerup SSP® HLA-B*44

Product number: 101.563-24/03 – including *Taq* pol.
101.563-24u/03u - without *Taq* pol.

Lot number: 11Y

Expiry date: 2017-September-01

Number of tests: 24 tests – Product No. 101.563-24/24u
3 tests – Product No. 101.563-03/03u

Number of wells per test: 63+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*44 LOT (14V):

| Well | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|---------------------------------------------------------------------------------------|
| 1 | Added | - | 5'-primer added for the B*44:02:34 allele. |
| 2 | - | Moved | 3'-primer modified for improved HLA-specific amplification. |
| 3 | Added | - | 5'-primer added for the B*44:02:34 allele. |
| 8 | Moved | Moved | Primer pair moved to well 62. |
| 16 | - | Added | 3'-primer added for the B*44:50:02 allele. |
| 22 | - | Added | 3'-primer added for the B*44:177 allele. |
| 35 | - | Added | 3'-primer added for the B*44:137 allele. |
| 47 | - | Added | 3'-primer added for the B*44:191 allele. |
| 48 | - | Added | 3'-primers added for the B*44:191 and B*44:195N. |
| 50 | - | Added | 3'-primer added for the B*44:137 allele. |
| 52 | Added | - | 5'-primer added for the B*44:173 allele. |
| 54 | - | Added | 3'-primer added for the B*44:195N allele. |
| 56 | - | Added | 3'-primer added for the B*44:198N allele. |
| 60 | Added | Added | Updated negative control moved to well 64, primer pair added for the B*44:187 allele. |
| 61 | New | New | New primer pair for the B*44:188 allele. |
| 62 | Added | Added | Primer pair added from well 8. |
| 63 | New | New | New primer pair for the B*44:173 allele. |
| 64 | - | - | Updated negative control added from well 60. |

THE NUMBER OF WELLS is increased from 60 to 64 wells.

ALLELE COVERAGE:

B*44:01 to B*44:214, i.e. all the currently recognized HLA-B*44 alleles, will be amplified by the primers in the HLA-B*44 subtyping kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2015-January-19, release 3.19.0.

The HLA-B*44 kit enables separation of the confirmed HLA-B*44 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

101.563-24/03 – including *Taq* polymerase
101.563-24u/03u – without *Taq* polymerase

Lot No.: **11Y**

The HLA-B*44 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The following HLA-B*44 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

| Alleles | Primer mix | Alleles | Primer mix |
|------------------|------------|--------------------|------------|
| B*44:32, 44:124 | 21 | B*44:78, 44:136 | 48 |
| B*44:58N, 44:139 | 47 | B*44:94, 44:98 | 31 |
| B*44:69, 44:85 | 26 | B*44:160Q, 44:198N | 56 |
| B*44:72, 44:93 | 47 | | |

The HLA-B*44 subtyping cannot distinguish the following silent mutations: the B*44:02:01:01, 44:02:01:03-44:02:09 and 44:02:11-44:02:36 alleles, the B*44:03:01-44:03:24 alleles, the B*44:05:01-44:05:04 alleles, the B*44:27:01-44:27:02 allele, the B*44:28:01-B*44:28:02 alleles, the B*44:37:01-44:37:02 allele, the B*44:41:01-44:41:02 allele, the B*44:43:01-44:43:02 alleles or the B*44:79:01-44:79:02 alleles.

¹ Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-B*44 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*44 SUBTYPINGS BY NON-HLA-B*44 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

In primer mix 2, a 3'-primer was modified for improved HLA-specific amplification.