

**101.852-12 – including *Taq* polymerase**

**101.852-12u – without *Taq* polymerase**

**Lot No.: 4E4**

## **Olerup SSP® HLA-B\*07 Add-on**

<b>Product number:</b>	<b>101.852-12 – including <i>Taq</i> polymerase</b> <b>101.852-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>4E4</b>
<b>Expiry date:</b>	<b>2019-04-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>6+1</b>

The HLA-B\*07 Add-on primer set is unchanged compared to the previous *Olerup SSP®* HLA-B\*07 Add-on (**Lot No. 4D6**).

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The HLA-B\*07:02, 07:05, 07:06, 07:61 and 07:161N alleles give different patterns in the B\*07 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-July-14, release 3.25.0.

The HLA-B\*07 Add-on kit cannot distinguish these silent mutations: the B\*07:02:01-07:02:05, B\*07:02:07-07:02:25, B\*07:02:27-07:02:49, the B\*07:05:01:01-07:05:08 or the B\*07:06:01-07:06:02 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF B\*07 ADD-ON SUBTYPINGS BY NON- HLA-B\*07 ALLELES:**

None of importance.