

HLA-C High Resolution for frequent alleles

101.707-12 – including *Taq* polymerase

101.707-12 – without *Taq* polymerase

Lot No.: 6H8

Olerup SSP[®] HLA-C high resolution for frequent alleles

Product number:	101.707-12 – including <i>Taq</i> polymerase 101.707-12u – without <i>Taq</i> polymerase
Lot number:	6H8
Expiry date:	2023-10-01
Number of tests:	12
Number of wells per test:	95 + 1

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] HLA-C HIGH RESOLUTION LOT (1F9):

The HLA-C high resolution primer set is unchanged compared to the previous Olerup SSP[®] HLA-C high resolution (Lot No. 1F9).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

C*01:02 to C*18:12, i.e. all the currently recognized HLA-C alleles, will be amplified by the primers in the HLA-C high resolution for frequent alleles subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2019-January-23, release 3.35.0.

The following HLA-C alleles can be distinguished by the different sizes of the specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
C*01:06, 01:11	3	C*06:15, 06:16N	78, 83
C*01:15:01, 01:20, 01:60	6	C*06:89, 06:161	73
C*01:113, 01:140	73	C*07:02:75, 07:679	8
C*02:51, 02:57	31	C*07:22, 07:35, 07:404	56
C*03:09, 03:23, 03:80:02, 03:287:01-03:287:02	18	C*07:25, 07:137:01, 07:138, 07:352	56
C*03:21-03:22Q, 03:116:01, 03:413	18	C*07:33N, 07:51	55
C*03:34, 03:44	80	C*07:37, 07:54	2
C*03:37:01-03:37:02, 03:46	11	C*07:47, 07:48	91
C*04:08, 04:09N	27	C*07:60, 07:109:01-07:109:02	17
C*04:16, 04:256	11	C*07:66, 07:665	90
C*04:23, 04:31, 04:277	88	C*07:130, 07:611	31
C*04:40, 04:223:02	31	C*07:140-07:141:02, 07:304, 07:512, 07:514, 07:671	5
C*05:06, 05:15	36	C*12:06,12:81	65
C*05:10, 05:55	21	C*12:07 ,12:46N	16
C*05:14, 05:16, 05:19	37	C*16:33, 16:109	73
C*06:08, 06:13	43		

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

HLA-C High Resolution for frequent alleles

101.707-12 – including *Taq* polymerase

101.707-12 – without *Taq* polymerase

Lot No.: 6H8

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.