

101.615-24/04 – including *Taq* polymerase  
 101.615-24u/04u – without *Taq* polymerase

Lot No.: **03V**

## **Olerup SSP® HLA-C\*07**

**Product number:** 101.615-24/04 – including *Taq* pol.  
 101.615-24u/04u - without *Taq* pol.  
**Lot number:** 03V  
**Expiry date:** 2016-May-01  
**Number of tests:** 24 tests – Product No. 101.615-24/24u  
 4 tests – Product No. 101.615-04/04u  
**Number of wells per test:** 76+1

**CHANGES COMPARED TO THE PREVIOUS HLA-C\*07 LOT (42S):**

Well	5'-primer	3'-primer	rationale
20	Added	-	5'-primer added for the C*07:48 allele, for improved HLA-specific amplification.
39	Added	-	5'-primer added for the C*07:308 allele.
74	-	Modified	3'-primer modified for decreased tendency of primer oligomer formation.
75	New	New	New primer pair for the C*07:116 allele
76	New	New	New primer pair for the C*07:102 allele.
77	New	New	Negative Control.

**THE NUMBER OF WELLS** is increased from 75 to 77.

**ALLELE COVERAGE:**

C\*07:01 to C\*07:315, i.e. all the currently recognized HLA-C\*07 alleles, will be amplified by the primers in the HLA-C\*07 SSP kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-July-25, release 3.13.1.

The HLA-C\*07 kit enables separation of the confirmed HLA-C\*07 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C\*07 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

Alleles	Primer mix
C*07:74, 07:308	39
C*07:51, 07:314	19
C*07:78, 02:235Q	40
C*07:93, 07:207	42
C*07:100, 07:161	49
C*07:104N, 07:119, 07:124	50
C*07:117, 07:126	51
C*07:118, 07:203	52
C*07:129, 07:153	54

101.615-24/04 – including *Taq* polymerase  
101.615-24u/04u – without *Taq* polymerase

Lot No.: **03V**

C*07:132, 07:179	55
C*07:147, 07:241	60
C*07:164N, 07:212	58
C*07:165, 07:180	60
C*07:227N, 07:254	30

The HLA-C\*07 subtyping kit cannot distinguish the following silent mutations: the C\*07:01:01:01-07:01:04, C\*07:01:06-07:01:24 and 07:01:26-07:01:34 alleles, the C\*07:02:01:01-07:02:43 alleles, the C\*07:04:01-07:04:08, the C\*07:17:01-07:17:02, the C\*07:27:01-07:27:02, the C\*07:31:01-07:31:02, the C\*07:38:01-07:38:02, the C\*07:56:01-07:56:02, the C\*07:96:01-07:96:02, the C\*07:137:01-07:137:02 or the C\*07:141:01-07:141:02 alleles.

<sup>1</sup>HLA-C alleles listed on the IMGT/HLA web page 2013-July-25, release 3.13.1, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

<sup>2</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

#### **INFLUENCE ON THE INTERPRETATION OF HLA-C\*07 SUBTYPINGS BY NON-HLA-C\*07 ALLELES:**

None frequently occurring.

#### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.