

Olerup SSP® HLA-C*12

Product number: 101.624-12 – including *Taq* polymerase
101.624-12u – without *Taq* polymerase
Lot number: 7K1
Expiry date: 2024-02-01
Number of tests: 12
Number of wells per test: 47+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*12 LOT (1H9):

Well	5'-primer	3'-primer	Rationale
10	Removed	-	5'-primer removed for improved HLA-specific amplification.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

C*12:02 to C*12:304, i.e. all the currently recognized HLA-C*12 alleles, will be amplified by the primers in the HLA-C*12 SSP kit²; www.ebi.ac.uk/imgt/hla, 2019-October-17¹, release 3.38.0.

The HLA-C*12 kit enables separation of the confirmed HLA-C*12 alleles as listed in the IMGT/HLA database 3.27.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C*12 kit also enables identification of null and alternatively expressed alleles.

The following HLA-C*12 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
C*12:03:09, 12:159	17	C*12:32, 12:34	30
C*12:06, 12:48, 12:81	6	C*12:39N, 12:167	25
C*12:11, 12:227	14	C*12:45, 12:50, 12:232N	28
C*12:23, 12:203	22	C*12:46N, 12:139	20
C*12:29, 12:38	29	C*12:86, 12:104N	29
C*12:30, 12:36	27	C*12:109, 12:171	34

101.624-12 – including *Taq* polymerase

101.624-12u – without *Taq* polymerase

Lot No.: 7K1

The following alleles give rise to identical amplification patterns with the HLA-C*12 primer set. These alleles can be distinguished by the HLA-C low resolution kit and the HLA-C*05, HLA-C*07, HLA-C*08 or HLA-C*16 kit, respectively.

Alleles

C*12:09, C*05:16, C*05:85, C*05:107, C*16:88

C*12:33, C*16:91

C*12:181, C*07:723

C*12:208, C*12:222, C*12:233, C*12:251, C*12:258, C*08:05, C*08:25, C*16:15:01

RESOLUTION IN HLA-C*12 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*12 SUBTYPINGS BY NON-HLA-C*12 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.