

## Olerup SSP® HLA-C\*14

**Product number:** 101.625-06 – including *Taq* polymerase  
101.625-06u – without *Taq* polymerase  
**Lot number:** 20L  
**Expiry date:** 2013-September-01  
**Number of tests:** 6  
**Number of wells per test:** 23

### CHANGES COMPARED TO THE PREVIOUS HLA-C\*14 LOT (50G):

Well	5'-primer	3'-primer	rationale
1	-	Added	Primer added for the C*14:27 allele.
2	-	Added	Primer added for improved yield of specific primer pair.
7	Modified	Modified	Modified primers for improved yield of specific primer pair.
10	-	Added	Primer added for the C*14:21N allele.
12	Added	Added	Primer pair added for the C*14:23 allele.
14	-	Added	Primer added for the C*14:15 allele.
17	New	New	New primer pairs for the C*14:22 and C*14:27 alleles.
18	New	New	New primer pairs for the C*14:17 and C*14:26 alleles.
19	New	New	New primer pairs for the C*14:18 and C*14:29 alleles.
20	New	New	New primer pairs for the C*14:19 and C*14:28 alleles.
21	New	New	New primer pair for the C*14:20 allele.
22	New	New	New primer pair for the C*14:24 allele.
23	New	New	New primer pair for the C*14:25 allele.

**THE NUMBER OF WELLS** has been increased from 17 to 23.

#### ALLELE COVERAGE:

C\*14:02 to C\*14:29, i.e. all the currently recognized HLA-C\*14 alleles, give rise to unique amplification patterns<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-January-14, release 3.3.0.

<sup>1</sup>The C\*14:18 and 14:29 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 19.

#### RESOLUTION IN HLA-C\*14 HOMO- AND HETEROZYGOTES:

Good.

101.625-6 – including *Taq* polymerase

101.625-6u – without *Taq* polymerase

Lot No.: **20L**

**INFLUENCE ON THE INTERPRETATION OF HLA-C\*14 SUBTYPINGS BY NON-HLA-C\*14 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.