

## Olerup SSP<sup>®</sup> HLA-C\*14

**Product number:** 101.625-06 – including *Taq* polymerase  
101.625-06u – without *Taq* polymerase  
**Lot number:** 12S  
**Expiry date:** 2015-October-01  
**Number of tests:** 6  
**Number of wells per test:** 24

### CHANGES COMPARED TO THE PREVIOUS HLA-C\*14 LOT (54N):

Well	5'-primer	3'-primer	rationale
10	-	Added	3'-primer added for the C*14:47N allele.
17	-	Added	3'-primer added for the C*14:47N allele.
22	Added	-	5'-primer added for the C*14:24:02 allele.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

C\*14:02 to C\*14:49, i.e. all the currently recognized HLA-C\*14 alleles, will be amplified by the primers in the HLA-C\*14 kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-January-12, release 3.11.0.

The HLA-C\*14 kit enables separation of the confirmed HLA-C\*14 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C\*14 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles

<sup>1</sup>The C\*14:18 and 14:29 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 19.

The C\*14:24:01-14:24:02 and 14:31 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 22.

The HLA-C\*14 subtyping kit cannot distinguish the following silent mutations: the C\*14:02:01 and 14:02:04-14:02:11 alleles.

### RESOLUTION IN HLA-C\*14 HOMO- AND HETEROZYGOTES:

Good.

### INFLUENCE ON THE INTERPRETATION OF HLA-C\*14 SUBTYPINGS BY NON-HLA-C\*14 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.