HLA-B*35 Product Insert Page 57 of 60

101.522-24/03 – including *Taq* polymerase IFU-01 Rev. No. 02 can be downloaded from

General "Instructions for Use"

Lot No.: 07L **Lot-specific Information** www.olerup-ssp.com

CERTIFICATE OF ANALYSIS

Olerup SSP® HLA-B*35 SSP

Product number: 101.522-24/03 - including *Taq* pol.

Lot number: 07L

2013-September -01 Expiry date:

Number of tests: 24 tests - Product No. 101.522-24

3 tests - Product No. 101.522-03

Number of wells per test: 72

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2009-667-01	25	2011-824-25	49	2011-824-49
2	2007-354-02	26	2009-667-26	50	2011-824-50
3	2007-354-03	27	2009-667-27	51	2011-824-51
4	2007-354-04	28	2007-354-28	52	2011-824-52
5	2011-824-05	29	2007-354-29	53	2011-824-53
6	2007-354-06	30	2007-354-30	54	2011-824-54
7	2009-667-07	31	2011-824-31	55	2011-824-55
8	2007-354-08	32	2007-354-32	56	2011-824-56
9	2008-491-09	33	2009-667-33	57	2011-824-57
10	2009-667-10	34	2007-354-34	58	2011-824-58
11	2007-354-11	35	2007-354-35	59	2011-824-59
12	2009-667-12	36	2008-491-36	60	2011-824-60
13	2007-354-13	37	2007-354-37	61	2011-824-61
14	2007-354-14	38	2011-824-38	62	2011-824-62
15	2007-354-15	39	2007-354-39	63	2011-824-63
16	2007-354-16	40	2007-354-40	64	2011-824-64

17	2007-354-17	41	2009-667-41	65	2011-824-65
18	2009-667-18	42	2009-667-42	66	2011-824-66
19	2007-354-19	43	2009-667-43	67	2011-824-67
20	2007-354-20	44	2008-491-44	68	2011-824-68
21	2007-354-21	45	2007-354-45	69	2011-824-69
22	2007-354-22	46	2011-824-46	70	2011-824-70
23	2007-354-23	47	2009-667-47	71	2011-824-71
24	2009-667-24	48	2007-354-48	72	2011-824-72

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 10, 33, 35, 36, 41 to 44 and 50 to 71 were available.

The specificities of the primers in primer solutions 10, 33, 35, 36, 41, 43, 44, 50, 58, 61, 63, 64, 68 and 69 were tested by separately adding one or two additional 5'-primers, respectively one or two additional 3'-primers.

In primer solutions 42, 52 to 55, 57, 62 and 67 it was only possible to test the 5'primers, the 3'-primers were not possible to test.

In primer solutions 51, 56, 59, 60, 65, 66, 70 and 71 it was only possible to test the 3'-primers, the 5'-primers were not possible to test.

IFU-01 Rev. No. 02 can be downloaded from

Lot No.: 07L Lot-specific Information www.olerup-ssp.com

In primer solutions 10, 16, 33, 35, 36, 43, 64 and 69 one or two of the 5'-primers were not possible to test, and in primer solutions 10, 12, 35, 41, 43, 44, 50, 58, 61, 64 and 69 one or two of the 3'-primers were not possible to test.

Additional primers in primer solutions 22, 26 and 31 were tested by separately adding one 5'-primer and/or one 3'-primer.

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2011-August-18

Approved by:

Asa Olaciss.

Quality Control, Supervisor

HLA-B*35 Product Insert Page 59 of 60 101.522-24/03 – including *Taq* polymerase General "Instructions for Use"

IFU-01 Rev. No. 02 can be downloaded from

Lot No.: 07L Lot-specific Information www.olerup-ssp.com

Declaration of Conformity

Product name:

Olerup SSP® HLA-B*35

Product number:

101.522-24/03

Lot number:

07L

Intended use:

HLA-B*35 high resolution histocompatibility testing

Manufacturer:

Olerup SSP AB Franzengatan 5

SE-112 51 Stockholm, Sweden

Phone: +46-8-717 88 27 **Fax:** +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

The Authorized Representative located within the Community is: Olerup SSP AB.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Stockholm, Sweden 2011-August-18

Ann-Cathrin Jareman

a-()

Head of QA and Regulatory Affairs