HLA-B\*51 Product Insert Page 45 of 48 101,561-24u/03u – without *Tag* polymerase General "Instructions for Use"

IEU-02 Rev. No. 02 can be downloaded from

Lot No.: 48M Lot-specific Information www.olerup-ssp.com

## CERTIFICATE OF ANALYSIS

Olerup SSP® HLA-B\*51 SSP

Product number: 101.561-24u/03u – without Tag pol.

Lot number: 48N

Expiry date: 2014-February-01

Number of tests: 24 tests - Product No. 101.561-24u

3 tests -- Product No. 101.561-03u

Number of wells per test: 64

Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2010-731-01	25	2010-731-25	49	2010-731-49
2	2010-731-02	26	2010-731-26	50	2010-731-50
3	2010-731-03	27	2010-731-27	51	2010-731-51
4	2010-731-04	28	2010-731-28	52	2010-731-52
5	2010-731-05	29	2010-731-29	53	2010-731-53
6	2010-731-06	30	2011-891-30	54	2010-731-54
7	2010-731-07	31	2010-731-31	55	2010-731-55
8	2010-731-08	32	2010-731-32	56	2010-731-56
9	2010-731-09	33	2010-731-33	57	2010-731-57
10	2010-731-09	34	2010-731-34	58	2010-731-58
11	2010-731-11	35	2010-731-35	59	2010-731-59
12	2010-731-12	36	2010-731-36	60	2010-731-60
13	2010-731-13	37	2010-731-37	61	2010-731-61
14	2010-731-14	38	2010-731-38	62	2010-731-62
15	2010-731-15	39	2010-731-39	63	2010-731-63
16	2010-731-16	40	2010-731-40	64	2010-731-64
17	2010-731-17	41	2010-731-41	-	
18	2010-731-18	42	2010-731-42	1	
19	2010-731-19	43	2010-731-43	1	
20	2010-731-20	44	2010-731-44	1	
21	2010-731-21	45	2010-731-45		
22	2010-731-22	46	2010-731-46		
23	2010-731-23	47	2010-731-47	1	
24	2010-731-24	48	2010-731-48	1	

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 5, 13, 17, 20, 21, 25, 27, 32, 33, 36, 37, 39, 42 to 46, 51 to 53, 55 and 57 to 64 were available. The specificities of the primers in primer solutions 5, 17, 27, 32, 37, 39, 42, 44, 46, 51, 57 and 61 were tested by separately adding one additional 5'-primer, respectively one additional 3'-primer. In primer solutions 13, 20, 21, 25, 33, 43, 45, 53, 55, 59, 60 and 62 it was only possible to test the 5'-primer, the 3'-primer was not possible to test. In primer solutions 36, 52, 58, 63 and 64

HLA-B\*51 Product Insert 101.561-24u/03u – without *Taq* polymerase

Page 46 of 48 General "Instructions for Use"

IFU-02 Rev. No. 02 can be downloaded from

Lot No.: 48M

Lot-specific Information

www.olerup-ssp.com

it was only possible to test the 3'-primer, the 5'-primer was not possible to test. In primer solutions 9, 16, 25, 39, 41, 42 and 44 one or two of the 5'-primers were not possible to test. In primer solution 15, 18, 26 to 28, 37, 42, 51, 57 and 61 one or two of the 3'-primers were not possible to test. Additional primers in primer solutions 18, 26, 28, 31, 41, 47 and 48 were tested by separately adding one 5'-primer or one 3'-primer.

Results:

No false positive or false negative amplifications were obtained.

Date of approval: 2011-August-18

Approved by:

Asa Olarku

Quality Control, Supervisor

HLA-B\*51 Product Insert Page 47 of 48

101.561-24u/03u - without Taq polymerase

General "Instructions for Use"

Lot No.: 48M Lot-specific Information www.olerup-ssp.com

**Declaration of Conformity** 

Product name:

Olerup SSP® HLA-B\*51

Product number:

101.561-24u/03u

Lot number:

48M

Intended use:

HLA-B\*51 histocompatibility testing

Manufacturer:

Olerup SSP AB Franzengatan 5

SE-112 51 Stockholm, Sweden

**Phone:** +46-8-717 88 27 **Fax:** +46-8-717 88 18

We, Olerup SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

The Authorized Representative located within the Community is: Olerup SSP AB.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Stockholm, Sweden 2011-August-18

Ann-Cathrin Jareman

Head of QA and Regulatory Affairs

Mu Cathreis Januar