

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01-01:02:17, 01:02:19-01:02:64, 01:05-01:08, 01:11-01:13, 01:15-01:01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79:01-01:89N, 01:91-01:100, 01:103-01:112, 01:114:01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:146-01:01:155, 01:157-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187	73																			92				
*01:02:18	73											84								92				
*01:03, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:176	73																			92	93			
*01:04	73											84		86						92				
*01:09:01-01:09:02	73													86						92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118, 01:186	73											83								92				
*01:21	73											84	85											
*01:43	73	74									82									92				
*01:60	73																		91	92				
*01:67, 01:122	73			76												88				92				
*01:90, 01:101-01:102, 01:113	73							80									89			92				
*01:116	73																89			92				
*01:125, 14:58 ³	73															88				92				
*01:131	73			76																92	93			
*01:136, 01:145N	73							80												92				
*01:156, 14:110, B*07:239, B*14:03 ³																				92				
*01:160	73													86				90		92				
*01:168-01:169:01																	87			92				
*01:169:02, 03:94 ³			75														87			92		94		
*01:177																	88			92		94		
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:31-02:02:37, 02:02:39, 02:02:41-02:02:43, 02:02:45-02:02:52, 02:02:54-02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186		74	75			78					83				87									
*02:02:13, 02:02:29		74	75								83				87							94		
*02:02:26, 02:21, 02:26:02-02:26:03, 02:39, 02:45		74	75								83				87									
*02:02:30, 02:02:40		74				78					83				87									
*02:02:38, 02:02:53		74	75			78					83				87							94		
*02:02:44		74	75			78					83	84			87									
*02:04		74	w			78					83				87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87		89							
*02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01		74	w			78					83													
*02:14:02, 02:164		74				78					83													
*02:22, 02:62, 02:118		74	75			78					83				87							93		
*02:26:01, 02:134, 02:152, 02:167		74									83				87									
*02:27:01-02:27:02, 02:65, 02:131		74	75			78					83													
*02:40:01-02:40:02		74	75			78					83				87				91					
*02:47		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87				90					
*02:58			75			78					83				87									
*02:68		74	75			78					83				87		89							
*02:87		74									82													
*02:94		74	75			78	79				83				87									

Negative Control

*02:101	74	75			w						83			87														
*02:107	74										83																	
*02:115	74	75			78						83	85				90												
*02:180	74	75			78						83			87		90	91											
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96				
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:02:14-03:02:19, 03:02:21-03:02:22, 03:05, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482,			75	76																							94	
*03:02:10, 03:264, B*44:274 ³				76																								
*03:02:13, 03:02:20, 03:16, 03:110, 03:332, B*46:77 ³			75																								94	
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:57, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-				76	77																						93	94
*03:03:21, 03:03:51, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195, 03:344:01-03:344:01:02, 03:454, 03:459, 03:504				76																							93	
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451					77																						93	94
																											Negative Control	

*07:01:01-01-07:01:02:08, 07:01:04-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:02:01-01- 07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:06:03, 07:08, 07:11- 07:25, 07:27-01-07:27:02, 07:29:01-01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50- 07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149- 07:151, 07:153, 07:155, 07:157-07:160, 07:162- 07:172:02, 07:174-07:176, 07:178-07:183, 07:185- 07:195, 07:197-07:209, 07:211-07:215, 07:217- 07:222, 07:224-07:237, 07:239-07:241, 07:243- 07:245, 07:248-07:266, 07:268-07:292, 07:297- 07:313, 07:318-07:321, 07:323-07:327, 07:329N- 07:334, 07:336-07:337, 07:339-07:350N, 07:352- 07:355, 07:357-07:366, 07:368:01-07:401, 07:404, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520- 07:530, 07:532-07:535, 07:537-07:562, 07:564- 07:567, 07:569-07:577, 07:580-07:582Q, 07:584- 07:597, 07:599-07:625, 07:627-07:713, 07:715- *07:01:41, 07:02:41, 07:223, 07:415-07:416												81							95						
*07:01:74			78									81			85				89	90				w	95
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:26:01-01- 07:26:03, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405, 07:449, 07:493, 07:517, 07:714, 07:724, 07:821N													81												
*07:02:09			78										81		85				90						95
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598													81	83						92					95
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568													81								93				95
*07:31:01-07:31:02, 07:514													81		86				90						95
*07:39-07:40, 07:335													81							92					
*07:53, 07:216													81	83											95
*07:64, 07:402													81			88									94
*07:96:01-07:96:02, 07:314:02, 07:578, 07:756			75										81												94
*07:101, 07:148, 07:161			74										81	82											95
*07:123, 07:173, 07:294, 07:626, 07:725, 07:775													81						89						95
*07:125, 07:356, 07:531							78						81	83	85	86									95
*07:154													81						90						95
*07:177													81		86				90		92				
*07:210, 07:238													81	83							92				
*07:242						76							81												95
*07:293, 07:519													81			88									95
*07:316, 07:338, 07:579			73										81												95
*07:563													81								92				95
*07:583			74				78						81	82	85				90						
*07:723							78						81		85										
*08:01:01-01-08:01:17, 08:01:19-08:01:29, 08:03:01- 08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01- 08:16:02, 08:20, 08:22:01-01-08:22:01:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117- 08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186- 08:187, 08:189-08:190, 08:192-08:194, 08:196- 08:197, 08:199, 08:203																			90			93			

Negative Control

*12:208								80					85					89								
*12:215	73												84	85	86											
*12:233														85									93			
*12:269								78					84	85	86			88								
*12:287													84	85				89	90							
*14:02:01:01-14:02:07, 14:02:09-14:02:20, 14:02:22-14:02:31, 14:02:33-14:03:03, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:24:02, 14:26-14:27, 14:29-14:44, 14:46:01-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:62, 14:65-14:76, 14:78-14:86, 14:88-14:89, 14:91, 14:93N-14:104, 14:106-14:107, 14:109, *14:02:08, 14:02:21																		88				92				
*14:02:32, 14:06, 14:08-14:09, 14:28:01-14:28:02, 14:53, 14:63, 14:105Q, A*29:86, A*30:62, A*33:63, B*08:249, B*18:64, B*35:183, B*35:252, B*39:114, B*39:140, B*40:138, B*41:59 ³																		88								
*14:04, 14:64													83					88				92				
*14:10, 14:116																		88		90		92				
*14:12, 14:49, 14:108													83									92				
*14:15, 14:87															86			88		90						
*14:16															86			88								
*14:17																		88		90	91	92				
*14:25	74																	88				92		94		
*14:45																		87	88			92				
*14:54																		88						93		
*14:77													83					88								
*14:90													83	85				88							95	
*14:92																		88	89							
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:20, 15:02:22-15:02:24, 15:02:26-15:02:44, 15:06:01-15:06:03, 15:08:01-15:09:01:02, 15:10:02-15:10:03, 15:13:01:01-15:13:01:02, 15:15, 15:17-15:19, 15:26-15:28, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44:01-15:45, 15:47-15:53, 15:56-15:57, 15:60, 15:64, 15:67-15:68, 15:73-15:74, 15:78:01-15:84Q, 15:86-15:89, 15:91-15:101, 15:103, 15:106, 15:109, 15:112-15:114, 15:119-15:124, 15:126-15:137, 15:139, 15:141-15:143, 15:145N-15:146, 15:150-15:151, 15:154-15:157, 15:159-15:161, 15:163, 15:166, 15:168, 15:170-15:176, 15:179-15:180, 15:182-15:184, 15:186-15:189N, 15:191-15:200, 15:202-15:203, 15:205-15:207			75					80				83					87		89				93			
*15:02:10, 15:02:17			75					80				83					87		89				93	94		
*15:02:21, 15:02:25, 15:10:01, 15:10:04, 15:61								80				83					87		89				93			
*15:03, 15:16			75					80				83	85				87		89				93			
*15:04:01:01-15:04:03, 15:24, 15:30, 15:40, 15:66, 15:76-15:77, 15:147, 15:177N			75									83					87		89				93			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		
*15:05:01:01-15:05:08, 15:05:10, 15:05:13-15:05:14, 15:05:16, 15:29, 15:46, 15:54, 15:59, 15:69-15:70, 15:90, 15:105Q, 15:108, 15:110, 15:115N, 15:117, 15:140, 15:148, 15:152:01-15:153, 15:164N-15:165, 15:167, 15:185N, 15:204			75					w				83				87		89				93				
*15:05:09, 15:05:15								w				83				87		89				93				
*15:05:11-15:05:12, 15:116			75					w								87		89				93				
*15:07:01:01-15:07:01:02			75					80										89				93				
*15:11, 15:118			75			78		80				83				87		89				93				
*15:12					77			80				83				87		89	90			93				
*15:21, 15:71			75					80								87		89				93				
*15:22, 15:125			75					w				83				87		89								
*15:23:01-15:23:02								w				83						89				93				
*15:25			75										85					90	91		93					
*15:31, 15:62			75					80				83				87		89	90		93					
*15:36						78		w				83						89				93				
*15:37, 15:55, 15:58, 15:72, 15:149			75					80				83				87		89								
*15:42	74		75					80				83				87		89				93				
*15:43			75					80										89				93	94			

Negative Control

A-B-C 2020 9K8 C

*15:63, 15:138, 15:158								80			83				89					93				
*15:65			75								83			87	89	90	91							
*15:75			75					80			83	85		87	89	90					93	95		
*15:85, 15:162			75					80			83				89						93			
*15:102			75								83			87	89					92				
*15:104	73		75					w			83			87	89						93			
*15:107							79	80			83			87	89						93			
*15:111			74	75				w			83			87	89						93			
*15:144						77		80							89	90					93			
*15:169				75				w			83			87	89	90					93			
*15:178				75			79	80			83			87	89						93			
*15:181				75				w			83				89						93			
*15:190				75							83			87	89									
*15:201			74	75				w			83			87	89		91				93			
*16:01:01:01-16:01:15, 16:01:18-16:01:36, 16:06-16:08, 16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98, 16:100, 16:105-16:106, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164											83	85			90							95		
*16:01:16								80			83	85			90							95		
*16:01:17, 16:15:02											83	84	85		90							95		
*16:02:01:01-16:02:12, 16:02:14-16:02:17, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163											83	85		87	90							95		
*16:02:13											83	84	85	87	90							95		
*16:04:01:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149													86		90		92							
*16:10, 16:67												85			90							95		
*16:20								80			83	85			89	90						95		
*16:29			74										86		90		92							
*16:34			75				78				83	85			90							95		
*16:35													86		90									
*16:40, 16:110							79							88										
*16:43											83	85										95		
*16:45												85			90		93					95		
*16:48											83		86	87	90									
*16:50			74								83	85			90							95		
*16:53							79							88					92					
*16:64, 16:87											83	85			89	90						95		
*16:68														90		92								
*16:70			75								83	85	87	89	90							95		
*16:88											83		87		90		92							
*16:91											83		86	87	90		92							
*16:103											83		87		90									
*16:109								80					86		89	90	92							
*16:113							79				83	85		88								95		
*16:121			75				78				83	85	87		90							95		
*16:139											82	83	85		90							95		
*16:140											83		87											
*16:150							78						86	88	90	92								
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:51:01:02			74								83			87			91				93			
*17:01:04			74								83	85	87				91				93			
*17:07, 17:17, 17:23, 17:33											83		87				91				93			
*17:20			74								83						91				93		95	
*17:22			74														91				93			

Negative Control

*18:01:01:01-18:02:04, 18:04-18:08, 18:10-18:13											83									92	93	95	Negative Control	
*18:03	74										83									92	93	94		95
*18:09											83		87							92	93			95
*18:14											83									92				95
A*24:64, A*24:96, A*24:106, B*08:16, B*08:123, B*08:145																								95
A*24:174											83													95
B*13:13:02, B*40:01:49-40:01:50																						94		
B*15:125:02, B*44:02:44, B*56:01:09, B*58:74	75																							
B*15:277, B*15:430, B*35:222, B*42:22, B*54:38																					87			
B*15:539, B*39:76, B*54:18	73																							
B*58:02:01:01																								80
B*73:01:01:01																						w		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94		95