

*25:30	-								8	9				13				17			21				
*25:36, 25:50, 25:52, 25:65	-								8	9															
*25:43	-			4					8	9															
*25:46	-													13									24		
*25:59	-	1		4								10									20				
*26:01:01-26:01:47, 26:01:49-26:01:62, 26:01:64-26:01:65, 26:08:01:01-26:08:03, 26:10-26:15, 26:17-26:18, 26:20, 26:23-26:28, 26:31-26:43:02, 26:45-26:48, 26:50-26:53, 26:55-26:56, 26:58-26:63, 26:66-26:71N, 26:73-26:74, 26:76-26:77, 26:79-26:88, 26:90, 26:93-26:110, 26:112-26:119, 26:121, 26:123-26:151, 26:153-26:176, 26:178-26:183, 26:185-26:192, 26:194-26:197, 26:199-26:203, 26:206N	A26, A10, Null, -								8	10			13												
*26:01:48, 26:64, 26:89	-											10		13											
*26:01:63, 26:57, 66:15 ^d	-								8					13											
*26:02:01-26:02:02, 26:29, 26:49	A26, -								8	10			w												
*26:03:01-26:03:01:02, 26:06, 26:21, 26:111	A26, -								8		11		13									23			
*26:04, 26:75, 26:205	A26, -								8	10															
*26:05	A26								8					13								23			
*26:07:01-26:07:02	A26											10		13								24			
*26:09, 26:91	A26, -								8	10	12														
*26:16	A26							7				10		13											
*26:19	-			4								10		13	14										
*26:22	A26								8	10				13	14						21				
*26:30	A26								8					13								23			
*26:54	-								8	10	12		13						20						
*26:65	-								8	10				13								24			
*26:72, 26:184	-			4								10		13											
*26:78	-								8		11											23			
*26:92, 66:42 ^d	-											11		13								24			
*26:120	-	1							8	10				13											
*26:122	-								8	10				13					20						
*26:152, 26:193	-								8	10	11			13											
*26:177	-											11		13				18		21		23			
*26:198	-								8	10				13	14										
*26:204	-			4					8	10															
*26:207	-											10		13				18							
*29:01:01-29:06, 29:08N-29:12, 29:15-29:27, 29:29-29:38, 29:40-29:47, 29:50-29:66, 29:68-29:79, 29:81-29:97, 29:99-29:104, 29:106-29:146	A29, Null, -													14											
*29:07, 29:49	A29, -					5	6							14											
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*29:13, 29:39, 29:98	-														14			17							
*29:14	-														14		16								
*29:28	-								8						14										
*29:48	-														14				18						
*29:67	-				4										14										
*29:80	-														14							21			
*29:105	-										10				14				18						
*30:01:01-30:04:03, 30:06-30:07, 30:09:01-30:20, 30:22-30:54, 30:56-30:88, 30:90-30:124, 30:126-30:139, 30:141-30:149, 30:151-30:162, 30:164-30:176, 30:178N-30:181, B*07:260, C*12:328 ^d	A30, Null, -															15									
*30:08:01-30:08:02, 30:150	A30, -				4											15									
*30:55	-													13		15									
*30:89, 30:177	-														15			17							
*30:125, 30:140	-											11			15										
*30:163	-							7							15										
*31:01:02:01-31:02:02, 31:05-31:07, 31:09-31:20, 31:22-31:23, 31:26-31:28, 31:30-31:34, 31:36-31:40, 31:42-31:61, 31:63-31:84, 31:86-31:108, 31:110-31:122, 31:124-31:175, 31:177-31:178, 31:180-31:184N	A31, Null, -																	16							

*66:02-66:03:01:03, 66:16, 66:21, 66:25-66:26Q, 66:28N, 66:34, 66:39N, 66:43	A66, A10, Null, -						8				12											24			
*66:06	A66						8		11		13											21			
*66:09	-						8		11		13	14										21			
*66:10	-						8		11		w														
*66:12	-						8				13											24			
*66:19	-						8		11	12	13														
*68:01:01:01-68:04:02, 68:06-68:08:01:02, 68:09-68:10, 68:12:01-68:12:02, 68:13:02-68:14, 68:16-68:18N, 68:21:01-68:24, 68:26-68:28, 68:30-68:40, 68:42-68:44, 68:46-68:65, 68:67-68:70, 68:72-68:83, 68:85-68:102, 68:104:01-68:104:02, 68:106-68:108, 68:110-68:116, 68:118-68:129, 68:131-68:133, 68:135, 68:137-68:156, 68:158-68:177, 68:179-68:217, 68:219-68:226, 68:228-68:229, 68:231-68:248, 68:251N	A68, A28, Null, -						8														20				
*68:05, 68:15, 68:20, 68:109, 68:136	A68, -						8														20	23			
*68:08:02, 68:105, 68:157	A68, -						8														20	21			
*68:11N, 68:41, 68:230	Null, -						8														20	22			
*68:13:01, 68:66:01-68:66:02, 68:134	A68, -				4		8														20				
*68:19, 68:130:02	A68, -						8				14										20				
*68:25	A68						8											19	20						
*68:29	A68						8										18	20	21						
*68:45, 68:117	-					7	8														20				
*68:71	-						8				13										20				
*68:84, 68:227	-						8		10												20				
*68:103:01-68:103:02	-				3		8														20				
*68:130:01	-						8				12										20				
*68:178	-				4		8														20				
*68:218	-						8																		
*68:249	-					5	8														20				
*68:250	-																18	20	21						
*69:01:01:01-69:01:03, 69:03-69:06, 69:08	A69, -						8															21			
*69:02	-						8		11													21			
*69:07	-						8				13											21			
*74:01:01:01-74:01:02, 74:01:04-74:03:02, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:29, 74:31-74:32N, 74:34, 74:36, 74:38-74:39	A74, Null, -																				19				
*74:01:03	-										12										19	24			
*74:04, 74:21	A74, -																	18	19						
*74:07, 74:33	A74, -													17							19				
*74:13	-										13		16								19				
*74:19	-				4																19				
*74:23	-				3																19				
*74:30	-									11											19				
*74:35	-																				19	21			
*74:37	-										13										19				
*80:01:01:01-80:01:02, 80:03, 80:05-80:07	A80, -					6																23			
*80:02	-					6				11															
*80:04	-					6												17				23			
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*49:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:01:38, B*57:60	-, B53								9																
B*18:27	-					6																			
B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569	-																					22			
C*03:365, C*12:131	-				4																				
C*06:187	-																					23			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24