

A-B-C 2021 OM2 C

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele <sup>1,2</sup>																								
*01:02:01-01:02:17, 01:02:19-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72, 01:05-01:08, 01:11-01:13, 01:15-01-01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01-01:89N, 01:91-01:100, 01:103-01:112, 01:114-01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:199	73																			92				
*01:02:18	73											84									92			
*01:02:49, 01:136, 01:145:01N, 01:146:01-01:146:02	73							80													92			
*01:02:68, 01:02:71, 01:168-01:169:01	73														87						92			
*01:03:01-01:03:02, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:176	73																				92	93		
*01:04	73											84		86							92			
*01:09:01-01:09:02	73													86							92			
*01:10	73	74																			92			
*01:14, 01:59, 01:118, 01:186	73											83									92			
*01:21	73											84	85											
*01:43	73	74										82									92			
*01:60	73																			91	92			
*01:67, 01:122	73			76													88				92			
*01:90, 01:101-01:102, 01:113	73							80										89			92			
*01:116	73																	89			92			
*01:125, 14:58 <sup>3</sup>	73																			88				
*01:131	73			76																	92	93		
*01:160	73													86					90		92			
*01:169:02, 03:94 <sup>3</sup>	73		75												87						92		94	
*01:177	73															88					92		94	
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:31-02:02:37, 02:02:39, 02:02:41-02:02:43, 02:02:45-02:02:52, 02:02:54-02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:110, 02:112-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186, 02:188, 02:190, 02:192N-02:196		74	75			78					83				87									
*02:02:13, 02:02:29		74	75								83				87								94	
*02:02:26, 02:21, 02:26:02-02:26:03, 02:39, 02:45, 02:189		74	75								83				87									
*02:02:30, 02:02:40		74				78					83				87									
*02:02:38, 02:02:53		74	75			78					83				87							94		
*02:02:44		74	75			78					83	84			87									
*02:04		74	w			78					83				87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87		89							
*02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01		74	w			78					83													
*02:14:02, 02:164		74				78					83													
*02:22, 02:62, 02:118		74	75			78					83				87							93		
*02:26:01, 02:134, 02:152, 02:167, 02:191		74									83				87									
*02:27:01-02:27:02, 02:65, 02:131, 02:187		74	75			78					83													
*02:40:01-02:40:02		74	75			78					83				87					91				
*02:47, 02:111		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87				90					
*02:58			75			78					83				87									
*02:68		74	75			78					83				87		89							
*02:87		74									82													
*02:94		74	75			78	79				83				87									
*02:101		74	75			w					83				87									
*02:107		74									83													
*02:115		74	75			78					83		85						90					
*02:180		74	75			78					83				87				90	91				
*02:197			75			78					83				87							93		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

A-B-C 2021 OM2 C

<p>*03:02:01-03:02:09, 03:02:11-03:02:25, 03:05, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505, 03:532</p>			75	76							94	
<p>*03:02:10, 03:264, B*15:96, B*44:274<sup>3</sup></p>			76									
<p>*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:59, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495, 03:503, 03:506-03:507, 03:509N, 03:514-03:517, 03:519-03:520, 03:524, 03:527-03:529, 03:533, 03:537, 03:541-03:542</p>			76	77					93	94		
<p>*03:03:21, 03:03:51, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195, 03:344:01:01-03:344:01:02, 03:454, 03:459, 03:504</p>			76							93		
<p>*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451, 03:534-03:535, 03:540</p>				77						93	94	

Negative Control







*06:118								80			83			86				90	91			95		
*06:132:02			75								83			86	87							94	95	
*06:204	73							80			83									92			95	
*07:01:01-07:01:02:08, 07:01:04-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:01:94-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:06:04, 07:08, 07:11-07:25, 07:27:01-07:27:02, 07:29:01:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:132, 07:133:02-07:138, 07:140:01-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:215, 07:217-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:401, 07:404, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:530, 07:532-07:535, 07:537-07:562, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:713Q, 07:715:01-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844, 07:846-07:849N, 07:851-07:856N, 07:858-07:869, 07:871-07:890, 07:893, 07:895-07:897								81												95				
*07:01:41, 07:02:41, 07:223, 07:415-07:416																						w		
*07:01:74								78			81			85				89	90			95		
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:26:01:01-07:26:03, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405, 07:449, 07:493, 07:517, 07:714, 07:724, 07:821N, 07:845, 07:870, 07:891, 07:894											81													
*07:02:09								78			81			85				90				95		
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598											81	83								92		95		
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568, 07:857, 07:892											81										93	95		
*07:31:01-07:31:02, 07:514											81			86				90				95		
*07:39-07:40, 07:335											81									92				
*07:53, 07:216											81	83										95		
*07:64, 07:402											81					88						94		
*07:96:01-07:96:02, 07:314:02, 07:578, 07:756			75								81											94		
*07:101, 07:148, 07:161			74								81	82										95		
*07:123, 07:173, 07:294, 07:626, 07:725, 07:775, 07:850											81							89				95		
*07:125, 07:356, 07:531								78			81	83		85	86							95		
*07:133:01, 07:242, 07:330:01, 07:748				76							81											95		
*07:154											81							90				95		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*07:177									81					86				90		92				
*07:210, 07:238									81		83									92				
*07:293, 07:519									81							88							95	
*07:316, 07:338, 07:579	73								81														95	
*07:563									81											92			95	
*07:583		74					78		81	82			85				90							
*07:723							78		81				85											
*08:01:01-08:01:17, 08:01:19-08:01:30, 08:03:01-08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22:01-08:22:01:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186-08:187, 08:189-08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203-08:204, 08:208N-08:210, 08:212, 08:217											82							90			93			

Negative Control



A-B-C 2021 OM2 C

*12:03:01-12:03:01:38, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:17-12:03:19, 12:03:21-12:03:22, 12:03:25, 12:03:27-12:03:29, 12:03:31-12:03:33, 12:03:35-12:03:46, 12:03:48-12:03:65, 12:06-12:07, 12:11-12:13:01:02, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57:01-12:57:02, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82, 12:87-12:91, 12:93-12:95, 12:97, 12:99:01-12:102, 12:107-12:111, 12:115-12:116, 12:120-12:122, 12:125, 12:129, 12:131, 12:133, 12:138-12:140:01, 12:141, 12:143, 12:149-12:150, 12:152, 12:156-12:160, 12:163, 12:165, 12:167, 12:170-12:174, 12:180, 12:182, 12:184, 12:186-12:187, 12:189-12:192, 12:197, 12:199, 12:201-12:202, 12:205-12:206, 12:209-12:211, 12:213, 12:216, 12:218, 12:220, 12:223, 12:225, 12:227, 12:229-12:230, 12:232N, 12:235, 12:237-12:238, 12:242, 12:244-12:246, 12:248-12:249, 12:253-12:254, 12:256-12:257, 12:259-12:260:02, 12:262, 12:264-12:267, 12:271-12:274:01N, 12:276-12:278, 12:283-12:284, 12:286, 12:288-12:293, 12:295N, 12:297, 12:300, 12:302, 12:305-12:306, 12:310-12:314, 12:316-12:324N, 12:326, 12:328					84	85	86									
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:03:30, 12:26, 12:59, 12:185								85	86							
*12:03:04, 12:98, 12:301						84		86								
*12:03:09			80			84		86								
*12:03:13, 12:03:24						84	w	86								
*12:03:23, 12:140:02	75					84		85	86							
*12:03:26, 12:03:66			80			84		85	86							
*12:03:34:01-12:03:34:02, B*15:510, B*35:475 <sup>4</sup>									86							
*12:04:01						83			86	87						
*12:04:02:01-12:04:02:02, 12:60						83	84	85	86	87						
*12:05:01-12:05:02						83	84		86	87						
*12:08							84	85			89					
*12:09						83				87				92		
*12:14:01								85			90	91				
*12:14:02							84	85			90	91				
*12:15, 12:113			80			84	85	86			89					
*12:16:01, 12:217						84	85									95
*12:16:02						84										95
*12:21			79			83	84			87						
*12:24							84	85						92		
*12:28, 12:58, 12:325							84	85	86		90					
*12:31:01-12:31:02							84	85	86						93	
*12:33			79			83			86	87						
*12:41						83	84	85		87						
*12:44							84	85							93	
*12:54, 12:194						83	84	85	86							
*12:55							84	85	86			88				
*12:72							84	85		87						
*12:81, 12:282							84	85	86		89					
*12:92							84	85	86			91				
*12:119	74						84	85	86							
*12:127						82		84	85							
*12:132, 12:169, 12:252							84	85			90					
*12:135							84	85	86	87	90					
*12:144								85	86						93	
*12:146						83	84	85		87	90					
*12:147, 12:195:02							84	85	86							95
*12:154, 12:178							84	85	86	87						
*12:176							84	85	86		90	91	92			
*12:181					81			85			90					
*12:188						83	84				89					
*12:195:01							84									95
*12:195:03									86							95
*12:198	74						84	85								
*12:203						82	84	85	86							
*12:208					80			85			89					
*12:215	73						84	85	86							
*12:233								85							93	
*12:269					78		84	85	86		88					
*12:287							84	85			89	90				
*12:298					80				85							

Negative Control



A-B-C 2021 OM2 C

*12:309					80				85				89				93							
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*14:02:01:01-14:02:07, 14:02:09-14:02:20, 14:02:22-14:02:31, 14:02:33-14:03:04, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:24:02, 14:27, 14:29-14:44, 14:46:01-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:62, 14:65-14:76, 14:78-14:86, 14:88-14:89, 14:91, 14:93N-14:104, 14:106-14:107, 14:109, 14:111-14:115, 14:117N-14:120, 14:122-14:123																88				92				
*14:02:08, 14:02:21, 14:03:05												84				88				92				
*14:02:32, 14:06, 14:08-14:09, 14:28:01-14:28:02, 14:53, 14:63, 14:105Q, A*29:86, A*30:62, A*33:63, B*08:249, B*18:64, B*35:183, B*35:252, B*39:114, B*39:140, B*40:138, B*41:59 <sup>4</sup>																88								
*14:04, 14:64											83					88					92			
*14:10, 14:116																88			90		92			
*14:12, 14:49, 14:108											83										92			
*14:15, 14:87														86		88			90					
*14:16														86		88								
*14:17																88			90	91	92			
*14:25		74														88					92		94	
*14:26								80								88					92			
*14:45		73														87	88				92			
*14:54																88						93		
*14:77											83					88								
*14:90											83		85			88								95
*14:92																88	89							
*14:110, B*07:239, B*14:03 <sup>4</sup>																					92			
*14:121									81							88					92			
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:20, 15:02:22-15:02:24, 15:02:26-15:02:45, 15:05:01:01-15:05:07, 15:05:10, 15:05:13-15:05:14, 15:05:16-15:06:03, 15:08:01-15:09:02, 15:10:02-15:10:03, 15:13:01:01-15:13:01:02, 15:15, 15:17-15:19, 15:26-15:29, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44:01-15:54, 15:56-15:57, 15:59-15:60, 15:64, 15:67-15:70, 15:73-15:74, 15:78:01-15:84Q, 15:86-15:101, 15:103, 15:105Q-15:106, 15:108-15:110, 15:112-15:115N, 15:117, 15:119-15:124, 15:126-15:137, 15:139-15:143, 15:145N-15:148, 15:150-15:157, 15:159-15:161, 15:163-15:168, 15:170-15:176, 15:179-15:180, 15:182-15:189N, 15:191-15:200, 15:202-15:207, 15:209-15:210, 15:212-15:217, 15:219-15:221, 15:223-15:224, 15:226-15:228			75					80			83					87	89					93		
*15:02:10, 15:02:17			75					80			83					87	89				93	94		
*15:02:21, 15:02:25, 15:05:09, 15:05:15, 15:10:01, 15:10:04, 15:61, 15:208								80			83					87	89				93			
*15:03, 15:16			75					80			83		85			87	89				93			
*15:04:01:01-15:04:05, 15:24, 15:30, 15:40, 15:66, 15:76-15:77, 15:177N, 15:211, 15:222			75								83					87	89				93			
*15:05:11-15:05:12, 15:21, 15:71, 15:116			75					80								87	89				93			
*15:07:01:01-15:07:01:02			75					80									89				93			
*15:11, 15:118			75				78	80			83					87	89				93			
*15:12						77		80			83					87	89	90			93			
*15:22, 15:37, 15:55, 15:58, 15:72, 15:125, 15:149, 15:225			75					80			83					87	89							
*15:23:01-15:23:02, 15:63, 15:138, 15:158								80			83						89					93		
*15:25			75										85					90	91		93			
*15:31, 15:62, 15:169			75					80			83					87	89	90			93			
*15:36						78		80			83						89				93			
*15:42, 15:111			74	75				80			83					87	89				93			
*15:43				75				80									89				93	94		
*15:65				75							83					87	89	90	91					
*15:75				75				80			83		85			87	89	90			93		95	
*15:85, 15:162, 15:181				75				80			83						89				93			
*15:102		73	75								83					87	89				92			
*15:104		73	75					80			83					87	89					93		
*15:107							79	80			83					87	89					93		
*15:144					77			80									89	90				93		
*15:178			75				79	80			83					87	89					93		
*15:190			75								83					87	89							
*15:201			74	75				80			83					87	89		91		93			
*15:218				75	76			80			83					87	89				93			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

A-B-C 2021 OM2 C

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*16:01:01:01-16:01:15, 16:01:18-16:01:39, 16:06-16:08, 16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98, 16:100, 16:105-16:106, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175											83		85					90						95
*16:01:16								80			83		85					90						95
*16:01:17, 16:15:02											83	84	85					90						95
*16:02:01:01-16:02:12, 16:02:14-16:02:18, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163, 16:166-16:167, 16:176											83		85		87			90						95
*16:02:13											83	84	85		87			90						95
*16:04:01:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149														86				90		92				
*16:10, 16:67													85					90						95
*16:20								80			83		85				89	90						95
*16:29		74													86			90		92				
*16:34			75			78					83		85					90						95
*16:35														86				90						
*16:40, 16:110							79									88								
*16:43											83		85											95
*16:45													85					90			93			95
*16:48											83			86	87			90						
*16:50		74									83		85					90						95
*16:53							79									88				92				
*16:64, 16:87											83		85				89	90						95
*16:68																		90		92				
*16:70			75								83		85		87		89	90						95
*16:88											83				87			90		92				
*16:91											83			86	87			90		92				
*16:103											83				87			90						
*16:109								80						86			89	90		92				
*16:113							79				83		85			88								95
*16:121			75			78					83		85		87			90						95
*16:139										82	83		85					90						95
*16:140											83				87									
*16:150						78								86		88		90		92				
*16:169								80					85				89	90						95
*16:170											83		85			88		90						95
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:57		74									83				87			91		93				
*17:01:04		74									83		85		87			91		93				
*17:07, 17:17, 17:23, 17:33											83				87			91		93				
*17:20		74									83							91		93				95
*17:22		74																91		93				
*18:01:01:01-18:02:05, 18:04-18:08, 18:10-18:13											83									92	93			95
*18:03		74									83									92	93	94		95
*18:09											83				87					92	93			
*18:14											83									92				95
A*02:42:01, B*15:539, B*39:76, B*54:18	73																							
A*24:64, A*24:96, A*24:106, B*08:16, B*08:123, B*08:145																								95
A*24:174, B*57:129											83													95
B*13:13:02, B*40:01:49-40:01:50																							94	
B*15:125:02, B*44:02:44, B*56:01:09, B*58:74			75																					
B*15:277, B*15:430, B*35:222, B*42:22, B*54:38																87								
B*46:84															86	87								
B*58:02:01:01-58:02:02								80																
B*73:01:01:01-73:02																					w			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control