

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01:01-01:02:17, 01:02:19-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72-01:02:92, 01:05-01:08, 01:11-01:13, 01:15-01-01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:47, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01-01:89N, 01:91-01:100, 01:103-01:112, 01:114-01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:214, 01:216-01:217:01:02, 01:219-01:232, 01:234-01:243, 01:245, 01:247-01:248	73																			92				
*01:02:18	73											84									92			
*01:02:49, 01:136, 01:145:01N, 01:146-01-01:146:02, 01:233	73							80													92			
*01:02:68, 01:02:71, 01:168-01:169:01	73														87						92			
*01:03:01-01:03:03, 01:24, 01:34, 01:49-01-01:50, 01:55, 01:78, 01:176	73																				92	93		
*01:04	73											84		86							92			
*01:09:01-01:09:02, 01:244	73													86							92			
*01:10	73	74																			92			
*01:14, 01:59, 01:118, 01:186	73											83									92			
*01:21	73											84	85											
*01:43	73	74										82									92			
*01:48	73		75																		92			
*01:60	73																				91	92		
*01:67, 01:122	73			76													88				92			
*01:90, 01:101-01:102, 01:113	73							80										89			92			
*01:116	73																	89			92			
*01:125, 01:218, 14:58 ³	73																88				92			
*01:131	73			76																	92	93		
*01:160	73													86				90			92			
*01:169:02, 03:94 ³	73	75													87						92	94		
*01:177	73															88					92	94		
*01:215	73												85								92		95	
*01:246	73																				92	94		
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:30-02:02:37, 02:02:39-02:02:43, 02:02:45-02:02:52, 02:02:54-02:04, 02:07-02:11, 02:15-02:16:02:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:110, 02:112-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186, 02:188, 02:190, 02:192N-02:196, 02:198-02:220		74	75			78					83				87									94
*02:02:13, 02:02:29		74	75								83				87									94
*02:02:26, 02:21, 02:26:01-02:26:03, 02:39, 02:45, 02:189, 02:191		74	75								83				87									94
*02:02:38, 02:02:53		74	75			78					83				87								94	
*02:02:44		74	75			78					83	84			87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87			89						
*02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01-02:14:02, 02:27:01-02:27:02, 02:65, 02:131, 02:164, 02:187		74	75			78					83													
*02:22, 02:62, 02:118		74	75			78					83				87							93		
*02:40:01-02:40:02		74	75			78					83				87				91					
*02:47, 02:111		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87			90						
*02:58			75			78					83				87									
*02:68		74	75			78					83				87			89						
*02:87		74									82													
*02:94		74	75			78	79				83				87									
*02:101		74	75			w					83				87									
*02:107		74									83													
*02:115		74	75			78					83		85						90					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*02:134, 02:152, 02:167		74									83				87									
*02:180		74	75			78					83				87			90	91					
*02:197			75			78					83				87							93		

Negative Control

*04:100						78	80				83	85				90			93					
*04:103:01-04:103:02											83	85	87			90			93					
*04:107		75				78					83	85	87						93					
*04:108						78	80				83	85				89	90							
*04:112						78					83	85				89	90			93				
*04:129						78	79				83	85					90			93				
*04:160:01-04:160:02		75				78					83	85	87			90	91			93				
*04:169						78					83					89				93				
*04:172						78					83	85		88						93				
*04:178, 04:422						78	80				83	85				89	90			93				
*04:201:02, 04:212			76			78					83	85					90			93				
*04:230, 04:263, 04:397, 04:432, 04:442						78					83	85					90	92		93				
*04:279N						78					83	w				90				93				
*04:299		75				78					83		87			90				93				
*04:335		75				78	80				83	85	87			90								
*04:357		74	75			78					83	85	87			90				93				
*04:360, 04:477, 04:496						78					83	85	86			90								
*04:383, 04:473			75			78					83	85				90				93				
*04:406		74				78					83	85				90				93				
*04:414						78					83									93				
*04:447		74				78					83	85				90	91		93	94				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*05:01:01:01-05:01:50, 05:01:52-05:01:55, 05:01:57, 05:01:59-05:01:74, 05:03, 05:05:01-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:59-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106:01-05:106:02, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:142, 05:144-05:145, 05:148-05:150, 05:152-05:155, 05:157-05:162, 05:164-05:190, 05:192-05:199, 05:201, 05:203-05:205, 05:208N-05:214, 05:216:01:01-05:237, 05:239N-05:240, 05:242-05:244N, 05:246-05:251N, 05:253N-05:271N, 05:273-05:280							79				83			87										93
*05:01:51, 05:01:58, 05:04:01-05:04:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147, 05:191, 05:207, 05:272							79				83			87										
*05:01:56							79	80			83			87										93
*05:16							79				83			87					92					93
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115							79				83			87			90							93
*05:20, 05:43, 05:113N, 05:156, 05:202Q							79				83													93
*05:22:01-05:22:02, 05:252							79	80			83			87										
*05:23, 05:62, 05:146, 05:151, 05:163							79				83			87			91							93
*05:26							78	79			83			87										93
*05:36							79				83			87		89								93
*05:42							79				83	85	87			90								93
*05:46:01:01-05:46:01:02							79				83	85	87											93
*05:58:01-05:58:04, 05:281			75				79				83			87										93
*05:78:01-05:78:02							78	79			83													93
*05:85	73						79				83			87					92					93
*05:105		74					79				83			87										93
*05:107	73						79				83			87					92					
*05:112							79				83	85	87											93
*05:129							79				83			87										93
*05:132, 05:215							79	80			83			87		89								
*05:134, 05:143							79				83			87			90	91						93
*05:200	73	74					79				83													93
*05:206, 17:05 ³		74					79				83			87			91							93
*05:238							79				83			87		89								
*05:241	73						79	80			83			87						92				
*06:02:01:01-06:02:01:94, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73-06:02:101, 06:02:103-06:02:108, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34:02, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:69, 06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:179, 06:181-06:196, 06:198-06:202, 06:205-06:216, 06:218-06:223, 06:225-06:235, 06:237-06:243, 06:245-06:247, 06:249-06:251, 06:253-06:263N, 06:265-06:270, 06:272-06:279, 06:281N-06:303, 06:306-06:316N, 06:318-06:328, 06:330-06:334, 06:336Q-06:347N, 06:349-06:357N, 06:360-06:363								80			83		86											95
*06:02:12, 06:02:65, 06:11, 06:52, 06:122, 06:124, 06:147, 06:217, 06:248, 06:252, 06:280								80					86											95
*06:02:16, 06:02:43, 06:04:01-06:04:02:02, 06:40, 06:81, 06:101, 06:123, 06:127:01:01-06:127:02, 06:136, 06:153:01-06:153:02, 06:197, 06:203, 06:264, 06:304, 06:317, 06:358								80			83													95
*06:02:38:01-06:02:38:02								80			83	84		86										95

Negative Control

*06:02:39, 06:02:102, 06:32, 06:94-06:95, 06:99, 06:180, 06:236, 06:244, 06:329:01-01-06:329:01:02, 06:348											83		86												95
*06:02:72											83	85	86				90								95
*06:03:01, 06:132:01		75									83	86	87												94
*06:03:02, 06:76:01											80	83	86	87											
*06:05							79	80			83						90			92					95
*06:06:01:01-06:06:01:02, 06:79N, 06:271								80			83										92				95
*06:08		74						80			83	86													95
*06:09:01-06:09:02								80			83	86												93	95
*06:14, 06:35, 06:72, 06:143-06:144								80			83													93	95
*06:17											83	86											92		95
*06:31								80			83	86					90			92					95
*06:67							79	80			83	86													95
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	
*06:70:01-06:70:02, 06:305			75					80			83			86											95
*06:76:02											83	84	85	86	87		90								
*06:77								80			83			86											
*06:89, 06:224								80			83			86		89									95
*06:118								80			83			86			90	91							95
*06:132:02			75								83			86	87										94
*06:204	73							80			83													92	95
*06:335											83	85	86												95
*06:359N							79	80			83			86	88										95
*07:01:01:01-07:01:02:14, 07:01:04-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:01:94-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:06:06, 07:08, 07:11-07:19, 07:21-07:25, 07:27:01-07:27:02, 07:29:01:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126, 07:128-07:132, 07:133:02-07:138, 07:140:01-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:215, 07:217-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:389, 07:391-07:401, 07:404, 07:405:02, 07:407-07:414, 07:417-07:448, 07:450:01-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:530:02, 07:532-07:535, 07:537-07:561, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:713Q, 07:715:01-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844, 07:846:01-07:849N, 07:851-07:856N, 07:858-07:869, 07:871-07:890, 07:893, 07:895-07:931, 07:933-07:935, 07:937-07:939, 07:941-07:969, 07:971-07:974Q, 07:976-07:986:02, 07:988, 07:990-07:992, 07:994-07:1009, 07:1011-07:1027, 07:1029-07:1044N, 07:1046-07:1055, 07:1057-07:1058, 07:1060, 07:1062N-07:1065										81															95
*07:01:41, 07:02:41, 07:223, 07:415-07:416, 07:1045											81													w	
*07:01:74						78					81			85			89	90							95
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:26:01:01-07:26:03, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405:01, 07:449, 07:493, 07:517, 07:714, 07:724, 07:821N, 07:845, 07:870, 07:891, 07:894, 07:936N, 07:970N, 07:989, 07:1056, 07:1061											81														
*07:02:09						78					81			85				90							95
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598											81	83												92	95
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568, 07:857, 07:892, 07:987, 07:1028											81												93		95
*07:20:01-07:20:02, 07:127:01-07:127:02, 07:263, 07:390, 07:562			75								81														95
*07:31:01-07:31:02, 07:514:01-07:514:02											81			86				90							95
*07:39-07:40, 07:335, 07:975, 07:993:01-07:993:02, 07:1059											81													92	
*07:53, 07:216, 07:932											81	83													95
*07:64, 07:402											81					88									94
*07:96:01-07:96:02, 07:314:02, 07:578, 07:756			75								81														94
*07:101, 07:148, 07:161			74								81	82													95
*07:123:01-07:123:02, 07:173, 07:294, 07:626, 07:725, 07:775, 07:850											81							89							95
*07:125:01, 07:356, 07:531						78					81	83	85	86											95
*07:125:02						78					81	83	85	86				90							95
*07:133:01, 07:242, 07:330:01, 07:748				76							81														95
*07:154											81							90							95
*07:177											81			86				90		92					
*07:210, 07:238											81	83												92	
*07:293, 07:519											81						88								95

Negative Control

*07:316, 07:338, 07:579	73								81										95					
*07:563									81									92	95					
*07:583		74				78			81	82			85					90						
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*07:723						78			81				85											
*07:940									81			84												95
*07:1010						78			81	83			86											95
*08:01:01-01-08:01:07, 08:01:09-08:01:17, 08:01:19-08:01:34, 08:03:01-01-08:03:05, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22:01-01-08:22:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186-08:187, 08:189-08:190, 08:192-08:194, 08:196-08:197, 08:199, 08:203-08:204, 08:208N-08:210, 08:212, 08:217, 08:219-08:221, 08:228, 08:232, 08:234-08:236N, 08:241, 08:244, 08:249-08:250, 08:253-08:255, 08:259, 08:261										82								90				93		
*08:01:18, 08:02:01-01-08:02:37, 08:02:39, 08:04:01-01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19:02, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62:01-08:63, 08:66-08:71, 08:73, 08:75-08:77, 08:90, 08:92-08:94, 08:100, 08:103:01-01-08:103:01:02, 08:107, 08:110-08:112, 08:114-08:116, 08:120, 08:123, 08:125-08:126, 08:132, 08:134, 08:138-08:140, 08:146, 08:150-08:151, 08:156, 08:158-08:161N, 08:165-08:167, 08:169-08:172, 08:177, 08:179-08:182, 08:184, 08:191, 08:195, 08:198, 08:200-08:202, 08:206, 08:207:02, 08:211, 08:213-08:216, 08:218, 08:222-08:227, 08:229-08:231, 08:233, 08:237-08:238, 08:240, 08:242-08:243, 08:245, 08:247-08:248, 08:256-08:258, 08:260, 08:262-08:264										82													93	
*08:01:35, 08:105			75							82								90					93	
*08:02:38, 08:48, 08:152										82														
*08:05, 08:25										82			85											93
*08:07, 08:104, 08:142										82										91				93
*08:09, 08:11, 08:205										82										90				
*08:10							79			82	83				87				90					93
*08:12:01-01-08:12:01:02										82													92	93
*08:21, 08:137										82			85							90				93
*08:31, 08:246		74								82														93
*08:47, 08:188										82										90	91			93
*08:74										82					88									93
*08:83								80		82										90				
*08:95, 08:251										82					88					90				93
*08:108									80	82														
*08:113									80	82							89							
*08:163						78				82					88					90				93
*08:183						78				82					88									93
*08:185		74								82											91			93
*08:207:01										80														93
*08:239			75							82														93
*08:252		74				78			81	82										90	92			93
*12:02:01-12:02:10, 12:02:12-12:02:13, 12:02:16, 12:02:18, 12:02:20-12:02:27, 12:02:31, 12:02:33, 12:02:36, 12:02:38-12:02:39, 12:02:41-12:02:49, 12:03:47, 12:10:01-12:10:03, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:85:01, 12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N, 12:151, 12:153, 12:155Q, 12:161-12:162, 12:164, 12:166, 12:168, 12:175, 12:177, 12:179, 12:183, 12:193, 12:196, 12:200, 12:204, 12:207, 12:212, 12:214, 12:219N, 12:221, 12:226, 12:228, 12:231, 12:234, 12:236N, 12:240-12:241, 12:243, 12:247, 12:250-12:251, 12:255, 12:258, 12:261, 12:263, 12:268, 12:270N, 12:275, 12:279-12:281, 12:285, 12:294, 12:296, 12:303-12:304, 12:307-12:308, 12:315, 12:327N, 12:329N, 12:335-12:336, 12:338, 12:351Q, 12:359, 12:363, 12:383												84	85											
*12:02:11								80				84												
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*12:02:14-12:02:15, 12:02:19, 12:02:30, 12:02:32, 12:02:34, 12:02:37, 12:02:40, 12:18:02, 12:85:02, 12:222, 12:224, B*67:02:01-01-67:02:01:02 ³													85											
*12:02:17, 12:239			75									84	85											
*12:02:28												84	w											
*12:02:29, 12:49, 12:123, 12:381, B*18:01:42, B*57:01:24 ³												84												

Negative Control

*15:201		74	75					80			83			87	89	91	93							
*15:218			75	76				80			83			87	89		93							
*16:01:01:01-16:01:15, 16:01:18-16:01:44, 16:06-16:08, 16:11, 16:13-16:14, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49:01-16:49:02, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98, 16:100, 16:105-16:106, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:180, 16:183, 16:185-16:187, 16:189-16:193, 16:196-16:198, 16:200-16:201, 16:204-16:206N										83	85			90				95						
*16:01:16								80			83	85			90				95					
*16:01:17, 16:15:02										83	84	85			90				95					
*16:02:01:01-16:02:12, 16:02:14-16:02:20, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163, 16:166-16:167, 16:176, 16:179, 16:181, 16:184, 16:194-16:195N, 16:202										83	85	87		90					95					
*16:02:13										83	84	85	87		90				95					
*16:04:01:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149, 16:182, 16:203												85		90	92									
*16:10, 16:67, 16:188												85		90					95					
*16:20								80		83	85			89	90				95					
*16:29		74										86		90	92									
*16:34			75			78				83	85			90					95					
*16:35												86		90										
*16:40, 16:110						79								88										
*16:43										83	85								95					
*16:45											85			90			93		95					
*16:48										83		86	87	90										
*16:50		74								83	85			90					95					
*16:53						79								88			92							
*16:64, 16:87										83	85			89	90				95					
*16:68														90	92									
*16:70			75							83	85	87		89	90				95					
*16:88										83		87		90	92									
*16:91										83		86	87	90	92									
*16:103										83		87		90										
*16:109								80				86		89	90	92								
*16:113						79				83	85		88						95					
*16:121			75			78				83	85	87		90					95					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*16:139										82	83	85					90						95	
*16:140											83			87										
*16:150						78							86	88		90	92							
*16:169								80				85			89	90							95	
*16:170										83	85			88		90							95	
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:61, 17:63-17:69		74									83			87			91	93						
*17:01:04		74									83	85	87				91	93						
*17:07, 17:17, 17:23, 17:33											83			87			91	93						
*17:20		74									83						91	93	95					
*17:22		74															91	93						
*17:62		74												87			91	93						
*17:70		74						80			83			87			91							
*17:71		74									83			87			91							
*18:01:01:01-18:02:05, 18:04-18:08, 18:10-18:13, 18:15-18:16, 18:18-18:19											83							92	93				95	
*18:03		74									83							92	93	94			95	
*18:09											83			87				92	93					
*18:14											83							92					95	
*18:17											83							92	93					
A*01:110, A*02:636, A*03:210, A*11:135, A*24:137, A*32:96, A*68:160, B*07:234, B*15:125:02, B*15:315, B*18:126, B*39:138, B*44:85:01-44:85:02, B*50:56, B*51:226, B*55:30, B*58:74, B*58:107			75																					
A*24:64, A*24:96, A*24:106, B*08:16, B*08:123, B*08:145, B*44:347, B*44:527, B*50:76																							95	
A*24:174, B*57:129											83												95	
B*13:13:02, B*40:01:49-40:01:50																							94	
B*15:96, B*35:537, B*44:274						76																		
B*15:277, B*15:430, B*35:222, B*35:525, B*42:22, B*47:01:06, B*54:38														87										
B*15:539, B*39:76, B*54:18	73																							

Negative Control

Negative Control

