

Well No.	HLA-A allele <sup>1,2</sup>	ser <sup>3</sup>	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
	*01:01:01:01:01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:55, 01:01:57-01:01:64, 01:01:67-01:01:78, 01:01:80-01:01:142, 01:02:01:01:01:04:01:02N, 01:06, 01:09:01:01:01:11N, 01:16N-01:18N, 01:20, 01:22N-01:25, 01:27N, 01:29, 01:31N-01:33, 01:35-01:42, 01:44-01:50, 01:52:01N-01:78, 01:80-01:94, 01:96-01:98, 01:100-01:120, 01:122-01:125, 01:127-01:135, 01:137-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:175, 01:177, 01:179N-01:191, 01:193, 01:195-01:199, 01:201-01:204, 01:206-01:207, 01:209-01:213, 01:215-01:227, 01:230-01:235, 01:237-01:243, 01:245-01:266:02, 01:268-01:270, 01:272-01:278, 01:280-01:288, 01:290N-01:297, 01:300, 01:302-01:313, 01:315-01:343, 01:345-01:359, 01:362-01:375, 01:377-01:379N, 01:382-01:390, 01:392-01:395, 01:397-01:402, 01:404-01:417N, 01:419-01:426	A1, Null, -				1																						
	*01:01:23, 01:01:48, 01:01:65-01:01:66, 01:01:79, 01:08, 01:14, 01:30, 01:79, 01:99, 01:145, 01:147Q, 01:149, 01:159Q, 01:167, 01:178N, 01:205, 01:208Q, 01:214, 01:228Q, 01:236, 01:267, 01:271, 01:279, 01:314, 01:376, 01:391, 01:396Q, 01:403, 01:418N	A1, Null, -				1																						
	*01:01:56, 01:361N	Null, -	1			4								12														
	*01:01:143, 01:15N	Null, -	1			w																						
	*01:07, 01:95, 01:289	A1, -	1			4													17									
	*01:12, 01:19, 01:200, 01:244	A1, -	1			3																						
	*01:13, 01:176	A1, -	1			4							11		13		15											
	*01:21, 01:126, 01:360	A1, -	1			3	4																					
	*01:26, 01:136, 01:192, 01:192, 01:380:01-01:380:01:02, 11:94, 11:112, 11:211, 11:290, 11:326 <sup>4</sup>	A1, -	1			4							11															
	*01:28, 01:229, 01:299, 11:282 <sup>4</sup>	-	1			4						11				15												
	*01:43, 11:271 <sup>4</sup>	-	1			4						10																
	*01:51	-	1			4				8		10																
	*01:121	-	1			4														19								
	*01:194	-	1			4									13		15											
	*01:298	-	1			4														17								
	*01:301Q	-	1			4	5		7																			
	*01:344, 01:381	-	1			4				8																		
	*01:427, 36:02 <sup>4</sup>	A36, -	1			3										14								22	23			
	*02:01:01:01:02:01:07, 02:01:09-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:16, 02:18-02:22:02, 02:24:01-02:33, 02:36-02:45, 02:47, 02:49:01-02:54, 02:57-02:61, 02:63-02:64:02, 02:66-02:69, 02:71-02:77, 02:79:01-02:80, 02:82N-02:97:02, 02:99:01-02:99:02, 02:101:01-02:102, 02:104-02:123, 02:125N-02:128, 02:130-02:134, 02:136-02:148, 02:150:01-02:168, 02:170-02:236, 02:238-02:242, 02:244-02:308, 02:310-02:406, 02:408-02:448, 02:450-02:453, 02:455-02:526, 02:528:01-02:570:02, 02:572-02:575, 02:577-02:579, 02:581Q, 02:583-02:643N, 02:645-02:650, 02:652-02:670, 02:672Q-02:681, 02:683-02:740, 02:742-02:803N, 02:805Q-02:808, 02:810-02:814, 02:816-02:828, 02:830:01-02:841, 02:843-02:906, 02:908-02:926, 02:928, 02:930-02:931, 02:933-02:946N, 02:948-02:1049, 02:1051-02:1089	A2, Low A2, A203, A210, Null, -				2																						
	*02:17:02:01-02:17:04	A2, -				2			w																			
	*02:34-02:35:01, 02:35:03, 02:56:01-02:56:02, 02:62, 02:103, 02:580, 02:907	A2, -				2										13						20		22				
	*02:35:02	A2													13							20		22				
	*02:46, 02:70, 02:1050	A2, -				2																20		22				
	*02:48, 02:129, 02:571	A2, -																				20		22				
	*02:55	A2, A28				2				8													21		23			
	*02:65, 02:407, 02:449, 02:932	A2, A19, -				2																	19					
	*02:78	-	w			4																20		22				
	*02:81, 02:124	A2				2					9																	
	*02:135	-				2										13											24	
	*02:149:01-02:149:02, 02:842, 02:947	-				2																	21					
	*02:169	-				2																						
	*02:237, 02:809, 02:927	-				2										14		16								23		
	*02:243:01-02:243:03	-				2																	18		21			
	*02:309	-				2								12		14								21		24		
	*02:454	-				2								12		14										24		
Well No.			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
*02:527, 02:582	-	2						8																	
*02:576, 02:682	-	2																					22		
*02:644	-							8	10												21				
*02:651	-	2																			20				
*02:671	-	2																					23		
*02:741	-	2						8					13								20	21	22		
*02:804	-	2					7																		
*02:815	-	2						8														21			
*02:829	-			4					9																
*02:929	-	2																18							
*03:01:01-03:01:18, 03:01:20-03:01:22, 03:01:24-03:01:37, 03:01:39-03:01:82, 03:01:84-03:01:93, 03:01:95-03:07:02, 03:09-03:11N, 03:13-03:14, 03:16-03:17:02, 03:20-03:23:02, 03:25-03:29, 03:31, 03:33-03:35, 03:37-03:40, 03:42, 03:44-01-03:49, 03:51-01-03:56:01:02, 03:58, 03:60-03:62:01:02, 03:64-03:71, 03:73-03:74, 03:76-03:81, 03:83-03:87, 03:90-03:94, 03:96-03:103:01, 03:104-03:106, 03:109-03:110, 03:112-03:123:01, 03:124-03:134, 03:136-03:141, 03:143-03:151, 03:153-03:170, 03:174-03:175, 03:179-03:185, 03:188-03:193, 03:195-03:197:02N, 03:201-03:202, 03:204, 03:206-03:210, 03:212-03:214, 03:216-03:218, 03:220-03:224, 03:226-03:230, 03:231-02-03:245, 03:247-03:248, 03:250-03:251, 03:253-03:259, 03:261-03:266N, 03:268-03:272, 03:274-03:281, 03:283N-03:321, 03:323N-03:329N, 03:331-03:342N, 03:344-01-03:346, 03:348-03:376, 03:378-03:381N, 03:383-03:394, 03:396-03:397, 03:399-03:412, 03:415-03:427, 03:430-03:458Q	A3, Null, -			3									13												
*03:01:19, 03:103:02	-			3								12	13											24	
*03:01:23, 03:01:83, 03:01:94, 03:08:01:01-03:08:01:02, 03:36N, 03:57, 03:59, 03:72, 03:89:02, 03:107-03:108, 03:111, 03:142, 03:172-03:173, 03:176, 03:178N, 03:198-03:199, 03:203, 03:205, 03:211, 03:252, 03:267, 03:330N, 03:398, 03:413	A3, Null, -			3																					
*03:01:38, 03:123:02, 03:171, 03:231:01, 03:428	-			3									13	14											
*03:12	A3			3	4																				
*03:15, 03:19, 03:30	A3			3		5							13												
*03:18, 03:135	-	1		4									13												
*03:24, 03:50	A3, -			3				8					13												
*03:32, 03:219, 03:429	-			3													17								
*03:41, 03:382	-			3																			23		
*03:43, 03:82, 03:186	-			3									13	15		17									
*03:63	-			3							11		13										23		
*03:75, 03:177, 03:377	-			3	4								13										23		
*03:88	-			3	4						11												23		
*03:89:01	-		2	3																					
*03:95, 03:343	-			3		5							13	14		16							23		
*03:152, 03:273	-			3		5											17								
*03:187	-			3									13										22		
*03:215, 03:225, 03:249, 03:322, 03:414, 25:77 <sup>4</sup>	-			3									13												
*03:246, 03:282	-			3									13						19						
*03:347	-			3		5	7						13										23		
*03:395	-	1		3									13	15											
*11:01:01-11:02:07, 11:04-11:07:01, 11:08-11:09, 11:12:01-11:16, 11:18-11:19, 11:21N-11:24:02, 11:26-11:27, 11:29-11:39, 11:41-11:52Q, 11:54-01-11:59, 11:61-11:93, 11:95-11:111, 11:114-11:115N, 11:117-11:129, 11:131-11:138, 11:141-11:154, 11:156-11:157, 11:159-11:161, 11:163-11:165, 11:167-11:174, 11:176-11:182Q, 11:184-11:190, 11:192-11:198, 11:200-11:208:02N, 11:210N, 11:212-11:221, 11:224-11:225, 11:227-11:268, 11:270, 11:272Q-11:281, 11:283-11:289, 11:291-11:305, 11:307-11:324, 11:327-11:347N, 11:349-11:365N, 11:367-11:371, 11:373-11:380, 11:382N-11:438	A11, Null, -			4								11													
*11:03:01-11:03:01:02, 11:20, 11:25:01-11:25:03, 11:60, 11:158, 11:175, 11:209, 11:348, 11:366, 11:381	A11, -			3	4							11													
*11:07:02, 11:372	-			4								11						19							
*11:10, 11:191, 11:269	A11, -			4				8				11													
*11:11, 34:19 <sup>4</sup>	-			4								11	12												
*11:17, 11:40, 11:223	-			4							10	11													
*11:113, 11:162	-			4								11			15										
*11:116, 11:140, 11:199:01-11:199:02, 11:222, 66:23 <sup>4</sup>	-			4								11	13												
*11:130	-			3	4							11		14									23		
*11:139	-			4		5	7					11													
*11:155	-			4								11											22		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24



*26:54	-						8	10	12	13						20										
*26:65	-						8	10		13									24							
*26:72, 26:184	-			4				10		13																
*26:78	-						8		11									23								
*26:92, 66:42 <sup>4</sup>	-								11	13									24							
*26:120	-		1				8	10		13																
*26:122	-						8	10		13						20										
*26:152, 26:193	-						8	10	11	13																
*26:177	-								11	13				18		21		23								
*26:198	-						8	10		13	14															
*26:204:01:01-26:204:01:02	-			4			8	10																		
*26:207	-							10		13				18												
*26:217	-						8	10		12	13															
*26:229	-						8		11	12	13					20			23							
*29:01:01:01-29:06, 29:08N-29:12, 29:15-29:27, 29:29-29:38, 29:40-29:47, 29:50-29:66, 29:68-29:79, 29:81-29:97, 29:99-29:104, 29:106-29:154, 29:157-29:163, 29:165-29:171	A29, Null, -										14															
*29:07, 29:49	A29, -				5	6					14															
*29:13, 29:39, 29:98	-										14			17												
*29:14, 29:172	-										14		16													
*29:28	-						8				14															
*29:48	-										14				18											
*29:67, 29:156	-			4							14															
*29:80, 29:155	-										14						21									
*29:105, 29:164	-								10		14				18											
*30:01:01:01-30:04:04, 30:06-30:07, 30:09:01:01-30:12:02, 30:14L-30:15, 30:17-30:20, 30:23-30:43, 30:45, 30:47-30:54, 30:56-30:88, 30:90-30:124, 30:126-30:139, 30:141-30:149, 30:151-30:162, 30:164-30:176, 30:178N-30:185, 30:187-30:206, B*07:260, C*12:328 <sup>4</sup>	A30, Null, -											15														
*30:08:01-30:08:02, 30:150	A30, -				4									15												
*30:13, 30:16, 30:22, 30:44, 30:46	A30, -		2											15												
*30:55, 30:186	-										13		15													
*30:89, 30:177	-												15		17											
*30:125, 30:140	-								11				15													
*30:163	-					7							15													
*31:01:02:01-31:02:02, 31:05-31:07, 31:09-31:20, 31:22-31:23, 31:26-31:28, 31:30-31:34, 31:36-31:40, 31:42-31:61, 31:63-31:84, 31:86-31:108, 31:110-31:122, 31:124-31:175, 31:177-31:178, 31:180-31:199, 31:201N-31:219	A31, Null, -													16												
*31:03	A31			w						12	14	16														
*31:04:01:01-31:04:02	A31, -			w						12	14	16							24							
*31:08, 33:53 <sup>4</sup>	A31, A24, -				5									16		18										
*31:21	-													16	17											
*31:24-31:25	-										13			16												
*31:35	-											15	16	17												
*31:41	-													16				21								
*31:62	-													16				22								
*31:85	-					6								16												
*31:109, 33:125, 33:131 <sup>4</sup>	-													16		18										
*31:123	-									12	14	16							24							
*31:176, 31:179	-				w									16												
*31:200	-			1										16												
*32:01:01:01-32:02, 32:07-32:14, 32:16-32:25, 32:27N-32:37, 32:39-32:41, 32:43-32:59, 32:61-32:63, 32:65-32:70, 32:73-32:77, 32:81-32:92N, 32:94-32:106:01:02, 32:108-32:127, 32:129-32:141, 32:143-32:151, 32:153-32:155, 32:157-32:169, 32:171-32:173	A32, Null, -							9								17										
*32:03:01:01-32:03:01:02, 32:38, 32:60, 32:71, 32:78, 32:80, 32:107, B*15:642 <sup>4</sup>	A32, -															17										
*32:04, 32:142	A32, A3, -			3				9								17										
*32:05, 32:79	A32, -				5									16	17	18										
*32:06	A32							9								17		21								
*32:15	A32						8	9								17	18									
*32:26:01-32:26:02, 32:152	-							9			13					17										
*32:42	-							9				14				17										
*32:64	-			4				9								17										
*32:72	-					6		9								17										
*32:93	-							9								17	18									
*32:128	-							9						16		17										
*32:156	-		1					9								17										
*32:170	-					5	6	9								17										
Well No.			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24



*74:13	-										13		16			19									
*74:19	-			4												19									
*74:23	-		3													19									
*74:30	-								11							19									
*74:35	-															19	21								
*74:37	-										13					19									
*74:44	-												16			19									
*80:01:01:01-80:01:02, 80:03, 80:05-80:09N	A80, Null, -					6													23						
*80:02	-					6			11																
*80:04	-					6								17					23						
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*49:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:01:38, B*57:60	B53, -								9																
B*07:439, B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569, C*07:990, C*07:994	-																	22							
B*18:27	-					6																			
B*35:01:60, B*58:01:24	-		w																						
C*03:365, C*12:131	-			4																					
C*06:187	-																		23						
Well No.	ser. <sup>3</sup>	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24