

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01:01-01:02:17, 01:02:19-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72-01:02:90, 01:05-01:08, 01:11-01:13, 01:15-01-01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:47, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01-01:89N, 01:91-01:100, 01:103-01:112, 01:114:01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:214, 01:216-01:217:01:02, 01:219-01:232, 01:234-01:243, 01:245, 01:247	73																			92				
*01:02:18	73											84								92				
*01:02:49, 01:136, 01:145:01N, 01:146:01-01:146:02, 01:233	73							80												92				
*01:02:68, 01:02:71, 01:168-01:169:01	73														87					92				
*01:03:01-01:03:03, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:176	73																			92	93			
*01:04	73											84		86						92				
*01:09:01-01:09:02, 01:244	73													86						92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118, 01:186	73											83								92				
*01:21	73											84	85											
*01:43	73	74										82								92				
*01:48	73		75																	92				
*01:60	73																			91	92			
*01:67, 01:122	73			76												88				92				
*01:90, 01:101-01:102, 01:113	73							80												89	92			
*01:116	73																			89	92			
*01:125, 01:218, 14:58 ³	73																			88	92			
*01:131	73			76																92	93			
*01:160	73													86						90	92			
*01:169:02, 03:94 ³	73		75												87					92	94			
*01:177	73															88				92	94			
*01:215	73												85							92	94		95	
*01:246	73																			92	94			
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:30-02:02:37, 02:02:39-02:02:43, 02:02:45-02:02:52, 02:02:54-02:04, 02:07-02:11, 02:15-02:16:02:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:110, 02:112-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186, 02:188, 02:190, 02:192N-02:196, 02:198-02:217		74	75			78						83								87				
*02:02:13, 02:02:29		74	75									83								87			94	
*02:02:26, 02:21, 02:26:01-02:26:03, 02:39, 02:45, 02:189, 02:191		74	75									83								87				
*02:02:38, 02:02:53		74	75			78						83								87			94	
*02:02:44		74	75			78						83	84							87				
*02:05:01-02:05:03, 02:17		74	75			78						83			86	87								
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80				83				87				89				
*02:12, 02:55:01-02:55:02		74	75			78						83		85		87								
*02:13, 02:18, 02:33, 02:75		74	75			78						83				87				90				
*02:14:01-02:14:02, 02:27:01-02:27:02, 02:65, 02:131, 02:164, 02:187		74	75			78						83												
*02:22, 02:62, 02:118		74	75			78						83				87						93		
*02:40:01-02:40:02		74	75			78						83				87				91				
*02:47, 02:111		74	75			78		80				83				87								
*02:49		74	75			78						83		85		87				90				
*02:58			75			78						83				87								
*02:68		74	75			78						83				87				89				
*02:87		74										82												
*02:94		74	75			78	79					83				87								
*02:101		74	75			w						83				87								
*02:107		74										83												
*02:115		74	75			78						83		85						90				
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*02:134, 02:152, 02:167		74									83				87									
*02:180		74	75			78					83				87			90	91					
*02:197			75			78					83				87						93			
*03:02:01-03:02:09, 03:02:11-03:02:28, 03:05:01:01-03:05:01:03, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505, 03:532, 03:544, 03:547, 03:549, 03:555, 03:561, 03:575Q, 03:584, 03:589, 03:596, 03:612-03:615, 03:618N-03:619N, 03:623, 03:625N-03:626			75	76																		94		
*03:02:10, 03:264, 03:610			75	76																				
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:63, 03:03:65-03:03:66, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495, 03:503, 03:506-03:507, 03:509N, 03:514-03:517, 03:519-03:520, 03:524, 03:527-03:529, 03:533, 03:537, 03:541-03:542, 03:545, 03:550-03:551, 03:556, 03:559-03:560N, 03:563, 03:567-03:568, 03:573-03:574, 03:578-03:579, 03:582-03:583, 03:585-03:587, 03:592, 03:594-03:595, 03:597-03:598, 03:606, 03:609, 03:611, 03:616-03:617, 03:622, 03:624					76	77														93	94			
*03:03:21, 03:03:51, 03:18:01, 03:39, 03:64:02, 03:144, 03:182, 03:195, 03:344:01:01-03:344:01:02, 03:504				76																		93		
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451, 03:534-03:535:01:02, 03:540, 03:590, 03:607						77																93	94	
*03:03:64				76	77														90			93		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

B*58:02:02									80													92				
B*73:01:01:01-73:03																							w			
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96		