

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
HLA-C allele ^{1,2}																								
*01:02:01-01:02:17, 01:02:19-01:02:48, 01:02:50-01:02:67, 01:02:69-01:02:70, 01:02:72-01:02:94, 01:05-01:08, 01:11-01:13, 01:15-01-01:20, 01:22-01:23, 01:25-01:33, 01:35-01:42, 01:44-01:47, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:79-01-01:89N, 01:91-01:100, 01:103-01:112, 01:114-01-01:115, 01:117N, 01:119-01:121Q, 01:123-01:124, 01:126-01:130, 01:132-01:135, 01:137N-01:144, 01:145:02N, 01:147-01-01:159, 01:161-01:167, 01:170-01:175, 01:178-01:185Q, 01:187-01:214, 01:216-01:217:01:02, 01:219-01:232, 01:234-01:243, 01:245, 01:247-01:279N	73																			92				
*01:02:18	73											84								92				
*01:02:49, 01:136, 01:145:01N, 01:146:01-01:146:02, 01:233	73							80												92				
*01:02:68, 01:02:71, 01:168-01:169:01	73														87					92				
*01:03:01-01:03:03, 01:24, 01:34, 01:49:01-01:50, 01:55, 01:78, 01:176	73																			92	93			
*01:04	73											84		86						92				
*01:09:01-01:09:02, 01:244	73													86						92				
*01:10	73	74																		92				
*01:14, 01:59, 01:118, 01:186	73											83								92				
*01:21	73											84	85											
*01:43	73	74																		92				
*01:48	73		75																	92				
*01:60	73																		91	92				
*01:67, 01:122	73			76												88				92				
*01:90, 01:101-01:102, 01:113	73							80									89			92				
*01:116	73																89			92				
*01:125, 01:218, 14:58 ³	73															88				92				
*01:131	73			76																92	93			
*01:160	73													86				90		92				
*01:169:02, 03:94:01-03:94:02 ³	73		75												87					92		94		
*01:177	73															88				92		94		
*01:215	73													85						92			95	
*01:246	73																			92		94		
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:02:30-02:02:37, 02:02:39-02:02:43, 02:02:45-02:02:52, 02:02:54-02:02:67, 02:02:69-02:04, 02:07-02:11, 02:15-02:16:02:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38:02N, 02:42-02:43:01, 02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61:02, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95-02:100, 02:102-02:106, 02:108-02:110, 02:112-02:114, 02:116-02:117, 02:119-02:130, 02:132-02:133, 02:135N-02:151, 02:153-02:163, 02:165N-02:166, 02:168-02:179, 02:181-02:186, 02:188, 02:190, 02:192N-02:196, 02:198-02:223N, 02:225-02:240Q, 02:242-02:250		74	75			78					83			87										
*02:02:13, 02:02:29		74	75								83				87							94		
*02:02:26, 02:21, 02:26:01-02:26:03, 02:39, 02:45, 02:189, 02:191, 02:224		74	75								83				87									
*02:02:38, 02:02:53		74	75			78					83				87							94		
*02:02:44		74	75			78					83	84			87									
*02:02:68, 02:12, 02:55:01-02:55:02		74	75			78					83		85		87									
*02:05:01-02:05:03, 02:17		74	75			78					83			86	87									
*02:06:01-02:06:02, 02:23, 02:36:01-02:36:02		74	75			78		80			83				87		89							
*02:13, 02:18, 02:33, 02:75		74	75			78					83				87			90						
*02:14:01-02:14:02, 02:27:01-02:27:02, 02:65, 02:131, 02:164, 02:187		74	75			78					83													
*02:22, 02:62, 02:118		74	75			78					83				87						93			
*02:40:01-02:40:02		74	75			78					83				87				91					
*02:43:02, 02:58			75			78					83				87									
*02:47, 02:111		74	75			78		80			83				87									
*02:49		74	75			78					83		85		87			90						
*02:68		74	75			78					83				87		89							
*02:87		74									82													
*02:94		74	75			78	79				83				87									

Negative Control

Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*02:101		74	75			w					83				87									
*02:107		74									83													
*02:115		74	75			78					83		85					90						
*02:134, 02:152, 02:167		74									83				87									
*02:180		74	75			78					83				87			90	91					
*02:197			75			78					83				87						93			
*02:241		74					79				83				87									
*03:02:01-03:02:09, 03:02:11-03:02:34, 03:05:01-01-03:05:01:03, 03:17:01-03:17:02, 03:25, 03:27, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:132:02, 03:135, 03:139:01-03:139:02, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:431, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505, 03:532, 03:544, 03:547, 03:549, 03:555, 03:561, 03:575Q, 03:584, 03:589, 03:596, 03:612-03:615, 03:618N-03:619N, 03:623, 03:625N-03:626, 03:637:01-03:638, 03:641, 03:643N, 03:647, 03:666Q, 03:668, 03:673, 03:676, 03:680N-03:681, 03:684, 03:686, 03:696-03:697N			75	76																	94			
*03:02:10, 03:264, 03:610			75	76																				
*03:03:01-01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:63, 03:03:65-03:03:69, 03:03:71, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30-03:31, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96:01-03:96:02, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495:01-03:495:02, 03:503, 03:506-03:507, 03:509N, 03:514-03:517, 03:519-03:520, 03:524, 03:527-03:529, 03:533, 03:537, 03:541-03:542, 03:545, 03:550-03:551, 03:556, 03:559-03:560N, 03:563, 03:567-03:568, 03:573-03:574, 03:578-03:579, 03:582-03:583, 03:585-03:587, 03:592, 03:594-03:595, 03:597-03:598, 03:606, 03:609, 03:611, 03:616-03:617, 03:622, 03:624, 03:632, 03:635N-03:636N, 03:639-03:640, 03:651, 03:659N, 03:661, 03:667, 03:670, 03:679, 03:685, 03:688-03:691, 03:693-03:694				76	77																93	94		
*03:03:21, 03:03:51, 03:18:01, 03:39, 03:64:02, 03:144, 03:182, 03:195, 03:344:01-01-03:344:01:02, 03:504, 03:628N				76																	93			
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451, 03:534-03:535:01:02, 03:540, 03:590, 03:607					77																93	94		
*03:03:64				76	77													90				93		
*03:03:70, 03:324, 03:479, 03:601, 03:634				76	77																	93		
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96

Negative Control

[illegible]

[illegible]

*05:01:01:01-05:01:50, 05:01:52-05:01:55, 05:01:57, 05:01:59-05:01:82, 05:03, 05:05:01-05:15, 05:18:01-05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44:01-05:45, 05:47-05:54, 05:56-05:57, 05:59-05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106:01-05:106:03, 05:108-05:111, 05:114, 05:116-05:122, 05:124-05:128N, 05:130-05:131, 05:133, 05:136-05:142, 05:144-05:145, 05:148-05:150, 05:152-05:155, 05:157-05:162, 05:164-05:190, 05:192-05:199, 05:201, 05:203-05:205, 05:208N-05:214, 05:216:01:01-05:237, 05:239N-05:240, 05:242-05:244N, 05:246-05:251N, 05:253N-05:271N, 05:273-05:280, 05:282Q-05:287, 05:290-05:294, 05:296N-05:302					79				83				87								93											
*05:01:51, 05:01:58, 05:04:01-05:04:02, 05:103:01-05:103:02, 05:123, 05:135, 05:147, 05:191, 05:207, 05:272, 05:295					79								83																			
*05:01:56					79				80				83				87								93							
*05:16								79					83				87						92	93								
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79, 05:115					79								83				87				90				93							
*05:20, 05:43, 05:113N, 05:156, 05:202Q, 05:288					79								83												93							
*05:22:01-05:22:02, 05:252, 05:289					79				80				83				87															
*05:23, 05:62, 05:146, 05:151, 05:163								79					83				87			91			93									
*05:26					78				79				83				87								93							
*05:36								79					83				87		89				93									
*05:42					79								83				85				87				90				93			
*05:46:01:01-05:46:01:02					79								83				85				87								93			
*05:58:01-05:58:04, 05:281	75								79								83				87								93			
*05:78:01-05:78:02					78				79				83																93			
*05:85	73							79					83				87					92	93									
*05:105		74						79					83				87						93									
*05:107	73							79					83				87					92										
*05:112								79					83			85	87						93	95								
*05:129					79								83				87								93							
*05:132, 05:215								79	80				83				87		89													
*05:134, 05:143					79								83				87				90				91				93			
*05:200	73	74						79					83				87						93									
*05:206, 17:05 ³	74								79								83				87				91				93			
*05:238								79					83				87		89													
*05:241	73							79	80				83				87					92										
*06:02:01:01-06:02:01:101Q, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:02:44-06:02:64, 06:02:66-06:02:71, 06:02:73-06:02:101, 06:02:103-06:02:109, 06:02:111-06:02:118, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30:01:02, 06:33-06:34:02, 06:36-06:39, 06:41-06:44:01, 06:45-06:51, 06:53:01-06:66, 06:68-06:69, 06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102:01-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:152N, 06:154-06:179, 06:181-06:196, 06:198-0																																

*06:132:02	75								83				86 87				94							
*06:204	73							80			83						92		95					
*06:335									83				85 86				95							
*06:359N							79	80			83		86		88				95					
*06:379									80			83		86			90		95					
*07:01:01:01-07:01:02:16, 07:01:04-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:01:94-07:02:07, 07:02:10:01-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:06:06, 07:08, 07:11-07:19, 07:21-07:25, 07:27:01-07:27:02, 07:29:01:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126, 07:128-07:132, 07:133:02-07:138, 07:140:01-07:147, 07:149-07:151, 07:153:01-07:153:02, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:215, 07:217-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:389, 07:391-07:401, 07:404, 07:405:02, 07:407-07:414, 07:417-07:448, 07:450:01-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:530:02, 07:532-07:535, 07:537-07:561, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:713Q, 07:715:01-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844, 07:846:01-07:849N, 07:851-07:856N, 07:858-07:869, 07:871-07:890, 07:893, 07:895-07:931, 07:933-07:935, 07:937-07:939, 07:940:02-07:969, 07:971-07:974Q, 07:976-07:986:02, 07:988, 07:990-07:992, 07:994-07:1009, 07:1011-07:1027, 07:1029-07:1044N, 07:1046-07:1055, 07:1057-07:1058, 07:1060, 07:1062N-07:1108, 07:1110-07:1118, 07:1120N, 07:1123, 07:1125-07:1134, 07:1137-07:1143, 07:1145N-07:1163N									81									95						
*07:01:41, 07:02:41, 07:223, 07:415-07:416, 07:1045										81									w					
*07:01:74						78				81			85			89	90		95					
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:26:01:01-07:26:03, 07:33N, 07:92, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:314:01, 07:314:03, 07:317, 07:322, 07:351, 07:405:01, 07:449, 07:493, 07:517, 07:714, 07:724, 07:821N, 07:845, 07:870, 07:891, 07:894, 07:936N, 07:970N, 07:989, 07:1056, 07:1061										81														
*07:02:09						78				81			85			90			95					
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598, 07:1124										81	83					92			95					
*07:10, 07:28, 07:41, 07:43:01-07:43:02, 07:184, 07:196, 07:367, 07:536, 07:568, 07:857, 07:892, 07:987, 07:1028:01-07:1028:02, 07:1109, 07:1135										81								93	95					
*07:20:01-07:20:02, 07:127:01-07:127:02, 07:263:01-07:263:02, 07:390, 07:562	75									81								95						
*07:31:01-07:31:02, 07:514:01-07:514:02										81			86			90			95					
*07:39-07:40, 07:335, 07:975:01-07:975:02, 07:993:01-07:993:02, 07:1059, 07:1136										81						92								
*07:53, 07:216, 07:932										81	83								95					
*07:64, 07:402:01-07:402:02										81				88				94						
*07:96:01-07:96:02, 07:314:02, 07:578, 07:756			75							81									94					
*07:101, 07:148, 07:161, 07:1122		74								81	82								95					
*07:123:01-07:123:02, 07:173, 07:294, 07:626, 07:725, 07:775, 07:850										81						89			95					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*07:125:01, 07:356, 07:531						78			81	83			85	86										95
*07:125:02						78			81	83			85	86			90							95
*07:133:01, 07:242, 07:330:01, 07:748				76					81															95
*07:154									81								90							95
*07:177									81				86				90		92					
*07:210, 07:238, 07:1119									81	83									92					
*07:293, 07:519, 07:1144									81						88									95
*07:316, 07:338, 07:579	73								81															95
*07:563									81										92					95
*07:583		74				78			81	82			85				90							
*07:723						78			81				85											
*07:940:01									81		84													95
*07:1010						78			81	83			86											95
*07:1121	73								81															

Negative Control

[illegible]

*12:309								80					85			89		93		2				
*12:333												84	85	86			90	92						
*12:349							78		81			84	85	86										
*12:353			75								83	84	85	86	87									
*12:362											83	84	85											
*12:420			75									84												
*14:02:01:01-14:02:07, 14:02:09-14:02:20, 14:02:22-14:02:42, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:21N, 14:23-14:24:02, 14:27, 14:29-14:34, 14:36-14:37, 14:39-14:40, 14:42-14:44, 14:46:01-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:60, 14:62, 14:65-14:69, 14:71-14:76, 14:78, 14:80-14:85, 14:88-14:89, 14:91, 14:93N-14:104, 14:106-14:107, 14:109, 14:111, 14:114-14:115, 14:117N, 14:119-14:120, 14:122-14:125, 14:127-14:139, 14:142-14:146, 14:148N-14:149, 14:151-14:155, 14:157-14:158Q, 14:160-14:162, 14:164-14:166																88		92						
*14:02:08, 14:02:21												84				88		92						
*14:03:01:01-14:03:04, 14:22, 14:35N, 14:38, 14:41, 14:61, 14:70, 14:79, 14:86, 14:112-14:113, 14:118, 14:126, 14:141N, 14:156N, 14:159			75													88		92						
*14:03:05			75									84				88		92						
*14:04, 14:64, 14:140												83				88		92						
*14:06, 14:08-14:09, 14:28:01-14:28:02, 14:63, 14:105Q, 14:150, 14:163, A*29:86, A*30:62, A*31:187, A*33:63, B*07:499, B*08:249, B*08:306, B*18:64, B*35:183, B*35:252, B*35:502, B*39:114, B*39:140, B*40:138, B*41:59 ³																88								
*14:10																88	90	92						
*14:12, 14:49, 14:108											83							92						
*14:15, 14:87														86		88	90							
*14:16														86		88								
*14:17																88	90	91	92					
*14:25			74													88		92	94					
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*14:26								80								88			92					
*14:45		73														87	88		92					
*14:53			75													88								
*14:54			75													88					93			
*14:77											83					88								
*14:90			75								83		85			88							95	
*14:92																88	89							
*14:110, 14:147, B*07:239, B*14:03, B*14:112 ³																			92					
*14:116			75													88	90		92					
*14:121									81							88			92					
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:51, 15:02:53-15:02:54, 15:02:56-15:02:66, 15:05:01:01-15:05:07, 15:05:09-15:05:10, 15:05:13-15:05:17, 15:05:19-15:06:03, 15:08:01-15:09:02, 15:10:02-15:10:03, 15:13:01:01-15:13:03, 15:15, 15:17-15:19, 15:26-15:29, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44:01-15:54, 15:56-15:57, 15:59-15:61, 15:64, 15:67-15:70, 15:73-15:74, 15:78:01-15:84Q, 15:86-15:101, 15:103, 15:105Q-15:106, 15:108-15:110, 15:112-15:115N, 15:117, 15:119-15:124, 15:126-15:137, 15:139-15:143, 15:145N-15:148, 15:150-15:157, 15:159-15:161, 15:163-15:168, 15:170-15:177N, 15:179-15:180, 15:182-15:189N, 15:191-15:200, 15:202-15:207, 15:209-15:210, 15:212-15:217, 15:219-15:221, 15:223-15:224, 15:226-15:235Q, 15:237-15:243Q, 15:245, 15:247N, 15:249-15:253Q, 15:255-15:259, 15:261-15:275, 15:277-15:283, 15:285-15:289			75					80			83				87		89				93			
*15:02:10, 15:02:17			75					80			83				87		89				93	94		
*15:02:52, 15:02:55, 15:05:18, 15:22, 15:37, 15:55, 15:58, 15:72, 15:125, 15:149, 15:225, 15:276			75					80			83				87		89							
*15:03, 15:16:01-15:16:02			75					80			83		85		87		89				93			
*15:04:01:01-15:04:06, 15:24, 15:30, 15:40, 15:66, 15:76-15:77, 15:211, 15:222, 15:236, 15:246, 15:254, 15:260			75								83				87		89				93			
*15:05:11-15:05:12, 15:21, 15:71, 15:116, 15:244, 15:248			75					80							87		89				93			
*15:07:01:01-15:07:01:02			75					80								89				93				
*15:10:01, 15:10:04, 15:208								80			83				87		89				93			
*15:11, 15:118			75			78		80			83				87		89				93			
*15:12					77			80			83				87		89	90			93			
*15:23:01-15:23:02, 15:85, 15:138, 15:158, 15:162, 15:181			75					80			83					89				93				
*15:25			75										85				90	91			93			
*15:31, 15:62, 15:169			75					80			83				87		89	90			93			
*15:36			75			78		80			83					89				93				
*15:42, 15:111			74	75				80			83				87		89				93			
*15:43			75					80								89				93	94			
*15:63								80			83					89				93				
*15:65			75								83				87		89	90	91					
*15:75			75					80			83		85		87		89	90			93		95	
*15:102	73		75								83				87		89			92				
*15:104	73		75					80			83				87		89				93			

Negative Control

*15:107								79	80			83			87			89				93		
*15:144							77		80									89	90			93		
*15:178			75					79	80			83			87			89				93		
*15:190, 15:284			75									83			87			89						
*15:201			74	75					80			83			87			89		91		93		
*15:218				75	76				80			83			87			89				93		
*16:01:01-16:01:15, 16:01:18-16:01:50, 16:06-16:08, 16:11, 16:13-16:14, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49:01-16:49:02, 16:51-16:52, 16:54, 16:56-16:59, 16:62, 16:65, 16:71-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:85-16:86, 16:92-16:98:02, 16:100, 16:105-16:106, 16:111-16:112, 16:114, 16:116-16:119, 16:122, 16:125-16:131, 16:134-16:135, 16:137-16:138, 16:141-16:142, 16:146-16:148, 16:151-16:152, 16:154, 16:157-16:162, 16:164-16:165, 16:168, 16:171-16:175, 16:177-16:178, 16:180, 16:183, 16:185-16:187, 16:189-16:193, 16:196-16:198, 16:200-16:201, 16:204-16:211, 16:213-16:214, 16:216-16:217, 16:219N, 16:221, 16:223												83		85						90			95	
*16:01:16									80			83		85				90					95	
*16:01:17, 16:15:02												83	84	85				90					95	
Well No.	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
*16:02:01:01-16:02:12, 16:02:14-16:02:24, 16:09, 16:12, 16:18-16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N, 16:84, 16:89N-16:90, 16:99, 16:101-16:102, 16:104, 16:107-16:108, 16:115, 16:120, 16:123N, 16:132N-16:133, 16:136, 16:143-16:145, 16:153, 16:155-16:156, 16:163, 16:166-16:167, 16:176, 16:179, 16:181, 16:184, 16:194-16:195N, 16:202, 16:215, 16:220N, 16:222												83		85		87				90			95	
*16:02:13												83	84	85		87			90				95	
*16:04:01:01-16:04:01:03, 16:04:03-16:04:05, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78, 16:82, 16:124, 16:149, 16:182, 16:203, 16:218															86			90		92				
*16:10, 16:67, 16:188									80			83		85				90					95	
*16:20												83		85				89	90				95	
*16:29		74												86				90		92				
*16:34			75				78					83		85				90					95	
*16:35														86				90						
*16:40, 16:110							79									88								
*16:43												83		85									95	
*16:45														85				90				93	95	
*16:48												83		86	87			90						
*16:50		74										83		85				90					95	
*16:53							79									88				92				
*16:64, 16:87												83		85				89	90				95	
*16:68																		90		92				
*16:70			75									83		85		87		89	90				95	
*16:88												83			87			90		92				
*16:91												83			86	87		90		92				
*16:103												83			87			90						
*16:109								80						86				89	90	92				
*16:113							79					83		85		88							95	
*16:121			75				78					83		85		87		90					95	
*16:139											82	83		85				90					95	
*16:140												83				87								
*16:150							78							86		88		90		92				
*16:169								80						85				89	90				95	
*16:170												83		85		88		90					95	
*16:212												83		85		87		89	90				95	
*17:01:01:02-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:32, 17:34-17:61, 17:63-17:69, 17:72-17:80		74										83			87				91		93			
*17:01:04		74										83		85		87			91		93			
*17:07, 17:17, 17:23, 17:33												83			87				91		93			95
*17:20		74																	91		93			
*17:22		74																	91		93			
*17:62		74													87				91		93			
*17:70		74						80				83			87				91					
*17:71		74										83			87				91					
*17:81		74										83							91		93			
*18:01:01:01-18:02:05, 18:04-18:08, 18:10-18:13, 18:15-18:16, 18:18-18:22												83								92	93		95	
*18:03		74										83								92	93	94	95	
*18:09												83			87					92	93			
*18:14												83								92			95	
*18:17												83								92	93			
A*01:110, A*02:636, A*02:1164, A*03:210, A*11:135, A*24:137, A*32:96, A*68:160, B*07:234, B*15:125:02, B*15:315, B*18:126, B*39:138, B*40:562, B*44:85:01-44:85:02, B*44:406, B*50:56, B*51:226, B*55:30, B*58:74, B*58:107			75																					

Negative Control

[illegible]