

*68:01:01-68:01:66, 68:01:68-68:04:02, 68:06-68:08:01:02, 68:09-68:10, 68:12:01-68:12:02, 68:13:02-68:14, 68:16-68:18N, 68:21:01-68:24, 68:26-68:28, 68:30-68:40, 68:42-68:44, 68:46-68:65, 68:67-68:70, 68:72-68:83, 68:85-68:102, 68:104:01-68:104:02, 68:106-68:108, 68:110-68:116, 68:118-68:129, 68:131-68:133, 68:135, 68:137-68:156, 68:158-68:177, 68:179-68:217, 68:219-68:226, 68:228-68:229, 68:231-68:248, 68:251N-68:254, 68:256-68:257, 68:259-68:261, 68:263Q-68:270, 68:272-68:275, 68:277-68:293Q, 68:295-68:321, 68:323-68:326	A68, A28, Null, -								8											20						
*68:01:67, 68:218, 68:258	-								8																	
*68:05, 68:15, 68:20, 68:109, 68:136, 68:255	A68, -								8											20			23			
*68:08:02, 68:105:01:01-68:105:01:02, 68:157, 68:276Q	A68, -								8											20	21					
*68:11N, 68:41, 68:230	Null, -								8											20		22				
*68:13:01, 68:66:01-68:66:02, 68:134	A68, -				4				8											20						
*68:19, 68:130:02	A68, -								8					14						20						
*68:25, 68:262	A68, -								8										19	20						
*68:29	A68								8										18	20	21					
*68:45, 68:117, 68:322	-						7	8												20						
Well No.			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*68:71	-									8					13					20						
*68:84, 68:227	-									8		10								20						
*68:103:01-68:103:02	-				3					8										20						
*68:130:01	-									8				12						20						
*68:178	-					4														20						
*68:249	-						5			8										20						
*68:250	-																		18	20	21					
*68:271	-									8			11							20						
*68:294	-									8				12						20						24
*69:01:01:01-69:01:06, 69:03-69:06, 69:08-69:09	A69, -									8												21				
*69:02	-									8			11									21				
*69:07	-									8					13							21				
*74:01:01:01-74:01:02, 74:01:04-74:03:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:29, 74:31-74:32N, 74:34, 74:36, 74:38-74:42, 74:45-74:47N, 74:49-74:52	A74, Null, -																			19						
*74:01:03	-													12						19						24
*74:04, 74:21, 74:43	A74, -																			18	19					
*74:07, 74:33, 74:48	A74, -																		17	19						
*74:13	-														13		16			19						
*74:19	-					4														19						
*74:23	-				3															19						
*74:30	-												11							19						
*74:35	-																			19		21				
*74:37	-														13					19						
*74:44	-																16			19						
*80:01:01:01-80:01:02, 80:03, 80:05-80:10	A80, Null, -							6																	23	
*80:02	-							6					11													
*80:04	-							6											17						23	
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*49:67, B*51:01:110, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:01:38, B*57:60	B53, -									9																
B*07:439, B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569, C*07:990, C*07:994	-																						22			
B*07:507, C*06:187	-																							23		
B*15:01:83, B*35:01:60, B*58:01:24	-				w																					
B*18:27	-							6																		
C*12:131	-					4																				
Well No.	ser. ³		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24