

*36:02	A36	1	3											14						22	23				
*36:04	A36	1		4																22					
*43:01	A43							8		10	12	13													
*66:01:01, 66:01:03, 66:04-66:05, 66:07-66:08, 66:11, 66:13-66:14, 66:17-66:20, 66:22, 66:24	A66, A26, -							8		11		13													
*66:01:02	-							8		11															
*66:02-66:03, 66:16, 66:21	A66, A10, -							8			12										24				
*66:06	A66							8		11		13								21					
*66:09	-							8		11		13	14							21					
*66:10	-							8		11		w													
*66:12	-							8				13									24				
*66:15	-							8				13													
*66:23	-			4				8		11		13													
*68:01:01:01-68:01:31, 68:02:01:01-68:04, 68:06-68:08:01, 68:09-68:10, 68:12, 68:13:02-68:14, 68:16-68:18N, 68:21:01-68:24, 68:26-68:28, 68:30-68:44, 68:46-68:65, 68:67-68:70, 68:72-68:83, 68:85-68:104:02, 68:106-68:108, 68:110-68:116, 68:118-68:129, 68:131-68:133, 68:135, 68:137-68:142N	A68, A28, Null, -							8												20					
*68:01:32, 68:11N ⁴	Null, -							8												20	22				
*68:05, 68:15, 68:20, 68:109, 68:136	A68, -							8												20	23				
*68:08:02, 68:105 ⁴	A68, -							8												20	21				
*68:13:01, 68:66, 68:134	A68, -				4			8												20					
*68:19	A68							8				14								20					
*68:25	A68							8											19	20					
*68:29	A68							8										18	20	21					
*68:45, 68:117	-						7	8												20					
*68:71	-							8				13								20					
*68:84	-							8		10										20					
*68:130	-							8			12									20					
*69:01, 69:03	A69, -							8													21				
*69:02	-							8		11											21				
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*74:01:01-74:01:02, 74:02:01:01-74:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:27	A74, Null, -																			19					
*74:01:03	-											12								19					24
*74:04, 74:21	A74, -																			18	19				
*74:07	A74																		17	19					
*74:13	-													13		16				19					
*74:19	-				4															19					
*74:23	-				3															19					
*80:01:01:01-80:01:01:02, 80:03	A80, -						6																		23
*80:02	-						6				11														
B*07:81, B*08:52, B*18:67, B*38:41, B*51:185, B*53:05, B*53:16, B*53:33	B53, -									9															
B*18:27	-						6																		
C*07:404	-										11														
C*07:481	-																								?
C*12:131	-				4																				
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24