

*15:01:17, 15:01:31, 15:05:01-15:06, 15:33, 15:39:01-15:40, 15:65, 15:78:01-15:78:03, 15:84, 15:107, 15:136, 15:141, 15:155, 15:179:01-15:179:02, 15:185, 15:195, 15:199, 15:218Q-15:219, 15:237, 15:249, 15:288, 15:295, 15:300, 15:322, 15:332, 15:348, 15:364	B62, B62/blank, -							31			36													62							71	
*15:01:19, 15:28, 15:120, 15:149N, 15:246N, 15:285	B62, Null, -							31						36																	69	71
*15:02:01-15:02:09, 15:88, 15:112, 15:121, 15:139, 15:170, 15:213-15:214, 15:223, 15:265, 15:283, 15:291, 15:301-15:302N, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:374, 15:378	B75, Null, -													31			41							57								71
*15:03:01:01-15:03:02, 15:03:04, 15:47:01-15:47:02, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220, 15:235, 15:242, 15:266, 15:274, 15:281-15:282, 15:369, 15:379	B72, B70/blank, -										32			35	36																	71
*15:03:03	B72							31	32				35	36																		71
*15:03:05, 15:64:01-15:64:02	B71, -								32				35	36																		71
*15:04:01-15:04:03	B62, -								31				36											w							69	71

Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72					
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72					
*27:01	B27					29										40													52			w															70							
*27:02:01-27:02:02, 27:53, 27:57, 27:65N, 27:83, 27:95, 27:102, 27:126, 27:134	B27, Null, -					29										40													52			56																	70					
*27:03-27:05:15, 27:05:17-27:05:19, 27:05:21-27:05:22, 27:05:24-27:06, 27:09-27:10, 27:13, 27:17, 27:20, 27:27, 27:31, 27:35-27:37, 27:45-27:46, 27:48-27:51, 27:54-27:56, 27:58, 27:60-27:61, 27:64N, 27:66N, 27:68-27:69, 27:72-27:74, 27:78-27:80, 27:82, 27:84-27:88, 27:90:01-27:91, 27:93-27:94N, 27:96:01-27:101, 27:103-27:106, 27:108, 27:110-27:114, 27:117-27:118, 27:120-27:124, 27:128, 27:131-27:133, 27:135-27:136, 27:141-27:146	B27, Null, -					29										40													52			56																	70					
*27:05:16, 27:05:20, 27:07:01-27:07:03, 27:11, 27:15, 27:21, 27:28, 27:34, 27:38, 27:43, 27:47, 27:63, 27:70-27:71, 27:76, 27:139	B27, -																40													52			56																	70				
*27:05:23, 27:52, 27:67, 27:140	-					29																								52			56																	70				
*27:07:04	-					29											40						46							52			56																	70				
*27:08, 27:26, 27:40, 27:42, 27:44, 27:89	B2708/B27/B7, -					29											40													52			56																			71		

*40:117	-			28																	49 ?						58	60		62									70													
*40:128, 40:183, 40:198	-	25		28																	? 49 ?									61	62												71									
*40:132, 40:178	-			28																	? 49 ?		53							61	62													71								
*40:136	-	25		28			32														? ? 51									61	62														71							
*40:137	-																				? 49				55					62														71	72							
*40:148, 40:161	-																				48 ?							59		61	62														71							
*40:157, 40:201	-																				? ?									62															70							
*40:158	-	25		28									36								? ?									61	62																71					
*40:159	-												36								48 ?									61	62																71					
*40:166	-			28																	? 49 ?									61																	71					
*40:170	-																				? 49									61	62																	71				
*40:174	-			27 28																	? ?									61	62																	71				
*40:184	-																				? ?									61	62																	71				
*40:185	-																				? 49									62																		71				
*40:208	-			28																	? 49 ?									62																			71			
*40:209	-																				? ?									62																				71		
*40:231, 40:279	-	25		28			32														? 49 ? 51				53					61	62																		71			
*40:237	-			28																	? ?									61	62																		71			
*40:265N	Null																				48 49									62																				71		
*40:308	-			28																	? 49 ?									62																				71		
*40:313	-																				48 ?									61																				71		
*40:322	-																				? ?	? ?								61 62																					71	
*41:01:01-41:01:04, 41:07, 41:12, 41:14, 41:16-41:17, 41:20-41:22, 41:25, 41:28, 41:33-41:35, 41:37, 41:44	B41, -			28			32														48			51		54				61	62																			71		
*41:02:01-41:02:06, 41:11, 41:13, 41:19, 41:23, 41:27, 41:31, 41:36, 41:38-41:40, 41:42-41:43, 41:45N	B41, Null, -			28			32														48			51						61	62																			71	72	
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72			
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72			

*44:03:01:01-44:03:06, 44:03:08-B44, 44:03:15, 44:03:17-44:03:18, 44:03:21-44:03:29, 44:03:32- 44:05:04, 44:13, 44:26, 44:28:01- 44:30, 44:32, 44:35-44:36, 44:38- 44:39, 44:42, 44:45, 44:50:02- 44:51, 44:56N, 44:58N, 44:61N, 44:69-44:70, 44:76-44:79:02, 44:82, 44:85, 44:88, 44:92, 44:94, 44:96, 44:98, 44:105, 44:107-44:109, 44:114-44:115, 44:120, 44:122-44:125, 44:128:01-44:128:02, 44:133, 44:136, 44:139, 44:141, 44:143, 44:147, 44:153-44:155, 44:157, 44:159-44:161, 44:163-44:165, 44:167, 44:174-44:175, 44:178, 44:180-44:181, 44:183, 44:186, 44:188-44:189, 44:198N-44:199, 44:202, 44:204-44:205:02, 44:207, 44:209, 44:215, 44:217N, 44:222-44:223, 44:226- 44:228, 44:233-44:234, 44:237N	Null, -									28			41				53			62			70	
*44:03:07, 44:03:16, 44:03:19- 44:03:20, 44:03:30-44:03:31, 44:07, 44:81, 44:103, 44:111, 44:134, 44:194, 44:224	B44, -												41				53			62			70	
*44:06	B44												41				53	54		60	63		68	70

*46:01:01-46:01:06, 46:01:08-46:01:21, 46:04-46:05, 46:07N, 46:09-46:10, 46:14-46:16, 46:20, 46:22-46:24, 46:26-46:31, 46:34-46:39, 46:42, 46:44, 46:46, 46:48, 46:51Q-46:53, 46:55N-46:57, 46:60, 46:63-46:65	B46, Null, -					31		36						52	55				69	71
*46:01:07	-					31		36			50	52	55						69	71
*46:02	B46					31		36				w	55						69	71
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25, 46:45, 46:54, 46:58-46:59, 46:62	B46, -					31		36				52	55							71
*46:06, 46:43	-					31		36			49	52	55							71
*46:11	B22					31		36				52	55			65				71
*46:12	-	27	28			31		36				52	55						69	71
*46:17	-					31		36	37			52	55						69	71
*46:18	-					31		36				52	55			64				71
*46:32	-					31		36				52	55	w					69	71
*46:33	-					31		36		41		50	52	55						71
*46:40:01	-					31		36					52			64			69	71
*46:40:02	-					31		36					52						69	71
*46:41N, 46:49	Null, -					31		36					55						69	71
*46:47	-					31		36					52						69	71
*46:50	-					31		36		43			52	55					69	71
*46:61	-					31		36					52	55					69	71
*47:01:01:01-47:01:02, 47:06-47:10	B47, -	28										50		56			62			70
*47:02	B47	28										50		56		61	62			71
*47:03	B47	28										50		56		w	62			70
*47:04-47:05	B17, -	28	29									50					62			70

*49:02	B49		28																59	62	65				70																								
*49:03	B49																		59	62	65				70																								
*49:04:01-49:05	B49, -		28																58 59	62					70																								
*49:07	-		28				34												58 59	62					70																								
*49:18	-		28			33													58 59	62					70																								
*49:20	-		28				35												58 59	62	65				70																								
*49:21	-		28																58 59	62	65																												
*49:25	-	26																	58 59	62	65				70																								
*49:26	-								40			46	47						59		65				70																								
*49:31	-																		58 59	62	65				70																								
*50:01:01-50:01:07, 50:04-50:05, 50:08, 50:10-50:11, 50:13, 50:16, 50:18-50:19, 50:32, 50:35-50:38, 50:40-50:42	B50, -		28																59	61	62			65		71																							
*50:02	B45		28				35												59	61	62			65		71																							
*50:06	-		28																59	62			65			71																							
*50:07, 50:34, 50:43	-																		59	61	62			65		71																							
*50:09, 50:20, 50:39	-		28																59	61	62					71																							
*50:12	-		28																	61	62			65		71																							
*50:14	-		28			33													59	61	62					71																							
*50:15	-		28																59	61				65		71																							
*50:17	-		28							41										61	62			65		71																							
*50:31	-		28				34												59	61	62					71																							
*50:33	-		28					36											59	61	62					71																							
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72

*51:104, 51:118N, 51:147	Null, -									33										54									59 60											70									
*51:106	-									33										54									59 60			65		68				70											
*51:112	-									33										54			58	59	60											70													
*51:148	-				27																													68		70													
*51:157	-																			54					60			64						68		70													
*51:176	-				26					33															59 60									68		70													
*51:186	-									33					41					54					59 60									68		70													
*52:01:01:01-52:01:01:03, 52:01:03, 52:01:06-52:01:08, 52:01:11-52:01:14, 52:01:16- 52:01:28, 52:02:02, 52:04-52:05, 52:06:03-52:08, 52:10:01- 52:10:02, 52:11-52:13, 52:15, 52:17-52:18, 52:22-52:24, 52:27- 52:31:01, 52:32-52:38, 52:40- 52:53, 52:55-52:57, 52:59	B52, Null, -									33															59 60		62									70													
*52:01:02, 52:01:04, 52:01:09, 52:02:01, 52:03, 52:06:01- 52:06:02, 52:09, 52:10:03, 52:31:02	B52, -									33										54					59 60		62									70													
*52:01:05, 52:01:15, 52:14, 52:58	B52, -																								59 60		62									70													
*52:01:10	-									33									46						59 60		62									70													
*52:16	-									33	35														59 60		62									70													
*52:19	-							32		33										54					59 60		62									70													
*52:20	-									33															59		62									70													
*52:21	-									33										54					59 60		62							69		70													
*52:25	-																								59 60		62		65							70													
*52:26	-						30			33															59 60		62									70													
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
Well No.		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72

*52:39	-					33	36								59	60	62				70				
*52:54	-					33						48				59	60	62				70			
*53:01:01-53:01:11, 53:05, 53:08:01-53:08:02, 53:10, 53:16, 53:18, 53:20-53:21, 53:23-53:27, 53:29, 53:32, 53:34-53:35, 53:37, 53:40-53:41	B53, -						36				41						60					68	70		
*53:02	B53						36				41			55			60					68	70		
*53:03, 53:09, 53:11-53:13	B53/B35						36				41											68	70		
*53:04, 53:07, 53:19, 53:33	B53, -										41											68	70		
*53:06	B53/B51						36				41			54	55			60					68	70	
*53:14	-			27	28		36										60					68	70		
*53:15	-						36				41						60					68	70		
*53:17:01-53:17:02	-						36				41						60			62		70	72		
*53:22	-						36	37			41			53			60					68	70		
*53:28	-					33	36				41				55			60			62		70		
*53:30, 57:45, 57:51, 57:69 ⁴	-				28		36				41					58	60					70			
*53:31, 53:36	-										41											68	70		
*53:38	-						36				41										62		70		
*53:39	-						36				41						58	60					70		
*54:01:01-54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21-54:25, 54:27-54:29, 54:31-54:32, 54:34	B54, Null, -													52			56					63	64	71	
*54:03	B22													52			56					63	65	71	
*54:04, 54:15	B54, -													52			56					63		71	
*54:06	-						36							52			56					63		69	71
*54:09, 54:14	-						36							52			56					63		71	
*54:11	B54													48			52					63		71	
*54:12	-													52			56			58		63	64	70	

*56:39	-										52		56	w					65							71		
*56:43	-										52		56						65						69	71		
*56:44	-										35								65							71		
*56:45	-												42						52						56	71		
*57:01:01-57:01:04, 57:01:06-57:01:19, 57:01:21-57:01:22, 57:06, 57:08, 57:10, 57:15-57:16, 57:18-57:20, 57:22-57:23, 57:25-57:27, 57:29-57:30, 57:33-57:38, 57:40-57:41, 57:43-57:44, 57:48-57:50, 57:52-57:56, 57:58, 57:60-57:62, 57:64-57:65, 57:67:02-57:68, 57:72-57:75, 57:77-57:79N, 57:81	B57, Null, -																								60	66	70	
*57:01:05, 57:13, 57:21, 57:31, 57:47	B57, -																								60	66	70	
*57:01:20	-																								60	66	70	
*57:02:01-57:03:02, 57:17, 57:28N, 57:32, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:76, 57:80	B57, Null, -																									60	66	70
*57:04:01	B57																									60	66	70
*57:04:02	-																									60	70	
*57:05	B57																									60	66	70
*57:07	B57																									60	66	70
*57:09	B57																									60	66	70
*57:11	B57																									60	66	70
*57:12	B57																									60	66	71
*57:14:01-57:14:02	-																									60	66	70
*57:24	-																									60	66	70
*57:59	-																									60	66	70

