





A-B-DR 2023 6R6 A

*02:929	-	2						18				
*03:01:01:01-03:01:18, 03:01:20-03:01:22, 03:01:24-03:01:37, 03:01:39-03:01:82, 03:01:84-03:01:93, 03:01:95-03:07:02, 03:09-03:11N, 03:13-03:14, 03:16-03:17:02, 03:20-03:23:02, 03:25-03:29, 03:31, 03:33-03:35, 03:37-03:40, 03:42, 03:44:01-03:49, 03:51:01-03:56:01:02, 03:58, 03:60-03:62:01:02, 03:64-03:71, 03:73-03:74, 03:76-03:81, 03:83-03:87, 03:90-03:94, 03:96-03:103:01, 03:104-03:106, 03:109-03:110, 03:112-03:123:01, 03:124-03:134, 03:136-03:141, 03:143-03:151, 03:153-03:170, 03:174-03:175, 03:179-03:185, 03:188-03:193, 03:195-03:197:02N, 03:201-03:202, 03:204, 03:206-03:210, 03:212-03:214, 03:216-03:218, 03:220-03:224, 03:226-03:230, 03:231:02-03:245, 03:247-03:248, 03:250-03:251, 03:253-03:259, 03:261-03:266N, 03:268-03:272, 03:274-03:281, 03:283N-03:321, 03:323N-03:329N, 03:331-03:342N, 03:344-03:346, 03:348-03:376, 03:378-03:381N, 03:383-03:394, 03:396-03:397, 03:399-03:412, 03:415-03:427, 03:430-03:453	A3, Null, -											
*03:01:19, 03:103:02	-	3				12	13					24
*03:01:23, 03:01:83, 03:01:94, 03:08:01:01-03:08:01:02, 03:36N, 03:57, 03:59, 03:72, 03:89:02, 03:107-03:108, 03:111, 03:142, 03:172-03:173, 03:176, 03:178N, 03:198-03:199, 03:203, 03:205, 03:211, 03:252, 03:267, 03:330N, 03:398, 03:413	A3, Null, -	3										
*03:01:38, 03:123:02, 03:171, 03:231:01, 03:428	-	3					13	14				
*03:12	A3	3	4									
*03:15, 03:19, 03:30	A3	3		5				13				
*03:18, 03:135	-, -	1		4				13				
*03:24, 03:50	A3, -	3			8			13				
*03:32, 03:219, 03:429	-	3							17			
*03:41, 03:382	-	3										23
*03:43, 03:82, 03:186	-	3						13	15	17		
*03:63	-	3				11		13				23
*03:75, 03:177, 03:377	-							13				23
*03:88	-	3	4			11						23
*03:89:01	-	2	3									
*03:95, 03:343	-							13	14	16		23
*03:152, 03:273	-	3		5						17		
*03:187	-							13				22
*03:215, 03:225, 03:249, 03:322, 03:414, 25:77 <sup>4</sup>	-							13				
*03:246, 03:282	-		3					13		19		
*03:347	-			5	7			13				23
*03:395	-	1	3					13	15			









A-B-DR 2023 6R6 A

*32:64	-		4			9				17				
*32:72	-			6		9				17				
*32:93	-					9				17	18			
*32:128	-					9			16	17				
*32:156	-	1				9				17				
*32:170	-			5	6	9				17				
*33:01:01:01-33:01:17, 33:03:01:01-33:12, 33:14-33:18:02, 33:20-33:21, 33:23, 33:25-33:37, 33:39-33:47, 33:50, 33:52, 33:54-33:60, 33:62-33:118, 33:120-33:124, 33:126-33:130, 33:132-33:143N, 33:145-33:178, 33:180-33:183, 33:185-33:186, 33:188-33:189, 33:191-33:200, 33:202-33:233	A33, Null, -										18			
*33:13, 33:48	-					10				18				
*33:19	-			5	7					18				
*33:22	-									18	21			
*33:24	-									18		23		
*33:49, 33:144, 33:179, 33:201	-		w							18				
*33:51	-					8				18				
*33:61, 33:190	-								13	18				
*33:119	-				7					18				
*33:184	-							11		18		23		
*33:187	-			5	6					18				
*34:01:01:01-34:01:07, 34:05-34:06, 34:11-34:12, 34:14, 34:16-34:18, 34:23, 34:27, 66:40 <sup>4</sup>	A34, A10, -					8		11	12					
*34:02:01:01-34:03, 34:07, 34:10N, 34:13, 34:15, 34:20-34:21, 34:24-34:26N, 34:28-34:30N	A34, Null, -		3			8		11	12					
*34:04	A34		3			8		11	12	14				
*34:08	A34		3			8		11	12	13				
*34:09	-		3			8		12					24	
*34:22	-		3			8		11	12				22	
*36:01:01:01-36:01:02, 36:03, 36:05-36:12	A36, -	1											22	
*36:02	A36	1	3							14			22	23
*36:04	A36	1		4									22	
*43:01-43:03	A43, Null, -							10	12	13				
*66:01:01:01-66:01:07, 66:04-66:05, 66:07-66:08, 66:11, 66:13-66:14, 66:17-66:18, 66:20, 66:22, 66:24, 66:27N, 66:29-66:33, 66:35-66:38, 66:41, 66:44-66:46	A66, A26, Null, -					8		11	13					
*66:02:01:01-66:03:01:03, 66:16, 66:21, 66:25-66:26Q, 66:28N, 66:34, 66:39N, 66:43	A66, A10, Null, -					8		12						24
*66:06	A66					8		11	13				21	
*66:09	-					8		11	13	14			21	
*66:10	-					8		11	w					
*66:12	-					8			13					24
*66:19	-					8		11	12	13				





A-B-DR 2023 6R6 A

*74:35	-																	19	21						
*74:37	-											13						19							
*80:01:01:01-80:01:02, 80:03, 80:05-80:09N	A80, Null, -					6															23				
*80:02	-					6				11															
*80:04	-					6										17					23				
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*49:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:01:38, B*57:60	B53, -									9															
B*07:439, B*35:491, B*40:359, B*55:96, B*57:65, C*04:31, C*06:137, C*07:569, C*07:990, C*07:994	-																				22				
B*18:27	-					6																			
B*35:01:60, B*58:01:24	-	w																							
C*03:365, C*12:131	-			4																					
C*06:187	-																				23				
Well No.	ser. <sup>3</sup>	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24