

Well No.		88	89	90	91	92	93	94	95	96
DQB1 allele <sup>1,2</sup>	ser. <sup>3</sup>									
*02:01:01-02:02:02, 02:04-02:33	DQ2, Null, -			90						
*02:03	DQ2			90			93			
*03:01:01:01-03:01:21, 03:04, 03:09-03:10:01, 03:13-03:14:02, 03:16, 03:19, 03:21-03:22, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:60, 03:69, 03:71, 03:73, 03:75-03:77, 03:80, 03:82-03:84N, 03:92-03:94	DQ7, DQ3, Null, -				91			94		
*03:02:01-03:02:12, 03:05:01-03:05:04, 03:07-03:08, 03:11, 03:18, 03:32, 03:37, 03:45, 03:61, 03:63-03:64, 03:66N-03:68, 03:70, 03:85	DQ8, Null, -					92		94		
*03:03:02:01-03:03:10, 03:06, 03:12, 03:15, 03:20, 03:25-03:26, 03:30-03:31, 03:33-03:34, 03:38-03:41, 03:43, 03:65, 03:74, 03:79, 03:86-03:91Q, 03:95N-03:99Q	DQ9, DQ3, Null, -						93	94		
*03:17:01-03:17:02, 03:23, 03:62, 03:72, 03:78, 03:81, 03:100	DQ9, -							94		
*04:01:01-04:01:02, 04:01:04-04:02:05, 04:04-04:17	DQ4, -								95	
*04:01:03	-							94	95	
*04:03:01-04:03:02	-						93		95	
*05:01:01:01-05:47	DQ5, Null, -	88								
*06:01:01-06:03:09, 06:04:01-06:28, 06:30-06:50, 06:51:02-06:65, 06:67-06:95, 06:97-06:119	DQ6, DQ5, DQ1, Null, -		89							
*06:03:10, 06:51:01, 06:66, 06:96	-		89				93			
*06:29	DQ6		89			92				
Well No.		88	89	90	91	92	93	94	95	96

Negative Control