









A-B-DR-DQ 2022 0R5 B

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
*14:08:01-14:08:02, 14:55, 39:01:35, 39:30, 39:32-39:33, 39:43, 39:47, 39:50, 39:74, 39:82, 39:102, 39:107, 39:112, 39:128, 39:136, 39:169, 39:171, 39:174 <sup>4</sup>	B16, B39, -					26																43	44																						64	
*14:10, 14:102	-					26									37							43	44				49																			64
*14:21	-				25																		44			47																			64	
*14:37	-						27															43	44								53					58								64		
*14:42:01	-				25																		44																						64	
*14:42:02	-				25	27																	44																						64	
*14:48	-				25	27																43	44												55									64		
*14:57	-				25	27																	44				47																		64	
*14:73	-				25	27				32													44																						64	
*14:93	-				25	27																	42	43																					63	
*15:01:01:01-15:01:01:59, 15:01:03, 15:01:06-15:01:16, 15:01:18, 15:01:20-15:01:30, 15:01:32-15:01:53, 15:01:55-15:01:70, 15:01:72-15:01:74, 15:26N-15:27:04, 15:30:01:01-15:30:01:03, 15:32:01, 15:34:01-15:35:01:02, 15:38:01-15:38:02, 15:50, 15:56, 15:58, 15:60, 15:63, 15:70-15:71, 15:75, 15:79N, 15:81-15:82, 15:92, 15:94N, 15:96-15:97, 15:102, 15:104-15:105, 15:109-15:111N, 15:113, 15:117-15:118, 15:122, 15:125:01-15:125:02, 15:128-15:129, 15:135, 15:140, 15:145-15:147, 15:159-15:160:01:02, 15:163-15:164, 15:166-15:167, 15:169, 15:171-15:172, 15:174-15:175, 15:178, 15:181N-15:182N, 15:184, 15:187, 15:190N, 15:192-15:193, 15:201, 15:203, 15:205-15:206, 15:211, 15:217, 15:225, 15:227-15:228, 15:231-15:232, 15:234, 15:241, 15:244-15:245:02Q, 15:247, 15:251:01-15:251:02, 15:257-15:258N, 15:260-15:262:02N, 15:264, 15:267, 15:269, 15:272N, 15:276, 15:279-15:280, 15:286, 15:296, 15:299, 15:309, 15:315-15:317, 15:320-15:321Q, 15:326-15:327, 15:333-15:334, 15:340, 15:342-15:344, 15:346, 15:349:01-15:349:02, 15:353, 15:355, 15:359, 15:363:01-15:363:02, 15:365, 15:368, 15:370-15:373, 15:375N, 15:381, 15:385-15:387, 15:390-15:391, 15:398, 15:400N-15:401, 15:406, 15:409-15:410, 15:413, 15:415-15:417, 15:419, 15:422, 15:432, 15:440-15:441, 15:443, 15:452-15:453, 15:456-15:457, 15:459, 15:461, 15:464-15:467, 15:469-15:474, 15:476-15:478, 15:480-15:484, 15:488-15:490, 15:492-15:493, 15:496N, 15:499, 15:501, 15:512, 15:514, 15:520Q-15:522, 15:528N-15:529, 15:533-15:536, 15:538, 15:540N-15:542, 15:547, 15:549N, 15:554, 15:556, 15:558, 15:561-15:562N, 15:565, 15:568, 15:572, 15:576, 15:580-15:581, 15:583, 15:585, 15:590-15:591, 15:593, 15:595, 15:597-15:598, 15:600-15:601, 15:604N-15:606, 15:608, 15:618, 15:621, 15:624	B62, Null, -						28										39																		56				62				64			
*15:01:02, 15:278	B62, -						28			30								39																			56					62			64	
*15:01:04, 15:01:54, 15:251:03	B62, -																	39																		56					62			64		
*15:01:17, 15:01:31, 15:05:01:01-15:06, 15:33, 15:39:01:01-15:40:02, 15:48, 15:65, 15:78:01:01-15:78:04, 15:84, 15:107, 15:136, 15:138, 15:141, 15:150, 15:155, 15:179:01-15:179:02, 15:185, 15:188, 15:195, 15:199, 15:218Q-15:219, 15:237, 15:248-15:249, 15:287-15:288, 15:295, 15:300, 15:322, 15:332, 15:348, 15:364, 15:395, 15:421, 15:497, 15:519, 15:577, 15:626	B62, -						28											39																		56									64	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				
*15:01:19, 15:28, 15:120, 15:149N, 15:246N, 15:285, 15:392, 15:428	B62, Null, -							28										39																								62	64					
*15:01:71, 15:377Q, 15:430	-							28										39																									62	64				
*15:02:01:01-15:02:12, 15:88, 15:112, 15:121, 15:139, 15:170, 15:194, 15:213-15:214, 15:223, 15:265, 15:283, 15:291, 15:301-15:302N, 15:319, 15:325, 15:328, 15:330, 15:341, 15:357-15:358, 15:374, 15:378, 15:425, 15:437-15:438, 15:454N, 15:460, 15:491, 15:495, 15:537, 15:557, 15:567, 15:569, 15:574, 15:596N, 15:627	B75, Null, -							28									37	39																										64				
*15:03:01:01-15:03:02, 15:03:04, 15:03:06-15:03:09, 15:03:11, 15:03:13, 15:47:01-15:47:02, 15:49, 15:61, 15:98, 15:103, 15:123, 15:127, 15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220:01:01-15:220:03, 15:235, 15:242:01-15:242:02, 15:253, 15:266, 15:274, 15:281-15:282, 15:369, 15:379, 15:389, 15:397, 15:444, 15:448, 15:475, 15:502, 15:527, 15:539, 15:560, 15:571, 15:582, 15:610, 15:614:01:01-15:614:01:02, 15:616Q-15:617, 15:628	B70, B72, -								29			32						39																										64				
*15:03:03	B72							28	29			32						39																											64			
*15:03:05, 15:64:01-15:64:02	B71, -							29				32						39																											64			
*15:03:10, 15:03:12, 15:74	B72, -							29										39																											64			
*15:04:01:01-15:04:04, 15:137, 15:531, 15:609	B62, -							28										39																										62	64			
*15:07:01:01-15:07:03, 15:45, 15:126, 15:207, 15:324, 15:331, 15:405, 15:431, 15:450	B62, -			24	25			28										39																											62	64		
*15:08:01:01-15:08:02, 15:11:01-15:11:03, 15:11:05-15:11:07, 15:15:01:01-15:15:01:02, 15:148, 15:191, 15:209N, 15:215, 15:284, 15:305, 15:350, 15:367, 15:479, 15:513, 15:622	B62, B75, Null, -							28										39																										62	64			
*15:09:01-15:09:02, 15:435	B70, -							29	30					34				39																												64		
*15:10:01:01-15:10:01:10, 15:10:03-15:10:08, 15:18:01:01-15:18:09, 15:37, 15:51-15:52, 15:72, 15:90, 15:93, 15:99, 15:108, 15:114, 15:119, 15:133-15:134, 15:153, 15:176, 15:197:01, 15:198, 15:200, 15:221, 15:226N, 15:229, 15:238, 15:252, 15:263, 15:275:01-15:275:02, 15:290, 15:292, 15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:335, 15:337-15:338, 15:351, 15:354, 15:380N, 15:382, 15:388, 15:414, 15:426-15:427, 15:433, 15:439, 15:445, 15:447, 15:451, 15:455, 15:487N, 15:503-15:504, 15:526, 15:530, 15:548, 15:563, 15:570, 15:607, 15:611-15:612, 15:620	B70, B71, Null, -							29						34				39																											64			
*15:10:02	B71							28	29					34				39																												64		
*15:11:04	B75																	39																											62	64		
*15:12:01-15:12:03, 15:14, 15:19, 15:270, 15:298, 15:304N, 15:551	B76, Null, -							28					33					39																												62	64	
*15:13:01-15:13:03	B77, -							28						34			37	39																												63	64	
*15:16:01:01-15:16:05, 15:67, 15:254, 15:362, 15:408, 15:550, 15:602, 15:613	B63, -										30			34				39																												63	64	
*15:17:01:01-15:17:07, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:230, 15:268, 15:273, 15:356, 15:361, 15:396, 15:403, 15:411, 15:423-15:424, 15:446, 15:500, 15:516, 15:523, 15:532, 15:555, 15:575N, 15:603, 15:619	B63, Null, -							28		30				34				39																													63	64
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				

A-B-DR-DQ 2022 0R5 B

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64							
*15:20, 15:25:01:01-15:25:04, 15:204, 15:240, 15:250, 15:271, 15:366, 15:393, 15:402, 15:407, 15:458, 15:463N, 15:468, 15:543-15:544N, 15:564, 15:566, 15:584N, 15:587, 15:589, 15:599, 15:625	B62, Null, -							28									37	39																													64				
*15:21:01:01-15:21:01:02, 15:44, 15:394, 15:404, 15:420, 15:615, 15:623	B75, -							28					34				37	39													51																64				
*15:23, 15:115	B70, B5, -							29					34					39																												63					
*15:24:01-15:24:02, 15:87, 15:157, 15:442, 15:586	B15, B62, -							28					34					39																												62	63				
*15:29:01:01-15:29:01:02, 15:553	B70, -							29										39														52				56										64					
*15:31:01:01-15:31:01:03, 15:384, 15:508	B75, -							28										39														51														64					
*15:32:02	-							28										39																											62	64					
*15:36	-			w				28									37	39																												63					
*15:42	B15							28										39																													64				
*15:43, 15:256	B15, -							28										39																													62	63			
*15:46	-				25			28	29									39																													62	64			
*15:53, 15:212	-				25			29										39																														62	64		
*15:54	B72							29				32						39																														62	64		
*15:55	B15			24	25			28										39																															64		
*15:57	B62							28										39											w																		62	64			
*15:62	-							29				32						37	39																													64			
*15:66	B62							28						34				39																														62	64		
*15:68	-			24	25								32														44																						62	64	
*15:69	B72							29				32						39			41																											64			
*15:73:01:01-15:73:01:02	B62, -							28										39																														62	64		
*15:76	-							28						34				39											47																				62	64	
*15:77	B62						27	28										37	39																														62	64	
*15:80	B70							29						34				37	39																														64		
*15:83:01:01-15:83:01:02	B62, -							28	29									39																															64		
*15:85, 15:154, 15:165	B62, -							28										37	39																														62	64	
*15:86, 15:224, 15:485	-							28										39																																64	
*15:89:01-15:89:02	B77, -				w			28										37	39																														63		
*15:91, 15:131	-							29				32	33					39																															64		
*15:95	B63									30				34				39																															63		
*15:101	-							28						34																																			62	64	
*15:106	-				25			28	29									37	39																															64	
*15:116, 15:352	B62, -							28										39																																64	
*15:124	B71							29						34				39																																64	
*15:142, 15:429	-			23				28										39																																62	64
*15:143	-							28	29									39																																62	64
*15:144	-							28										37	39																															62	64
*15:152	-							28										39																																62	64
*15:161	-							29					33	34				39																																64	
*15:180	-			23																																														62	64
*15:183	-							28										39																																62	64
*15:186	-							29						34				39			41																													64	
*15:189:01-15:189:02	-						27	28						34				39																																62	64
*15:197:02, 15:329	-				24	25		29						34				39																																64	
*15:202, 15:239	-							28										38	39																															62	64
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64							

A-B-DR-DQ 2022 0R5 B

Well No.	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			
*15:222	-								30				34				39												52												63					
*15:233	-					27	28										39																								62	64				
*15:236	-						28																																		62	64				
*15:243	-							29	30				34				39																													
*15:255	-							29					34													47																				
*15:259	-							28	29			32					39																									62	64			
*15:277	-							28		30							39																									62				
*15:289	-							28								37	39											49															64			
*15:293	-								29				34				39																	53									64			
*15:297	-	22						28								37	39											49		51												64				
*15:303	-							28									39													51	w							59			62	64				
*15:308	-							28								37	38	39																						61			64			
*15:310	-							28									39																53										64			
*15:318, 15:579	-							28									39											49														62	64			
*15:336, 15:559	-				25												39																									62	64			
*15:339	-							28									39																											63		
*15:345	-							28								37	39														51				55								64			
*15:347	-				25												39																									62	64			
*15:360, 15:383	-			24	25			28									39																										64			
*15:376	-												34				38	39																										64		
*15:399	-				25			28								37	39																											64		
*15:412, 15:507	-							28									39									45																62	64			
*15:418	-												34			37	39																										63			
*15:434, 15:436, 15:573	-												34				39																											64		
*15:449, 15:494	-							28									39				41																						64			
*15:462	-							28		30				34			39																									60	63			
*15:486	-							28									39																										62	64		
*15:498	-							29					34				39																										62	64		
*15:505	-							28								37	39																											63		
*15:506	-																																											64		
*15:509	-							29					34				39										47																	64		
*15:510	-				25			28									39																										62	64		
*15:511	-							29				32					39																											64		
*15:515	-																																											64		
*15:517	-							29				32					39																											64		
*15:518	-											33					39																											64		
*15:524	-				24			28									39																											64		
*15:525	-							28	29								39																											62	64	
*15:545	-							28						35			39																											62	64	
*15:546Q	-																																											63		
*15:552	-							28									39																											64		
*15:578	-							28									39												48															64		
*15:588	-				25			29				32				36	39																												64	
*15:592	-							28				33					37	39																											64	
*15:594	-							28									37																												62	63
Well No.	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			









A-B-DR-DQ 2022 0R5 B

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
*35:14:01-35:14:03, 35:62, 35:163, 35:213, 35:343, 35:390N, 35:413, 35:512	B35, Null, -																37	38																							61	62	64							
*35:16, 35:318	B35, -			24													37	38																		56								64						
*35:17:01-35:17:04, 35:30:01-35:30:01:02, 35:113, 35:125, 35:164, 35:242, 35:346	B35, -			24													37	38																							61			64						
*35:19, 35:274	B35, -																37																			55								64						
*35:21, 35:24:01-35:24:02, 35:81, 35:96, 35:109, 35:157, 35:188, 35:190, 35:233, 35:287	B35, B78, -																37	38											49												61			64						
*35:25, 35:49:01-35:49:02, 35:198, 35:210, 35:254, 35:271	B35, -																37																											64						
*35:26	B35																37																											64						
*35:35, 35:468	B35, -						27										37	38																			54						61		64					
*35:36	B35																37	38																								61		63	64					
*35:38, 35:115, 35:169, 35:309, 35:452, 35:513N	B35, Null, -																37	38											48												61				64					
*35:43:01, 35:43:03-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:476	B35, -																37	38																									61	62		64				
*35:45, 35:71	B35												33				37	38											48												61				64					
*35:46	B35				25												37																												64					
*35:47, 35:154, 35:217, 35:450	B35, -				25												37																												64					
*35:60	B35																37	38																										61		64				
*35:63, 35:317, 40:129 <sup>f</sup>	-				25												37																												64					
*35:66	B35	22			25												37	38																										61		64				
*35:76	B35, 22																37	38										47																w		64				
*35:82, 35:85, 35:294	-																37	38																										61		64				
*35:84, 35:90, 35:150:01-35:150:02, 35:184, 35:201, 35:441	-																37	38																											61		64			
*35:87	-	22							29								37	38										46																61		64				
*35:114	-			24													37	38																											61		64			
*35:135	-																37	38																										61	62		64			
*35:137, 35:147, 35:182, 35:278, 35:540	-				25												37	38																										61		64				
*35:185	-																37	38																											61	62		64		
*35:186	-																37	38											47																	61		64		
*35:207	-				25												37																													64				
*35:218, 35:256	-			23													37	38																											61		64			
*35:229	-																37	38											47																61		64			
*35:230	-																37	38									45																		61		64			
*35:251:01	-																37	38																											61		64			
*35:260	-																37	38											48																	61		64		
*35:265	-																37	38																												61	62		64	
*35:276	-																37	38																												61		64		
*35:289, 35:487	-																37	38																												61		64		
*35:311	-																37																														64			
*35:329	-																37	38																												61		63		64
*35:331	-																37	38																													61		64	
*35:361	-																37																															64		
*35:368	-			24	25												37																														61		64	
*35:432, 35:446	-																37	38																														61		64
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						





















A-B-DR-DQ 2022 0R5 B

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64
*46:01:01:01-46:01:36, 46:04-46:05, 46:07N, 46:09, 46:14-46:16, 46:20, 46:22-46:24, 46:26-46:31, 46:34-46:39, 46:42, 46:44, 46:46, 46:48, 46:51Q-46:53, 46:55N-46:57, 46:60-46:61, 46:63-46:66, 46:68-46:70, 46:73-46:76, 46:78-46:82, 46:85, 46:87-46:88, 46:90-46:92, 46:94-46:95N	B46, Null, -							28																		47	49														62	64		
*46:02	B46							28																		w	49														62	64		
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25, 46:45, 46:54, 46:58-46:59, 46:62, 46:83, 46:93	B46, -							28																			47	49																64
*46:06, 46:43	-							28																	45	47	49																64	
*46:10	B46							28																			47	49					55								62	64		
*46:11	B22							28																			47	49									59					64		
*46:12	-			24	25			28																			47	49														62	64	
*46:17	-							28				33															47	49														62	64	
*46:18	-							28																			47	49									58					64		
*46:32:01-46:32:02	-							28																			47	49				w									62	64		
*46:33	-							28									37										47	49														64		
*46:40:01	-							28										38									47											58				62	64	
*46:40:02	-							28										38									47															62	64	
*46:41N, 46:49	Null, -							28																				49														62	64	
*46:47, 46:72, 46:86	-							28																			47															62	64	
*46:50	-							28																			47	49															62	64
*46:67	-							28									37										47	49															62	64
*46:71	-							28																			47	49															62	64
*46:77	-							28																			47																62	64
*46:84	-				25			28																			47	49															62	64
*46:89	-																							44		47	49																62	64
*47:01:01:02-47:01:06, 47:06-47:12	B47, -				25																																						63	
*47:02	B47				25																								50														64	
*47:03	B47				25																								50						w	56						63	w	
*47:04-47:05	B17, -				25	26																																					63	
*48:01:01:01-48:01:10, 48:07, 48:09, 48:15, 48:19-48:20, 48:22, 48:27, 48:30, 48:32-48:38, 48:42-48:50, 48:52	B48, -	22										32												44						51														64
*48:02:01-48:02:03, 48:25	B48, B70, -											32					37							44																			64	
*48:03:01:01-48:03:02, 48:23, 48:40	B48, -											32												44	45																		64	
*48:04:01:01-48:04:02, 48:24, 48:51, 48:53	B48, -											32																		51													64	
*48:05, 48:08, 48:10	B48, -	22										32																															64	
*48:06	-	22										32																															64	
*48:11	B48	22										32															44		48		51												64	
*48:12, 48:14	B48	22			25							32																															64	
*48:13	-											32																															64	
*48:16, 48:39	B48, -	22										32																		51													64	
*48:17	-	22										32													44	45			49														64	
*48:18	B48	22										32																		51	52											63		
*48:21, 48:26	-					26						32													44																		64	
*48:28	-	22										32																		51													64	
*48:29	-	22			25							32																		51													64	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64









Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			
*54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10, 54:13-54:19, 54:21-54:25, 54:27-54:29, 54:31-54:32, 54:34, 54:36-54:37, 54:39, 54:41-54:44	B54, Null, -																										47			50							57	58							64		
*54:03	B22																										47			50							57		59						64		
*54:06	-																										47			50							57						62		64		
*54:09	-																										47			50							57								64		
*54:11	B54																								44			47			50						57	58							64		
*54:12	-																										47			50		52					57	58						63		64	
*54:20	-																										47			50		53					57	58							64		
*54:26	-										31																47			50							57								64		
*54:30	-																											47			50						57	58							64		
*54:33	-																											47			50		w				57	59			62			64			
*54:35	-																											47			50				56	57	58								64		
*54:38	-				25																						47			50							57	58							64		
*54:40	-																										47			50							57	59			62			64			
*55:01:01:01-55:01:06, 55:01:08-55:01:28, 55:02:11, 55:03, 55:05, 55:11, 55:15, 55:17, 55:25, 55:28:01-55:29, 55:31, 55:33, 55:36, 55:38, 55:44-55:45, 55:52-55:55N, 55:58-55:60, 55:64, 55:66, 55:68, 55:73-55:76, 55:79, 55:84-55:85, 55:87, 55:91-55:95, 55:97N-55:102, 55:104-55:106, 55:108-55:109, 55:111, 55:113-55:116, 55:118N-55:119	B55, Null, -																										47			50		53							58							64	
*55:01:07, 55:02:01:01-55:02:10, 55:02:12-55:02:16, 55:07, 55:10, 55:12-55:13, 55:16, 55:19, 55:23, 55:26-55:27, 55:30, 55:32, 55:35, 55:37, 55:39, 55:41-55:43, 55:46-55:48, 55:50, 55:57, 55:61-55:63, 55:65, 55:67, 55:69-55:72, 55:77, 55:80-55:83N, 55:88-55:89N, 55:96, 55:107, 55:112, 55:117N, 55:120, 56:10, 56:66 <sup>4</sup>	B54, B55, Null, -																										47			50									58							64	
*55:04, 55:08, 55:49, 55:86, 56:01:05, 56:01:13, 56:15, 56:18-56:19N, 56:22, 56:31-56:32, 56:50 <sup>4</sup>	B55, B56, Null, -																											47			50																64
*55:09, 55:24	B22, B55										31																47			50		53						58								64	
*55:14	-																37										47			50							58									64	
*55:18	-											32																47			50					56		58								64	
*55:20	B22							29																			46	47		50							58									64	
*55:21	B55																										47			50		w					58				62				64		
*55:22	B22										31																47			50																64	
*55:34	-																										47			50					56		58									64	
*55:40	-																										47			50		53			56		58									64	
*55:51	-			24																							47			50							58									64	
*55:56	-							29																			47			50							58									64	
*55:78	-																														53						58									64	
*55:90, 55:121	-																										47			50		53					58						63			64	
*55:103	-																										47			50		52	53				58							63			64
*55:110	-				25																						47			50																	64
*55:122	-											32															47			50		53						58									64
*55:123	-																										47			50								58									64
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64			

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				
*56:01:01:01-56:01:04, 56:01:06-56:01:12, 56:01:14-56:02:03, 56:04:01:01-56:04:04, 56:08, 56:14, 56:16, 56:20:01, 56:24, 56:26-56:30, 56:34-56:35, 56:38N, 56:40-56:41, 56:47-56:49, 56:51-56:53, 56:55:01:02-56:59, 56:61-56:64, 56:69-56:73, 56:75-56:79, 56:81-56:86	B56, Null, -																										47			50										59						64		
*56:03	B22, 15																										47			50													62		64			
*56:05:01	B56									30																	47			50		53						58						64				
*56:05:02	B56									30																	47			50		53													64			
*56:06	B78									30								38									47			50		53						58		w				64				
*56:07, 56:46	B56, -																										47			50									59					63				
*56:09, 56:11-56:12	B56, B22																37										47			50															64			
*56:13, 82:04 <sup>4</sup>	B56, -								29																		47			50											59				64			
*56:17, 56:33	B56, -																																												64			
*56:20:02	-																																												64			
*56:21	-									30																		47			50														63			
*56:23	-																																												64			
*56:25, 56:37	-																											47			50		53												64			
*56:36	-																											47			50					55									64			
*56:39	-																											47			50	w													64			
*56:42	-																																												64			
*56:43	-																											47			50														62	64		
*56:44	-											32																																	64			
*56:45	-																	38										47			50												w			64		
*56:54	-																																												64			
*56:60	-										31																	47			50															64		
*56:65	-																											47			50		52													63		
*56:67	-				25																																									64		
*56:68	-																																													64		
*56:74	-											32																47																		64		
*56:80	-																	37											47			50															64	
*57:01:01:01-57:01:04, 57:01:06-57:01:22, 57:01:24-57:01:31, 57:01:33-57:01:44, 57:02:01:01-57:03:05, 57:05-57:06, 57:08, 57:10, 57:15-57:20, 57:22-57:23, 57:25, 57:27-57:30, 57:32-57:35, 57:37-57:44, 57:46, 57:48-57:50, 57:52-57:58, 57:60-57:68, 57:70, 57:72-57:81, 57:84-57:87, 57:90-57:92, 57:94-57:108, 57:110-57:117, 57:120-57:127, 57:129-57:144, 57:146, 57:148-57:156, 58:36 <sup>4</sup>	B57, Null, -																																														63	
*57:01:05, 57:13, 57:21, 57:31:01-57:31:02, 57:47, 57:82-57:83, 57:118, 57:145	B57, -																																														63	
*57:01:23, 57:01:45, 57:36, 57:89, 57:109, 58:01:01:01-58:01:02, 58:01:04-58:01:15, 58:01:17-58:01:42, 58:04-58:05, 58:10N-58:15, 58:19, 58:21-58:24, 58:28:01-58:29, 58:31N-58:35, 58:37, 58:39N-58:42, 58:45:01-58:45:02, 58:47-58:59:02, 58:62-58:63, 58:65-58:72N, 58:74-58:75, 58:77, 58:79-58:83, 58:85-58:93N, 58:95-58:99, 58:101-58:103, 58:105-58:107, 58:109-58:110, 58:112-58:113, 58:115-58:120, 58:123-58:126, 58:129-58:131, 58:133N-58:135 <sup>4</sup>	B58, Null, -																																															63
*57:01:32	-									30																																					63	
*57:04:01	B57						27																																								63	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				

Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64										
*57:04:02, 58:114 <sup>4</sup>	-						27										37																	54										63										
*57:04:03	-						27											37																											60		63							
*57:07, 57:24	B57, -																	37																												60		63						
*57:09	B57								29									37																												60		63						
*57:11	B57				25																																										60		63					
*57:12	B57																																														60		64					
*57:14:01-57:14:02	-																																															60		63				
*57:26	-																																															60		63				
*57:45, 57:51, 57:69	-				25																																											60		63				
*57:59	-																																															60		63				
*57:71, 57:88	-				25																																											60		63				
*57:93, 57:119	-																																															60		63				
*57:147	-																																																60		62	63		
*58:01:03, 58:01:16, 58:16:01-58:16:02, 58:26-58:27, 58:46, 58:61, 58:121-58:122, 58:127	B58, -																																																			63		
*58:02:01:01-58:02:02, 58:06, 58:25, 58:38, 58:43, 58:60, 58:104, 58:132	B58, -				25																																														63			
*58:07	B58				25								33																																					63				
*58:09, 58:76	B17, -																																																		63			
*58:17N, 58:108	Null, -																																																			63		
*58:18	-				24	25																																													63			
*58:20	-																																																		63			
*58:44	-												32																																					63				
*58:73	-																																																	63				
*58:78	-																																																	62	63			
*58:84	-																																																	63				
*58:94N	Null																																																	63				
*58:128N	Null																																																		63			
*59:01:01:01-59:03, 59:05-59:07, 59:10N-59:11	B59, Null, -																																																		63			
*59:04	-																																																		63			
*59:08	-																																																		63			
*59:09	-																																																		63			
*59:12	-																																																		63			
*67:01:01, 67:03-67:04, 67:06-67:07	B67, -						26																																													64		
*67:01:02-67:02:01:02, 67:05	B67, -																																																				64	
*73:01:01:01-73:03	B73, -																																																				62	64
*78:01:01:02-78:01:02, 78:02:02, 78:07-78:10	B78, -																																																				64	
*78:02:01, 78:04	B78, B35																																																		64			
*78:03	B78																																																		64			
*78:05-78:06	-																																																		64			
*81:01:01:01-81:01:05, 81:05, 81:09-81:10	B81, -				22																																														64			
*81:03-81:04N, 81:06-81:08	Null, -				22																																														64			
*82:01:01:01-82:03	B82, -																																																			64		
*83:01	-																																																			64		
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64										

