

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	
DRB1 allele ^{1,2}	ser ³																								
*01:01:01:01:01:02:13, 01:04-01:38, 01:40N-01:41, 01:43-01:60, 01:62:01N-01:88, 01:90-01:98, 01:100-01:101, 01:103-01:111, 01:113-01:126	DR1, Null, -	65																							
*01:03:01-01:03:05, 01:39N, 01:99, 01:102, 01:112	DR1, DR103, Null, -		66																						
*01:42, 01:61, 01:89	-	65	66																						
*03:01:01:01-03:01:05, 03:01:07-03:01:08, 03:01:10-03:01:23, 03:01:25-03:01:37, 03:04:01-03:06, 03:09, 03:11:01, 03:11:03-03:16, 03:18-03:20, 03:22-03:23, 03:25:01-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43-03:45, 03:47-03:48, 03:50-03:52:02, 03:54-03:63, 03:66-03:68N, 03:70-03:73, 03:75, 03:77-03:86, 03:89, 03:91, 03:93, 03:95-03:96, 03:98-03:100:02, 03:104, 03:106-03:110, 03:112-03:114, 03:116-03:118, 03:121-03:124, 03:127-03:129, 03:132-03:139, 03:141-03:149, 03:151-03:153, 03:158-03:161, 03:163-03:166, 03:168-03:170, 03:172, 03:175, 03:180-03:184, 03:186-03:187, 03:191-03:192	DR3, DR17, Null, -					69	70											81							
*03:01:06, 03:01:09, 03:46, 03:64, 03:92, 03:111, 03:130, 03:156N, 03:162, 03:189N	DR17, Null, -					69	70																		
*03:01:24, 03:07:01-01-03:07:02, 03:17, 03:21, 03:24, 03:32, 03:35, 03:39-03:41:02, 03:49, 03:94, 03:97, 03:101, 03:105, 03:120, 03:150, 03:155, 03:171, 03:178, 03:185, 03:188	DR3, -					69													81						
*03:02:01:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74, 03:88, 03:90, 03:102-03:103, 03:115, 03:119, 03:131, 03:154, 03:167, 03:176, 03:179	DR3, DR18, -					69	71												81						
*03:08, 03:65, 03:140	DR3, -					69	70				75		77												
*03:10	DR3					69	70													83	84				
*03:42, 03:69, 03:87, 03:173, 03:190	-					69																			
*03:76, 03:174N	Null, -					69	70					76			79			81							
*03:125	-					69	70					76			79	80		81							
*03:126, 13:193 ⁴	-						70	71										81		83					
*03:157	-					69	70											81			84				
*03:177	-					69	70	71										81		83					
*04:01:01:01-04:05:11, 04:05:13-04:61, 04:63-04:68, 04:70-04:72:02, 04:74:01-04:104, 04:106-04:121, 04:123-04:145, 04:147-04:161, 04:163-04:210, 04:212N-04:250, 04:252-04:283, 04:285-04:336	DR3, DR4, Null, -								72																
*04:62, 04:69, 04:73:01-04:73:02, 04:105:01-04:105:02, 04:122, 04:146, 04:162, 04:251, 04:284:01-04:284:02	DR4, -								72										82						
*04:211	-								72											83					
*07:01:01:01-07:01:01:02, 07:01:01:04-07:01:28, 07:03-07:05, 07:07-07:127, 07:129N-07:133	DR7, Null, -									73															
*07:06	DR7									73	74														
*07:128	-									73									82						
*08:01:01-08:01:02, 08:01:04-08:01:05, 08:01:07-08:02:04, 08:04:01:01-08:07, 08:11:01:01-08:11:01:02, 08:16-08:17, 08:22, 08:24:01-08:24:02, 08:26, 08:28, 08:39, 08:42-08:44, 08:50, 08:52, 08:54-08:55, 08:59, 08:64, 08:72-08:73, 08:77-08:79, 08:82-08:86, 08:88-08:90, 08:93-08:94, 08:96, 08:99-08:101, 08:103, 08:105-08:107, 08:109, 08:111	DR8, Null, -										74							80							
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
*08:01:06, 08:03:02:01-08:03:10, 08:03:12, 08:10, 08:12-08:15, 08:18-08:19, 08:23, 08:25, 08:27, 08:29-08:30:03, 08:33-08:34, 08:36:01, 08:37-08:38, 08:45:01-08:49, 08:51, 08:56, 08:58, 08:60N-08:62, 08:65-08:66, 08:68:02-08:69, 08:71, 08:74, 08:76, 08:81, 08:91, 08:95, 08:97, 08:102, 08:104, 08:110	DR8, Null, -										74													
*08:03:11, 08:68:01	-										74											w		
*08:08	DR8										74						80		82					
*08:09:01	DR8										74						80				84			
*08:09:02, 14:15 ⁴	DR8, -										74						80	81			84			
*08:20, 13:18, 13:47, 13:55, 13:158:01-13:158:02, 13:164, 13:232, 13:286 ⁴	DR13, -															79	80	81			84			
*08:21	DR8										74					79	80				84			
*08:31, 11:05, 11:209 ⁴	DR8, DR11, -													77		80								
*08:32, 08:53	-										74				78						84			
*08:35, 08:36:02	-										74										84			
*08:40	-						70				74													
*08:41, 08:75, 11:67, 11:193:01-11:193:02 ⁴	-										74			77		80								
*08:57, 08:80, 08:98	-																80							
*08:63	-										74			77										
*08:67	-										73						80							
*08:70	-										74						80				w			
*08:87	-										74				78	80					w			
*08:92	-						70				73	74					80							
*08:108	-										73	74					80							
*09:01:02:01-09:01:05, 09:01:07-09:02:02, 09:04-09:22, 09:24-09:28, 09:30-09:51	DR9, Null, -												75							83				
*09:01:06, 09:03, 09:23, 09:29	DR9, -												75											
*10:01:01:01, 10:01:01:03-10:10, 10:12-10:18, 10:20-10:43	DR10, -													76										
*10:11	-										73			76										
*10:19	-													76									86	
*11:01:01:01, 11:01:01:03-11:01:17, 11:01:20-11:01:27, 11:01:29-11:01:45, 11:04:01:01-11:04:16, 11:04:18-11:04:20, 11:04:22, 11:06:01-11:06:03, 11:09:01:01-11:10:02, 11:12:01-11:12:03, 11:15:01:01-11:15:02, 11:24:01-11:24:02, 11:27:01, 11:28:01:01-11:29:02, 11:32-11:33, 11:35, 11:37:01-11:39, 11:43:01:01-11:44, 11:46:01-11:47, 11:49:01-11:51, 11:54:01-11:54:02, 11:56, 11:58:01-11:58:02, 11:60-11:62:02, 11:66:01-11:66:02, 11:72, 11:74:01-11:75, 11:77-11:78, 11:81, 11:84:01, 11:90-11:92, 11:94-11:95, 11:97, 11:99-11:102:02, 11:106, 11:108-11:117:02, 11:120-11:121, 11:123, 11:126, 11:128-11:129, 11:133-11:134, 11:137, 11:140-11:141, 11:144-11:147:02, 11:150, 11:152, 11:154-11:158, 11:160, 11:162-11:167, 11:169N-11:170, 11:172, 11:174-11:175, 11:177-11:181, 11:183, 11:185-11:190, 11:192, 11:194:01-11:201, 11:203-11:208, 11:210-11:212, 11:214-11:228, 11:230, 11:232, 11:234-11:236, 11:238-11:240, 11:243-11:247, 11:249-11:250N, 11:252-11:255, 11:257-11:259, 11:262-11:263, 11:265, 11:269, 11:271-11:275, 11:277-11:279, 11:281-11:282, 11:284-11:285, 11:287N, 11:289, 11:291-11:292, 11:294N-11:298	DR11, Null, -														77	79	80							
*11:01:18-11:01:19, 11:04:21, 11:08:01-11:08:03, 11:18-11:19:03, 11:42:01-11:42:02, 11:57, 11:130, 11:191, 11:248Q, 11:268, 11:280, 11:286	DR11, -														77	79								
*11:01:28	-													77		? 80								

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
----------	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	
*11:02:01:01-11:03:03, 11:11:01, 11:11:03, 11:14:01:01-11:14:02, 11:16:01:01-11:16:01:02, 11:20-11:21, 11:36, 11:40-11:41, 11:48, 11:63:01-11:63:02, 11:65:01-11:65:02, 11:68, 11:70, 11:76, 11:85-11:86, 11:93, 11:118:01-11:118:02, 11:122, 11:124:01, 11:127, 11:138-11:139, 11:151, 11:161, 11:168, 11:171, 11:176, 11:213, 11:231, 11:233, 11:237, 11:241-11:242, 11:251, 11:256, 11:260, 11:276, 11:290	DR11, DR13, -					70								77	79	80									
*11:03:04	-					70								77	79	80	81								
*11:04:17, 11:149	-													77	79	80	81								
*11:07:01-11:07:02, 11:103:01-11:103:02, 11:105, 11:107, 11:125, 11:173	DR11, -					69					75			77											
*11:13:01-11:13:02	DR11						w							77			81	83	84						
*11:17, 11:202	DR14, -													77			81	83	84						
*11:22, 11:53, 11:98, 11:104, 11:143, 11:261	DR11, -													77											
*11:23:01-11:23:02, 11:25, 11:96, 11:148	DR11, -													77	79	80				84					
*11:26, 11:34, 11:267	DR11, -						71							77											
*11:27:02-11:27:03, 11:84:02-11:84:03	DR11, -					69								77	79	80									
*11:30, 11:288	DR11, -				68									77	79	80									
*11:31, 11:45, 11:64, 11:119	DR11, -													77	79								84		
*11:52	DR11													77						83	84				
*11:55	DR11													77		80						84			
*11:59, 11:80, 11:83, 11:87, 11:135, 11:142, 11:182	DR11, -					70						76		77	79	80									
*11:69, 11:82	-													77	80	82									
*11:73, 11:79, 11:131-11:132, 11:153, 11:184	-					70								77	79										
*11:88	-													77	79	80		82							
*11:89, 11:159, 11:266, 11:283	-													77							84				
*11:124:02	-					69	70							77	79	80									
*11:136	-					69								77											
*11:229, 11:270	-					70				73				77		80									
*11:264, 13:05:01:01-13:05:03, 13:07:01-13:07:02, 13:11:01-13:11:02, 13:14:01-13:14:03, 13:21:01:01-13:21:02, 13:42:01-13:42:02, 13:46, 13:49:01-13:50:03, 13:62, 13:100, 13:108, 13:132, 13:136, 13:150, 13:169, 13:189, 13:199, 13:203, 13:229, 13:246, 13:284, 13:302 ⁴	DR11, DR13, DR6, -															79	80	81							
*11:293, 13:45 ⁴	DR13, -					70										79	80		82	84					
*12:01:01:01-12:01:11, 12:03:02-12:04, 12:06-12:12, 12:14, 12:17, 12:24N-12:25, 12:28-12:30, 12:34-12:36, 12:38, 12:40-12:41, 12:46-12:48, 12:53-12:54, 12:58-12:59, 12:62-12:63, 12:68, 12:70, 12:72N, 12:77, 12:79, 12:82-12:84, 12:86N-12:87, 12:89-12:91	DR12, Null, -														78										
*12:02:01:01-12:02:11, 12:13, 12:18-12:19, 12:21, 12:23, 12:26-12:27, 12:31N-12:33, 12:37, 12:42-12:45, 12:49-12:52, 12:55-12:56, 12:60N-12:61, 12:64-12:65, 12:67, 12:69, 12:71, 12:74N-12:76, 12:78, 12:80-12:81, 12:85, 12:92-12:94	DR12, Null, -														78	80									
*12:05	DR12														78								85		
*12:15	DR12														78	80							85		
*12:16:01-12:16:03	-									74					78	80									
*12:20	-														78	80	81								
*12:22	-									73	74				78										
*12:39:01-12:39:02	-									74					78										
*12:57, DRB3*02:111 ⁴	-				68										79	80							85		
*12:66	-														78					83					
*12:73	-					70				73					78	80									
*12:88	-									73					78										
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	
*13:01:01:01-13:02:01:11, 13:02:03-13:02:18, 13:04, 13:08, 13:16:01:01-13:16:01:02, 13:20:01:01-13:20:01:03, 13:22:01:01-13:24:01:03, 13:28:01-13:29:01:02, 13:31-13:32, 13:34-13:36, 13:38-13:40, 13:48, 13:51-13:52, 13:54, 13:59, 13:61:01, 13:63-13:65, 13:68-13:70, 13:72-13:76, 13:78-13:80, 13:83-13:84, 13:87, 13:91-13:93, 13:96:02-13:99, 13:102, 13:105-13:107, 13:109, 13:111-13:114, 13:117:01:01-13:117:01:02, 13:121, 13:123-13:128, 13:130-13:131, 13:135, 13:138-13:143, 13:145, 13:147-13:149, 13:153, 13:155, 13:162, 13:166, 13:168, 13:171:02, 13:173, 13:177, 13:182, 13:184-13:187, 13:190, 13:200N-13:202, 13:204-13:205, 13:207-13:209, 13:212-13:213, 13:215, 13:218, 13:221, 13:223, 13:225-13:226, 13:230, 13:233-13:234, 13:236-13:245, 13:247, 13:249N-13:252N, 13:254-13:276, 13:278Q-13:283, 13:285, 13:287-13:288, 13:290-13:294, 13:296-13:298N, 13:300-13:301, 13:303-13:313, 13:315-13:316, 13:318, 13:320	DR11, DR13, DR14, Null, -						70										79	80	81						
*13:02:02, 13:289N	DR13, Null						?											81							
*13:03:01:01-13:03:11, 13:10, 13:33:02-13:33:03, 13:37, 13:66:01-13:66:02, 13:81, 13:88-13:90, 13:94:02-13:95, 13:101, 13:115, 13:120, 13:133, 13:137N, 13:151-13:152, 13:170, 13:172, 13:174, 13:178, 13:180, 13:188, 13:194, 13:196, 13:216, 13:219, 13:227, 13:253, 13:317, 13:319N	DR13, Null, -						70										79		81						
*13:06, 13:12:01-13:12:05, 13:25, 13:30, 13:56, 13:58, 13:60, 13:77, 13:82, 13:118, 13:134, 13:157, 13:161, 13:163, 13:181, 13:220, 13:224, 13:231, 13:295N	DR11, DR13, DR6, Null, -																79		81						
*13:09, 13:192, 13:210-13:211, 13:222	DR13, -																80	81							
*13:13, 13:119, 13:154, 13:156, 14:84, 14:116, 14:144 ⁴	DR13, -															79		81			84				
*13:15, 13:19, 13:53, 13:57, 13:104, 13:198	DR13, -						70	71									79	80	81						
*13:17, 13:116, 13:175, 13:214	DR13, -						70			73							80								
*13:26:01-13:26:02	DR14, -							71									79	80	81						
*13:27:01-13:27:02, 13:41, 13:71, 13:277	DR13, -						70						76				79	80	81						
*13:33:01, 13:94:01	DR13, -					69	70										79		81						
*13:43, 13:159, 13:179, 13:191, 13:217	DR13, -						70										79	80			83	84			
*13:44, 13:86, 13:206, 14:183 ⁴	-							71											81						
*13:61:02, 13:96:01	-					69	70										79	80	81						
*13:67, 13:103, 13:160, 13:195	DR13, -																79	80							
*13:85	-						70	71									79		81						
*13:110, 13:167	-						?										79		81						
*13:122	-						70										79								
*13:129	-						70						76					80	81						
*13:144	-						70										79	80	81			84			
*13:146	-																	80	81			84			
*13:165	-						70										79	80							
*13:171:01	-						70			74							79	80			83	84			
*13:176	-						70						76				79		81						
*13:183	-						70											80	81						
*13:197	-									74							79	80			82	84			
*13:228	-																79								
*13:235, 14:98 ⁴	-							71									79		81						
*13:248	-						70											80							
*13:299	-						70																		
*13:314, 14:37 ⁴	DR14, -																	80	81	82					
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
*14:01:01, 14:01:05, 14:04:01:01-14:04:01:03, 14:04:04-14:04:06, 14:04:09, 14:07:01-14:07:02, 14:10:01:01-14:10:01:02, 14:54:01:01-14:54:01:12, 14:54:03, 14:54:05-14:54:06, 14:54:08-14:54:11, 14:57, 14:70, 14:137N, 14:152N, 14:191-14:192, 14:196-14:197N, 14:199, 14:201-14:202, 14:204-14:205:01, 14:207, 14:210Q, 14:213-14:216, 14:219, 14:224, 14:227-14:228, 14:230-14:234, 14:237-14:238, 14:240	DR14, DR1404, DR4, DR6, Null, -																	81	82	83	84			
*14:01:02, 14:01:04, 14:04:02-14:04:03, 14:26, 14:28, 14:31, 14:32:03, 14:35, 14:38:01, 14:39, 14:55, 14:60-14:62, 14:71, 14:75-14:76, 14:79, 14:86-14:88, 14:90, 14:99, 14:101, 14:107, 14:110, 14:112-14:114, 14:117-14:118, 14:124-14:125, 14:129, 14:131, 14:138-14:140, 14:142-14:143, 14:145, 14:147, 14:149-14:151, 14:153, 14:157-14:158, 14:163-14:164, 14:166N, 14:172, 14:180, 14:182, 14:184, 14:189-14:190, 14:203, 14:220, 14:225, 14:229	DR14, Null, -																	?	82	83	84			
*14:01:03, 14:23:02, 14:34, 14:72, 14:92N, 14:97, 14:126:01-14:126:02, 14:130, 14:133, 14:148, 14:175, 14:185	DR14, Null, -																	?		83	84			
*14:02:01:01-14:02:09, 14:06:01-14:06:06, 14:09, 14:13, 14:17, 14:20, 14:29-14:30, 14:33, 14:41, 14:47-14:48, 14:51, 14:80, 14:83, 14:94, 14:106, 14:108, 14:121, 14:159, 14:165, 14:170, 14:176-14:177, 14:194-14:195N, 14:198, 14:209, 14:212, 14:221	DR14, DR6, Null, -						71											81		83				
*14:03:01-14:03:02, 14:12:01-14:12:02, 14:40, 14:63, 14:67, 14:77-14:78, 14:85, 14:102, 14:115, 14:135, 14:156, 14:174, 14:181, 14:218, 14:223, 14:239	DR14, DR1403, DR6, -						71								79		81				84			
*14:04:07-14:04:08, 14:54:07	-																		82	83	84			
*14:05:01:01-14:05:05, 14:08, 14:11:01:01-14:11:01:02, 14:14, 14:23:01:01-14:23:01:03, 14:23:03-14:23:04, 14:36, 14:43-14:44:03, 14:54:02, 14:56, 14:59, 14:64, 14:91, 14:96, 14:100, 14:103, 14:123, 14:127:01-14:127:02, 14:134, 14:155, 14:162, 14:167, 14:169, 14:187-14:188N, 14:193, 14:205:02-14:206, 14:211, 14:236	DR14, Null, -																	81			83	84		
*14:16	DR6					70								79	80		81	82	83	84				
*14:18, 14:81, 14:154, 14:160	DR14, -						71											81		83	84			
*14:19, 14:21, 14:109	DR14, -					70	71								79			81		w				
*14:22, 14:105	DR14, -														79	80		?	82	83	84			
*14:24, 14:226	DR14, -						71										80	81						
*14:25:01-14:25:02, 14:128	DR6, -														79	80		?	82		84			
*14:27:01-14:27:02	DR14, -						71								79	80		81			84			
*14:32:01-14:32:02	DR14						w											?	82	83	84			
*14:38:02, 14:104	-				69													?	82	83	84			
*14:42, 14:136, 14:161	-																	81			84			
*14:45, 14:95	DR14, -					70												81		83	84			
*14:46	DR14																			83				
*14:49, 14:119, 14:146	DR14, -						71											?	82	83	84			
*14:50	DR14									73								?	82	83	84			
*14:52	DR14														78			?		83				
*14:53	DR13, 14															79	80				82	84		
*14:54:04	-																	81	82	w	84			
*14:58, 14:186	DR14, -																	?	82		84			
*14:65	DR6						w											81			83	84		
*14:68:01-14:68:02, 14:93	DR14, -									74								?	82	83	84			
*14:69	-															79		?	82		84			
*14:73	-																80	?	82	83	84			
*14:74, 14:173	-															79		?	82	83	84			
*14:82	-					70												81	82	83	84			
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87

A-B-DR-DQ 2022 0R5 DR

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
*14:89	-							71										81		84				
*14:111, 14:171	-					69												81	82	83	84			
*14:120, 14:122, 14:178	-						70											?	82	83	84			
*14:132	-						70											?		83	84			
*14:141	-																	81				85		
*14:168	-																	?	w	83	84			
*14:179	-						70	71					76					81		83				
*14:200	-							71												83	84			
*14:208	-																			82	84			
*14:217	-																	81	82	83	84		86	
*14:222N	Null							71									79	?			84			
*14:235	-																79	81	82	83	84			
*15:01:01:01-15:20, 15:22-15:24, 15:26, 15:28-15:33, 15:35-15:36, 15:37:02-15:53, 15:55-15:65, 15:67-15:99, 15:101-15:103, 15:104:02-15:121, 15:123-15:140, 15:142-15:157, 15:159N-15:166, 15:168-15:186, 15:188-15:198	DR15, DR2, Null, -				67																			
*15:21	DR2				67																w			
*15:25, 15:37:01, 15:100, 15:104:01	-				67	69																		
*15:27, 15:34, 15:54, 15:66:01-15:66:02, 15:122	-				67															83				
*15:141	-				67																		86	
*15:158, 15:187	-				67																83	w		
*15:167	-				67															82				
*16:01:01-16:03, 16:05:01-16:05:02, 16:07-16:17, 16:19-16:65, 16:67-16:71	DR16, Null, -				68																			
*16:04:01-16:04:02, 16:18	DR16, -				68																	w		
*16:66	-				67	68																		
DRB3*01:01:02:01-01:41, DRB3*01:43:01:99, DRB3*02:01:02:26, DRB3*02:28:02:110, DRB3*02:112-02:130, DRB3*02:132:02:159, DRB3*02:161-02:175, DRB3*03:01:01:01-03:21, DRB3*03:23-03:55	DR52, Null, -																						85	
DRB3*01:42	-											75											85	
DRB3*02:27	-															79							85	
DRB3*02:131	-															79					84			
DRB3*02:160	-															79	80						85	
DRB4*01:01:01:01, DRB4*01:01:01:03-01:03:01:01, DRB4*01:03:01:03-01:03:01:12, DRB4*01:03:01:14-01:13, DRB4*01:15-01:30, DRB4*01:32-01:154, DRB4*02:01N, DRB4*03:01N	DR53, Null, -																							86
DRB4*01:03:01:02N, DRB4*01:03:01:13N, DRB4*01:14N	Null																			82				86
DRB4*01:31	-																				84			86
DRB5*01:01:01:01-01:117, DRB5*02:02:01-02:30	DR51, Null, -																							87
Well No.	ser. ³	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87