















*15:517	-							29			32							39							53		56							64																		
*15:518	-						28				33						39									56							64																			
*15:524	-			24			28										39									56							64																			
*15:525	-						28	29									39									55	56					62	64																			
*15:545	-						28					35					39									56						62	64																			
*15:546Q	-																						52									63																				
<b>Well No.</b>		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64								
*15:552	-						28										39																				56		59													
*15:578	-						28										39										48									56	57															
*15:588	-				25			29			32				36		39																			56																
*15:592	-						28				33					37	39											48		51																						
*15:594	-						28									37																			56								62	63								
*15:634	-										33												44												56								62	64								
*18:01:01:01-18:01:36, 18:01:38-18:03:01:03, 18:04:02-18:08, 18:10-18:11, 18:13, 18:15, 18:17N-18:18:01:02, 18:20-18:21, 18:23N-18:28, 18:30-18:32, 18:34-18:36, 18:38-18:43, 18:45-18:47, 18:49-18:53, 18:55, 18:57:01-18:65, 18:68, 18:70-18:71, 18:73-18:79, 18:81-18:91, 18:93-18:98:01, 18:99-18:101, 18:103-18:104, 18:106Q, 18:109-18:124, 18:126-18:131:01:02, 18:133-18:135:01:02, 18:137-18:148, 18:151-18:163, 18:167, 18:170-18:174, 18:176-18:177, 18:180-18:190, 18:192-18:193, 18:196-18:203, 18:206-18:224	B18, Null, -													35		38																									61					64						
*18:01:37	-						26							35		38																									61					64						
*18:03:02, 18:66	-													35																																64						
*18:04:01, 18:98:02, 18:191	B18, -													35		38																							58		61					64						
*18:09	B18													35																											61			63								
*18:12:01-18:12:02	B18, -													35		38																			56											64						
*18:14, 18:44:01-18:44:02	B18, -				25									35		38																									61					64						
*18:19, 18:164	B18, -													35		38																									61	62			64							
*18:22, 18:69, 18:105, 18:149-18:150, 18:204	B18, -													35	37	38													49												61					64						
*18:29, 18:72:01-18:72:03, 18:92, 18:175, 35:01:55, 35:01:62, 35:31-35:32:03, 35:37, 35:53N, 35:64:01-35:64:02, 35:68:01-35:68:02, 35:75, 35:99, 35:118-35:119, 35:127, 35:151:01-35:151:02, 35:174, 35:205:01-35:205:03, 35:234-35:235, 35:251:02, 35:273, 35:292, 35:367, 35:369, 35:409, 35:433, 35:521, 35:523, 35:535 <sup>4</sup>	B35, Null, -															38																												61					64			
*18:33	-													35		38																			54						61					64						
*18:37:01-18:37:02	-										32			35																																64						
*18:48	-													35																					55											64						
*18:54	-										33			35						40																					61			63								
*18:56, 18:108, 18:132, 18:169	-										33			35		38																									61								64			
*18:67, 18:136, 18:194:01:01-18:194:01:02	-													35															54												61			63								
*18:80	-													35		38																						55			61								64			
*18:102	-													38																								55			61								64			
*18:107	-				25											38																						55			61								64			
*18:125	-										32			35																								55			58								64			
*18:165	-					26								35		38																			50						61								64			
*18:166	-						27							35		38																									61								64			
*18:168	-													35		38			41																						61								64			
*18:178	-													35																								56								64						
*18:179	-													35		38																									61								64			
*18:195	-										33			35					40																			56						61			63					
*18:205	-									30				35		38																									61								64			
<b>Well No.</b>		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64								
*27:01	B27				26									36														47		w						56											63					







*35:35, 35:468	B35, -																	37	38		41														61		64											
*35:36	B35																		37	38															61		63	64										
*35:38, 35:115, 35:169, 35:309, 35:452, 35:513N	B35, Null, -																		37	38																61		64										
*35:43:01, 35:43:03-35:44, 35:67, 35:79, 35:86, 35:102, 35:117, 35:476, 35:549	B35, -																			38																61	62	64										
*35:45, 35:71	B35																			37	38															61		64										
*35:46	B35																			37																	64											
*35:47, 35:154, 35:217, 35:450	B35, -																			37																	64											
*35:60	B35																			38																	58	61	64									
*35:63, 35:317, 40:129 <sup>4</sup>	-																			37																	55	56	64									
*35:66	B35																				38																	61	64									
*35:76, 35:561	B35, 22, -																				37	38																w	64									
*35:82, 35:85, 35:294	-																				37	38																61	64									
*35:84, 35:90, 35:150:01-35:150:02, 35:184, 35:201, 35:441	-																				37	38																61	64									
*35:87	-																				38																		61	64								
*35:135	-																				38																		61	62	64							
*35:137, 35:147, 35:182, 35:278, 35:540	-																				37	38																	61	64								
*35:185	-																				38																		58	61	62	64						
*35:186	-																				37	38																	64									
*35:207	-																				37																		56	64								
*35:218, 35:256	-																				37	38																	61	64								
*35:229	-																				37	38																	61	64								
*35:230	-																				37	38																	61	64								
*35:251:01	-																				38																		61	64								
*35:260	-																				38																			61	64							
*35:265	-																				38																		61	62	64							
*35:276	-																				38																			58	61	64						
*35:289, 35:487	-																				35	37	38																61	64								
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64				
*35:311	-																				37																									64		
*35:329	-																				37	38																							61	63	64	
*35:331	-																				38																									64		
*35:361	-																				37																									64		
*35:368	-																				38																								64			
*35:432, 35:446	-																				37	38																							64			
*35:471	-																					41																								61	63	64
*35:475	-																				38																									61	62	64
*35:516	-																				37	38																								64		
*35:525	-																				37	38																								61	64	
*35:530	-																				37	38																								64		
*35:534	-																				38																									64		
*35:541	-																				38																									64		
*35:545	-																				37	38																									64	
*37:01:01-01-37:01:07, 37:01:09, 37:01:11, 37:01:13-37:01:29, 37:03N, 37:06:01-37:06:02, 37:13, 37:15-37:21, 37:23-37:27, 37:30N-37:33N, 37:35-37:36, 37:38, 37:41-37:43, 37:45-37:50, 37:52-37:56, 37:58-37:59, 37:61, 37:63-37:65, 37:68-37:80, 37:82N, 37:85-37:93, 37:98-37:99, 37:101-37:105	B37, Null, -																				33																										63	
*37:01:08, 37:01:10, 37:09, 37:44, 37:51	B37, -																				33																										63	
*37:01:12, 37:10, 37:22:01-37:22:03, 37:95, 37:100	B37, -																				37																										63	
*37:02	B37, B27																				26																									63		
*37:04:01-37:04:02	B37, -																				33																									63		
*37:05	-																				33	37	38																								64	
*37:07	-																				22																										63	
*37:08	B37																				33																									61	63	64















*44:06	B44																					37								48					54		57			61	63									
*44:08	B44		w	25																			37							48						56	57				63									
*44:09, 44:254	B45, -			25																			37							48					55	56	57				64									
*44:10, 44:140, 44:340	B44, -			25							32												37							48						56			59			63								
*44:14	B44			25				29															37							48						56					63									
*44:15:01:01-44:15:01:02	B12, -			25			29			32													37							48						56			59			63								
*44:16, 44:91, 44:132	B44, B47, -			25		27																	37							48						56	57					63								
*44:17, 44:43:01-44:43:02, 44:144, 44:256, 44:308	B44, -			25							33												37							48						56						63								
*44:18:01:01-44:18:01:02, 44:355	B44, -			25				29		32													37							48					52		56		59			63								
*44:20, 44:100	B44, -			25				29															37							48						56	57					63								
*44:25, 44:50:01, 44:50:03	B44, -			25																			37							48					52	54	56					63								
*44:31, 44:182	B44, -			25																			37							48						45	48					63								
*44:37:01-44:37:03, 44:64:01-44:64:02, 44:268, 44:347	B44, -			25		27																	37							48							56					63								
*44:40, 44:210:01-44:210:02	B44, -			25							32												36	37						48												63								
*44:44	B44			25							32												36	37					48									56	57				63							
*44:46:01-44:46:02, 44:75, 44:129, 44:277, 44:438N	B12, Null, -			25																			37							48							55	56					64							
*44:47, 44:65, 44:192:01-44:193, 44:197, 44:305, 44:440, 44:491, 44:506	B44, -			25																			37							48													63							
*44:49	B44		23																				37							48								56	57					63						
Well No.			22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64					
Well No.			22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64					
*44:54	B44			24	25																									48																63				
*44:57	-		w	25																			37							48															63					
*44:60, 44:110	-			24	25																		37							48															63					
*44:62	-			25						29	30												37							48																63				
*44:90	-			25																			37							47	48						55	56	57						64					
*44:95	-			25																			37							48						52			56	57					63					
*44:97	-			25																			37							47	48					50									63					
*44:106	-			24	25																		37							48															63					
*44:117, 44:251	-			25										32									37							48															63					
*44:130	-													32									36	37						48															63					
*44:131	-																						37							48								55	56	57					64					
*44:146, 44:299	-				25																		37							48														62	63					
*44:150	-		22		25			27															37							48																63				
*44:156	-			23										32									37							48																63				
*44:166	-		22		25					29													37								48															63				
*44:169	-				25																		37								48															63				
*44:190	-				25					27													37								48	49														63				
*44:203:01-44:203:02	-				25																		37								48	49														63				
*44:211	-				24	25																	37								48														63					
*44:213	-				24	25				29													37									48														63				
*44:230	-				25																		37									48														63				
*44:231	-				w																		37								48														63					
*44:236	-				25																		37								48															63				
*44:263	-				25																		37								48							54								63				
*44:275	-				25					27													37									48													62	63				
*44:276	-				25																		37									48														63				
*44:290	-				25					27													37									49														63				
*44:339	-				25																		37									48														63				
*44:342	-			22		25																	37								46	48														63				
*44:345N, 44:458	Null, -				25																		37								48															63				
*44:356	-				25																		37								47	48														64				
*44:445	-				25																		37								48															63				
*44:459	-				25								30										37																							63				
*45:01:01:01-45:01:08, 45:01:11-45:01:12, 45:05, 45:07, 45:11-45:15, 45:17-45:21, 45:23, 45:25, 45:27-45:28N, 45:30	B45, Null, -			25						29																						44													55	56	59			64

*45:01:09, 45:16	-			25			29			32																		56		59																														
*45:01:10	-			25			29			32																		55	56		59			64																										
*45:02-45:03, 45:29	B45, -			25			29			32																		55	56					64																										
*45:04	B21			25			29																					55	56		59			64																										
*45:06	B45			25			29			32																		55	56		59			64																										
*45:08	-			25			29																					55	56		58			64																										
*45:09	-			25			29			33																		55	56					64																										
*45:10	-			25			29		31																			55	56					64																										
*45:22	-			25			29			32					38													55			59	61		64																										
*45:24	-			25			29																					55	56		w			64																										
Well No.				22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64														
Well No.				22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64														
*46:01:01:01-46:01:39, 46:04-46:05, 46:07N, 46:09, 46:14-46:16, 46:20, 46:22-46:24, 46:26-46:31, 46:34-46:39, 46:42, 46:44, 46:46, 46:48, 46:51Q-46:53, 46:55N-46:57, 46:60-46:61, 46:63-46:66, 46:68-46:70, 46:73-46:76, 46:78-46:82, 46:85, 46:87-46:88, 46:90-46:92, 46:94-46:96	B46, Null, -							28																					47	49															62		64													
*46:02	B46							28																				w	49																62		64													
*46:03, 46:08, 46:13:01-46:13:03, 46:19, 46:21:01-46:21:02, 46:25, 46:45, 46:54, 46:58-46:59, 46:62, 46:83, 46:93	B46, -							28																					47	49																			64											
*46:06, 46:43	-							28																				45	47	49																		64												
*46:10	B46							28																					47	49								55							62		64													
*46:11	B22							28																					47	49																	64													
*46:12	-			24	25			28																						47	49											59					62		64											
*46:17	-							28							33														47	49																62		64												
*46:18	-							28																					47	49																		64												
*46:32:01-46:32:02	-							28																					47	49					w											62		64												
*46:33	-							28																					47	49																		64												
*46:40:01	-							28																					38	47																		62		64										
*46:40:02	-							28																					38	47																		62		64										
*46:41N, 46:49	Null, -							28																																									62		64									
*46:47, 46:72, 46:86	-							28																					47																				62		64									
*46:50	-							28																					47	49																			62		64									
*46:67	-							28								37													47	49																			62		64									
*46:71	-							28																					47	49																			62		64									
*46:77	-							28																					47	49																			62		64									
*46:84	-						25		28																				47	49																			62		64									
*46:89	-																										44		47	49																				62		64								
*47:01:01:02-47:01:06, 47:06-47:12	B47, -						25																																												50		63							
*47:02	B47						25																																											50		64								
*47:03	B47						25																																												50		w	63	w					
*47:04-47:05	B17, -						25	26																																										50		56		63						
*48:01:01:01-48:01:10, 48:07, 48:09, 48:15, 48:19-48:20, 48:22, 48:27, 48:30, 48:32-48:38, 48:42-48:50, 48:52, 48:54	B48, -			22										32																																					51		56			64				
*48:02:01-48:02:03, 48:25	B48, B70, -													32																																						37		56		64				
*48:03:01:01-48:03:02, 48:23, 48:40	B48, -													32																																							44	45		56		64		
*48:04:01:01-48:04:02, 48:24, 48:51, 48:53	B48, -													32																																							44		51		56		64	
*48:05, 48:08, 48:10	B48, -			22										32																																							44		56		64			
*48:06	-			22																																																	44		51		56		64	
*48:11	B48			22										32																																						44		56		64				
*48:12, 48:14	B48			22			25							32																																						44		56		64				
*48:13	-																																																				44		56		64			
*48:16, 48:39	B48, -			22										32																																							44		51		56		64	
*48:17	-			22																										44	45																							49		56		64		
*48:18	B48			22																																																			51	52		56		63
*48:21	-						26																																																44		56		64	

*48:26	-				26				32					44			50		56				64																											
*48:28	-	22							32	35				44			51		56				64																											
*48:29	-	22	25						32					44			51		56				64																											
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
*48:31	-	22															44													51			56												64					
*48:41	-	22															44													51			55	56												64				
*49:01:01:01-49:01:01:18, 49:01:03-49:01:07, 49:01:09-49:01:13, 49:01:15-49:01:23, 49:06, 49:08-49:10, 49:12-49:14, 49:16-49:17, 49:19N, 49:22-49:24, 49:27-49:30, 49:32-49:44, 49:46, 49:48-49:49, 49:51-49:60N, 49:64-49:66, 49:69, 49:71-49:72, 49:74-49:75, 49:78	B49, Null, -				25																										52	53		56		59					63									
*49:01:02, 49:11, 49:45, 49:47:01-49:47:02	B49, -				25																										52		56		59											63				
*49:01:08, 49:04:01-49:05:02, 49:18:02, 49:50	B49, -				25																									52	53		56													63				
*49:01:14	-				25												44													52	53		56		59												63			
*49:02, 49:70	B49, -				25																									52	53		56		59												63			
*49:03, 49:31, 49:68, 49:77	B49, -				25																									52	53		56		59													63		
*49:07	-				25						31																			52	53		56														63			
*49:18:01	-				25					30																				52	53		56															63		
*49:20	-				25						32																			52	53		56		59													63		
*49:21	-				25																									52	53		56		59														63	
*49:25	-	23																												52	53		56		59														63	
*49:26	-				25											36		42	43											52	53		56		59														63	
*49:61	-				25	27											37													52		56																	63	
*49:62	-				25																								52	53	54	56		59															63	
*49:63, 49:67	-				25																								52	53		56		59															63	
*49:73	-				25						32				36														52	53		56		59															63	
*49:76	-				25																								52	53		56		58															63	
*50:01:01:01-50:01:17, 50:01:19-50:01:20, 50:01:22-50:01:23, 50:04:01-50:08, 50:10-50:11, 50:13, 50:16, 50:18-50:19, 50:32, 50:35-50:38, 50:40-50:42, 50:44, 50:46, 50:48, 50:50, 50:52-50:55, 50:57-50:65, 50:67-50:69, 50:72N-50:74, 50:76-50:79	B50, Null, -				25																		44							53	55	56		59														64		
*50:01:18	-				25																									53	55	56		59															64	
*50:01:21	-				25					30								44												53	55	56		59																64
*50:02:01:01-50:02:01:02, 50:21, 50:56	B45, -				25						32							44											53	55	56		59																64	
*50:09, 50:20, 50:33, 50:39	-				25													44											53	55	56		59																64	
*50:12, 50:51, 50:66, 50:70	-				25													44											55	56		59																	64	
*50:14	-				25						30							44											53	55	56		59																64	
*50:15	-				25													44											53	55		59																	64	
*50:17	-				25												37	44											55	56		59																	64	
*50:31	-				25						31							44											53	55	56		59																64	
*50:34, 50:43, 50:71	-				25													44										53	55	56		59																	64	
*50:45	-				25													44										53	55	56	58		59																64	
*50:47	-				25													44										w	55	56		59		62														64		
*50:49	-				25													44					41						55	56	58		59																64	
*50:75	-				25							32			36			44											53	55	56		59																64	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64						



*51:274	-					27	29	30							48		53	54		57			61	63																					
*51:310, 52:88, 53:28 <sup>4</sup>	-							30			37				49		52	54	56				61	63																					
*51:341	-							30															61	63																					
*51:367	-					28		30		34		39						53					61	63																					
*51:371	-						29	30										53					61	63																					
*52:01:01-01-52:01:04, 52:01:06-52:01:09, 52:01:11-52:01:14, 52:01:16-52:01:28, 52:01:30-52:01:49, 52:01:51-52:13, 52:14:02-52:15, 52:17-52:18, 52:22-52:24, 52:27-52:36, 52:38-52:53, 52:55-52:57, 52:59-52:66, 52:68-52:70, 52:72, 52:74-52:76, 52:78-52:83, 52:85-52:87, 52:89N, 52:91-52:96N, 52:98-52:105, 52:107-52:109	B52, Null, -							30									52	53	54	56				63																					
*52:01:05, 52:01:15, 52:14:01, 52:58, 52:67	B52, -											42					52	53	54	56				63																					
*52:01:10	-							30									52	53	54	56				63																					
*52:01:29, 52:20	-							30										53		56				63																					
*52:01:50, 52:54	-							30				44					52	53	54	56				63																					
*52:16	-							30	32								52	53	54	56				63																					
*52:19	-						29	30									52	53	54	56				63																					
*52:21:01-52:21:02	-							30									52	53	54	56			62	63																					
*52:25:01	-																52	53	54	56		59		63																					
*52:25:02	-							30									52	53	54	56		59		63																					
*52:26	-						27	30									52	53	54	56				63																					
*52:37	-							30									52	53	54	56		58		63																					
*52:71, 52:84	-							30				41					52	53	54	56		58		63																					
*52:77	-			23				30									52	53	54	56				63																					
*52:90	-							30									52	53		56				63																					
*52:97	-																52	w	54	56			62	63																					
*52:106	-				25			30									52	53	54	56				63																					
*53:01:01:01-53:01:17, 53:01:19-53:01:29, 53:04-53:05:02, 53:07-53:08:02, 53:10, 53:16, 53:18-53:21, 53:23-53:27, 53:29, 53:32-53:35, 53:37, 53:40-53:43:02, 53:45-53:46, 53:48N-53:54, 53:56-53:61, 53:63-53:66, 53:68-53:69	B53, Null, -										37								54				61	63																					
*53:01:18, 53:03, 53:09, 53:11-53:13, 53:31, 53:36, 53:55	B53, B35, -										37												61	63																					
*53:02, 53:06	B53, B51										37				49				54				61	63																					
*53:14	-			24	25														54				61	63																					
*53:15	-			22															54				61	63																					
*53:17:01-53:17:02	-																52	54	56					63																					
*53:22	-									33						48			54				61	63																					
*53:30, 58:111 <sup>4</sup>	-				25												52	54						63																					
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
Well No.		22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	
*53:38	-																37																			56								63	
*53:39, 57:128, 58:100 <sup>4</sup>	-																37															52	54											63	
*53:67	-				24	25																																					61	63	
*53:70	-																37																	54							61	62	63		
*54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10, 54:13-54:19, 54:21-54:25, 54:27-54:29, 54:31-54:32, 54:34, 54:36-54:37, 54:39, 54:41-54:46	B54, Null, -																	47		50																		57	58						64
*54:03	B22																	47		50																	57		59					64	
*54:06	-																	47		50																	57					62			64
*54:09	-																	47		50																	57								64
*54:11	B54																			50																57	58							64	
*54:12	-																	47		50		52													57	58							63		
*54:20	-																47		50			53												57	58									64	
*54:26	-													31				47		50															57									64	
*54:30	-																	47		50														57	58										64
*54:33	-																	47		50		w												57			59			62			64		
*54:35	-																	47		50															56	57	58								64









