

Well No.	1	2
HLA-B allele ¹		
*2701-270402, 270502-270508, 270510, 2708, 2710, 2713, 2715, 2717, 2725, 2728, 2736-2740, 2742, 2744, 2745	1	2
*2706, 2707, 2709, 2711, 2714, 2719-2721, 2724, 2727, 2730, 2732-2735, 2741, 2743	1	
*270509, 2712, 2716, 2718, 2723, 2726, 2729, 2731, 3702, 4704, 4705		2
Length of spec. PCR product	145	95
Well No.	1	2
Specific Reaction		
1st match: B*		
2nd match: B*		
Control Reaction		
Failed Controls		
False Positive		
False Negative		