

|   |   |   |
|---|---|---|
| Well No.  | 1 | 2 |
| HLA-B allele <sup>1,2</sup>   |   |   |
| *27:05:09, 27:05:23, 27:12, 27:16, 27:18, 27:23, 27:26, 27:29, 27:31, 27:75, 27:77, 27:85, 27:92, 27:101, 27:119, 27:129, 27:140, B*07:197, B*15:342, B*37:02, B*47:04-47:05, B*48:26   |   | 2 |
| *27:01-27:05:08, 27:05:10-27:05:15, 27:05:17-27:05:22, 27:05:24-27:05:31, 27:08, 27:10, 27:13, 27:15, 27:17, 27:25, 27:28, 27:36-27:40, 27:42, 27:44-27:45, 27:47-27:69, 27:71-27:74, 27:79-27:80, 27:82-27:84, 27:86-27:90:03, 27:93-27:100, 27:103-27:105, 27:108-27:118, 27:120-27:124, 27:126, 27:128, 27:131-27:135, 27:137, 27:139, 27:141-27:142 | 1 | 2 |
| *27:05:16, 27:06-27:07:04, 27:09, 27:11, 27:14, 27:19-27:21, 27:24, 27:27, 27:30, 27:32-27:35, 27:41, 27:43, 27:46, 27:70, 27:76, 27:78, 27:81, 27:91, 27:102, 27:106-27:107, 27:125, 27:127, 27:130, 27:136, 27:138, B*44:97   | 1 |   |
| Well No.  | 1 | 2 |
|   |   |   |