

B27 bulk 2018 2G7

Well No.	1	2
HLA-B allele ^{1,2}		
*27:05:09, 27:05:23, 27:18, 27:23, 27:26, 27:29, 27:31, 27:75, 27:77, 27:85, 27:101, 27:129, 27:140, 27:157, B*07:197, B*15:342, B*37:02, B*47:04-47:05, B*48:26, A*32:89		2
*27:01-27:05:08, 27:05:10-27:05:15, 27:05:17-27:05:22, 27:05:24-27:05:32, 27:08, 27:10, 27:12:01:01-27:13, 27:15-27:17, 27:25, 27:28, 27:36-27:40, 27:42, 27:44-27:45, 27:47-27:69, 27:71-27:74, 27:79-27:80, 27:82-27:84, 27:86-27:90:04, 27:92-27:100, 27:103-27:105, 27:108-27:124, 27:126, 27:128, 27:131-27:135, 27:137, 27:139, 27:141-27:149, 27:151-27:152, 27:155-27:156, 27:158-27:163	1	2
*27:05:16, 27:06-27:07:06, 27:09, 27:11, 27:14, 27:19-27:21, 27:24, 27:27, 27:30, 27:32-27:35, 27:41, 27:43, 27:46, 27:70, 27:76, 27:78, 27:81, 27:91, 27:102, 27:106-27:107, 27:125, 27:127, 27:130, 27:136, 27:138, 27:150, 27:153-27:154, 27:164, B*38:22, B*40:75, B*40:322, B*44:97, B*44:263, B*49:26	1	
Well No.	1	2