

B27 unit dose 2019 3H1

Well No.	1	2
HLA-B allele <sup>1,2</sup>		
*27:01-27:05:08, 27:05:10-27:05:15, 27:05:17-27:05:22, 27:05:24-27:05:36, 27:08, 27:10, 27:12:01:01-27:13, 27:15-27:17, 27:25, 27:28, 27:36-27:40, 27:42, 27:44-27:45, 27:47-27:69, 27:71-27:74, 27:79-27:80, 27:82-27:84, 27:86-27:90:04, 27:92-27:100, 27:103-27:105, 27:108-27:124, 27:126, 27:128, 27:131-27:135, 27:137, 27:139, 27:141-27:149, 27:151-27:152, 27:155-27:156, 27:158-27:163, 27:165-27:167, 27:169-27:179, 27:181-27:182	1	2
*27:05:09, 27:05:23, 27:18, 27:23, 27:26, 27:29, 27:31, 27:75, 27:77, 27:85, 27:101, 27:129, 27:140, 27:157, <i>B*07:197</i> , <i>B*15:342</i> , <i>B*37:02</i> , <i>B*47:04-47:05</i> , <i>B*48:26</i> , <i>A*32:89</i>		2
*27:05:16, 27:06-27:07:06, 27:09, 27:11, 27:14, 27:19-27:21, 27:24, 27:27, 27:30, 27:32-27:35, 27:41, 27:43, 27:46, 27:70, 27:76, 27:78, 27:81, 27:91, 27:102, 27:106-27:107, 27:125, 27:127, 27:130, 27:136, 27:138, 27:150, 27:153-27:154, 27:164, 27:168, 27:180, <i>B*38:22</i> , <i>B*40:75</i> , <i>B*40:322</i> , <i>B*44:97</i> , <i>B*44:263</i> , <i>B*49:26</i>	1	
Well No.	1	2