

INTEPRETATION TABLE								
DPA1 SSP typing								
Amplification patterns of the DPA1 alleles								
	Well ⁵							
	1	2	3	4	5	6	7	8
Length of spec.	85	255	205	115	105	155	100	100
PCR product(s)						195		
Length of int.	515	515	430	430	430	515	430	430
pos. control ¹								
5'-primer(s) ²	15 (138) 5'-ACg ^{3'}	11 (125) 5'-CgC ^{3'}	28 (177) 5'-gAA ^{3'}	4 (103) 5'-Cgg ^{3'}	84 (345) 5'-AAT ^{3'}	31 (185) 5'-gCA ^{3'}	11 (125) 5'-CgC ^{3'}	11 (125) 5'-CAT ^{3'}
						43 (222) 5'-TgT ^{3'}		
3'-primer(s) ³	31 (184) 5'-CAT ^{3'}	83 (340) 5'-ggT ^{3'}	83 (340) 5'-ggT ^{3'}	28 (177) 5'-TCg ^{3'}	2 nd I 5'-ggC ^{3'}	69 (298) 5'-gTC ^{3'}	31 (184) 5'-CTg ^{3'}	31 (184) 5'-CTg ^{3'}
						83 (340) 5'-ggT ^{3'}		
Well No.	1	2	3	4	5	6	7	8
DPA1 allele ⁴								
*01:03:01-01:03:02, 01:03:04	1	2	3					
*01:03:03		2	3					
*01:04	1	2		4				
*01:05	1				5			
*01:06:01-01:06:02		2	3			6	7	
*01:07	1	2	3					
*01:08	1	2		4				
*01:09	1	2	3					
*01:10	1	2	3			6		
*02:01:01-02:01:07					5		7	
*02:02:01-02:02:04					5			8
*02:03					5			
*02:04					5	6		8
*03:01			3					
*03:02			3					
*03:03				4				
*04:01	1				5			
DPA1 allele ⁴								
Well No.	1	2	3	4	5	6	7	8

INTERPRETATION TABLE								
DPA1 SSP typing								
Amplification patterns of the DPA1 alleles								
Well ⁵								
9	10	11	12	13	14	15	16	
205	85	90	205	135	140	245	220	Length of spec. PCR product(s)
430	515	515	430	430	430	430	430	Length of int. pos. control ¹
11 (125)	15 (138)	66 (290)	18 (145)	51 (244)	50 (242)	15 (138)	23 (161)	5'-primer(s) ²
5'-CAT 3'	5'-ACC 3'	5'-ATC 3'	5'-gAA 3'	5'-AAA 3'	5'-CCg 3'	5'-ACC 3'	5'-TAC 3'	
66 (290)	31 (184)	83 (340)	73 (310)	83 (340)	83 (340)	83 (340)	83 (340)	3'-primer(s) ³
5'-TCA 3'	5'-CAT 3'	5'-ggT 3'	5'-AgC 3'	5'-ggT 3'	5'-ggT 3'	5'-ggT 3'	5'-ggT 3'	
9	10	11	12	13	14	15	16	Well No. DPA1 allele ⁴
	10					15		*01:03:01-01:03:02, 01:03:04 *01:03:03 *01:04 *01:05 *01:06:01-01:06:02 *01:07 *01:08 *01:09 *01:10
				13	14		16	
9								*02:01:01-02:01:07 *02:02:01-02:02:04 *02:03 *02:04
	10	11				15		*03:01 *03:02 *03:03
			12					*04:01
9	10	11	12	13	14	15	16	DPA1 allele ⁴ Well No.