

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
DPA1 allele <sup>1,2,3</sup>																						
*01:03:01:01-01:03:02, 01:03:04-01:03:09, 01:03:11-01:03:13, 01:03:15-01:03:20, 01:18-01:26, 01:28-01:32N	1	2	3																			
*01:03:03		2	3						10					15								
*01:03:10	1	2	3														17					
*01:03:14	1		3																			
*01:04	1	2		4																		
*01:05	1				5																	
*01:06:01-01:06:02		2	3			6	7															
*01:07	1	2	3							11				14								
*01:08	1	2		4										14								
*01:09	1	2	3										13									
*01:10	1	2	3			6																
*01:11	1	2	3											14								
*01:12	1	2	3							11												
*01:13	1		3			6																
*01:14	1	2	3																		19	
*01:15	1	2	3															18				
*01:16	1	2	3				7															
*01:17	1	2	3																		20	
*01:27		2	3																			
*02:01:01:01-02:01:02:02, 02:01:04, 02:01:08-01-02:01:13, 02:13N, 02:18-02:19, 02:24, 02:26, 02:29					5		7											17				
*02:01:03, 02:01:05-02:01:07					5		7											?				
*02:02:02:01-02:02:02:05, 02:02:02:07, 02:02:07-02:02:09, 02:14, 02:17					5			8	9									17				
*02:02:03-02:02:06, 02:25					5			8	9									?				
*02:03					5				10									?				
*02:04					5	6		8	9									?				
*02:05					5			8	9					14				?				
*02:06, 02:15					5			8	9				13					17				
*02:07:01:01-02:07:02, 02:12					5			8	9							16		17				
*02:08					5		7											17		19		
*02:09					5		7											17	18			
*02:10					5			8	9									17				21
*02:11			3		5		7											?				
*02:16					5		7						?					?				
*02:20								8	9									17				
*02:21		2	3		5	6	7							14	15			17				
*02:22					5				9	10								17				
*02:23, 02:28					5													17				
*02:27			3		5	6		8	9					14	15			17				
*03:01:01:01-03:01:01:06, 03:05Q			3							10	11				15			17				
*03:01:02			3							10	11				15							
*03:02			3						9	10					15			?				
*03:03				4						10	11				15			?				
*03:04			3							10	11			14	15			17				
*03:06			3			6		8			11				15							
*04:01:01:01-04:01:01:03	1				5							12										
*04:02:01:01-04:02:01:02	1				5							12						17				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22

Negative Control