

DPA1 2022 7N5

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DPA1 allele <sup>1,2,3</sup>																							
*01:03:01:01-01:03:02, 01:03:04-01:03:09, 01:03:11-01:03:13, 01:03:15-01:03:26, 01:03:28-01:03:31, 01:03:33-01:03:38, 01:18-01:26, 01:28-01:32N, 01:34-01:38, 01:41:01-01:45, 01:48-01:50, 01:54-01:57, 01:59-01:61, 01:63, 01:65, 01:67, 01:69-01:94	1	2	3																				
*01:03:03		2	3						10					15									
*01:03:10, 01:33:01-01:33:02	1	2	3														17						
*01:03:14, 01:39	1	2	3																				
*01:03:27	1	2	3		5																		
*01:03:32, 01:27		2	3																				
*01:04:01, 01:04:03	1	2		4																			
*01:04:02	1	2		4														?					
*01:05	1				5																		
*01:06:01-01:06:02		2	3			6	7																
*01:07	1	2	3							11				14									
*01:08	1	2		4										14									
*01:09	1	2	3										13										
*01:10	1	2	3			6																	
*01:11:01-01:11:02	1	2	3											14									
*01:12	1	2	3							11													
*01:13	1		3			6																	
*01:14, 01:68	1	2	3																			19	
*01:15	1	2	3															18					
*01:16	1	2	3				7																
*01:17	1	2	3																			20	
*01:40, 01:46-01:47, 01:51-01:53, 01:64, 01:66N	1	2	3																				
*01:58:01:01-01:58:01:02	1				5													17					
*01:62					5													17	18				
*02:01:01:01-02:01:01:04, 02:01:01:06-02:01:02:03, 02:01:04, 02:01:07-02:01:19, 02:13N, 02:18-02:19, 02:24, 02:26:01:01-02:26:01:02, 02:29, 02:31-02:32N, 02:34:01:01-02:34:01:02, 02:36-02:37, 02:39-02:40, 02:45:01-02:46, 02:49, 02:55, 02:57, 02:61					5		7											17					
*02:01:03, 02:01:05-02:01:06, 02:53					5		7																
*02:02:02:01-02:02:02:05, 02:02:02:07-02:02:02:13, 02:02:07-02:02:11, 02:14, 02:17, 02:35, 02:38Q, 02:42, 02:44, 02:48, 02:51, 02:54, 02:56Q, 02:58, 02:62-02:64					5			8	9									17					
*02:02:03-02:02:06, 02:25					5			8	9														
*02:03:01					5					10													
*02:03:02, 03:07:02					5					10													
*02:04					5	6		8	9														
*02:05					5			8	9					14									
*02:06, 02:15, 02:30					5			8	9				13										
*02:07:01:01-02:07:02, 02:52N					5			8	9								16						
*02:08					5		7											17		19			
*02:09:01:01-02:09:01:04					5		7											17	18				
*02:10					5			8	9									17				21	
*02:11																							
*02:12:01			3		5		7																
*02:12:02					5			8	9								16						22
					5			8	9									17					22
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	22

Negative Control

DPA1 2022 7N5

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	22
*02:16					5		7						?				?						
*02:20								8	9								17						
*02:21:01		2	3		5	6	7							14	15		17						
*02:21:02		2	3			6	7							14			?						
*02:22					5				9	10							17						
*02:23, 02:28:01-02:28:02					5												17						
*02:27:01			3		5	6		8	9					14	15		17						
*02:27:02			3			6		8	9					14	15	16	?						
*02:33					5			8	9														
*02:41N					5			8	9				?				?						
*02:43							7										17						
*02:47					5			8	9								?	18					
*02:50		2	3			6	7							14	15								
*02:59					5		7																
*02:60					5		7									16	17		19				22
*03:01:01:01-03:01:01:15, 03:01:03, 03:05:01:01Q-03:05:01:02Q, 03:09			3							10	11				15		17						
*03:01:02:01-03:01:02:04			3							10	11				15								
*03:02			3						9	10					15		?						
*03:03				4						10	11				15		?						
*03:04			3							10	11				14	15	17						
*03:06:01			3			6		8			11				15								
*03:06:02			3			6		8			11				15	16	?						
*03:07:01										10							?						
*03:08					5					10							17	18					
*04:01:01:01-04:01:01:03	1				5							12											
*04:02:01:01-04:03	1				5							12					17						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	22

Negative Control