

*04:01:01:01-04:01:02, 04:01:04-04:02:14, 04:03:01-04:03:02, 04:06-04:09, 04:11-04:21, 04:23, 04:25N-04:32, 04:34-04:37, 04:39-04:42	DQ4, Null, -								8	9									
*04:01:03	-								8	9			12						
*04:03:03	-									9			12						
*04:04-04:05, 04:33	-								8										
*04:10	-	2							8	9								13	
*04:22, 04:24, 04:38	-								8	w									
*05:01:01:01-05:01:15, 05:01:17-05:03:09, 05:03:11-05:03:17, 05:05:01-05:33, 05:35-05:43:02, 05:45-05:51, 05:53, 05:55-05:59, 05:61-05:71, 05:73-05:76, 05:78-05:81, 05:83-05:97, 05:99-05:104, 05:106-05:113, 05:115, 05:117-05:127, 05:129, 05:131, 05:133-05:145, 05:147-05:154	DQ5, Null, -	1											10						
*05:01:16, 05:60, 05:82, 05:98, 05:116	-												10						
*05:03:10, 05:04, 05:34, 05:44, 05:52, 05:54, 05:72, 05:77, 05:105, 05:114, 05:128N, 05:130, 05:132Q, 05:146	DQ5, Null, -	1																	
*06:01:01-06:01:06, 06:01:08-06:01:15, 06:43, 06:54N-06:58, 06:98, 06:100-06:105, 06:108, 06:120, 06:132, 06:140, 06:153, 06:157, 06:167, 06:177, 06:181, 06:194, 06:205, 06:214, 06:229, 06:239	DQ6, Null, -	2																	15
*06:01:07, 06:03:11, 06:04:01-06:05:01, 06:07:01-06:07:02, 06:09:01:01-06:09:07, 06:12, 06:17, 06:21-06:22:03, 06:25, 06:28, 06:34-06:36, 06:38-06:39, 06:42, 06:45, 06:52-06:53:02, 06:69:01-06:69:02, 06:85-06:86, 06:88-06:89, 06:92-06:94, 06:118:01-06:119, 06:121, 06:129, 06:135, 06:142, 06:144N, 06:149, 06:155, 06:158N, 06:160, 06:164, 06:171, 06:180, 06:186, 06:189, 06:193N, 06:199, 06:202-06:204, 06:207, 06:210, 06:212, 06:217, 06:231, 06:241	DQ6, DQ1, Null, -	2																	
*06:02:01:01-06:02:12, 06:02:14-06:02:28, 06:14:01-06:14:03, 06:16, 06:19:01-06:20, 06:23-06:24, 06:33, 06:37, 06:47-06:50, 06:51:02, 06:68, 06:70-06:78, 06:80-06:84, 06:95, 06:97, 06:107, 06:109, 06:111-06:117, 06:122, 06:124-06:127, 06:136-06:138, 06:147, 06:150-06:152, 06:159, 06:161, 06:163, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:188, 06:192, 06:197-06:198, 06:200-06:201, 06:208, 06:211, 06:213, 06:215-06:216N, 06:219, 06:224-06:228, 06:232, 06:236-06:237, 06:240, 06:242	DQ6, Null, -	2																13	14
*06:02:13, 06:15:01-06:15:02, 06:46, 06:79:01-06:79:02, 06:146:01-06:146:02, 06:235	DQ6, -	2																13	
*06:03:01:01-06:03:09, 06:03:12-06:03:27, 06:08:01-06:08:03, 06:10-06:11:03, 06:13:01-06:13:02, 06:18:01-06:18:02, 06:26N-06:27:02, 06:30-06:32:02, 06:40-06:41, 06:44, 06:59-06:65, 06:67, 06:87, 06:90-06:91, 06:99:01, 06:106, 06:110, 06:128, 06:130-06:131, 06:133-06:134, 06:141, 06:143, 06:145, 06:148, 06:154, 06:165, 06:170, 06:176, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196, 06:206:01-06:206:02, 06:218, 06:221-06:223, 06:230, 06:233-06:234, 06:238	DQ6, DQ1, Null, -	2																	14
*06:03:10, 06:51:01, 06:96	-	2						6											14
*06:05:02	DQ6	2																	?
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15			
*06:06	DQ5	2																?	?
*06:29, 06:123, 06:139	DQ6, -	2				5													14
*06:66, 06:172	-	2					6												
*06:99:02	-	2																	14 15
*06:156, 06:162	-										10							13	14
*06:168	-	2					6												15
*06:169	-										10								14
*06:209	-	2			4														15
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15			