

DQ-DR 2021 8L5 DQ

*03:03:11	-					6	7			11	12						
*03:03:26, 03:62, 03:78, 03:81, 03:110, 03:151, 03:225, 03:244, 03:272, 03:308, 03:314, 03:413	-						7				12						
*03:05:01, 03:05:05, 03:61, 03:181, 03:250, 03:346	DQ8, -					5	7			11							
*03:05:03	DQ8					5	7			11	12						
*03:06, 03:25:01-03:25:02, 03:411N	DQ3, Null, -						6				12						
*03:08, 03:107	DQ8, -					5	w				12						
*03:10:01-03:10:02:02, 03:14:01-03:14:02, 03:195	DQ3, -			4			7				12			15			
*03:12	DQ9						6	7			12			15			
*03:17:01	DQ9						7			11	12						
*03:17:02	-						7										
*03:23:01	-		2				7										
*03:23:02, 03:217, 03:259	-		2				7				12						
*03:65, 03:137, 03:283	-						6	w			12						
*03:70, 03:179	-					5	7				12			15			
*03:72, 03:100	-						7			11							
*03:104	-					5	6	7									
*03:132	-					5	7	8									
*03:147, 03:418	-			4							12						
*03:183	-			4			7							15			
*03:262	-					5	w										
*03:355	-		2	4			7				12						
*04:01:01:01-04:01:02, 04:01:04-04:02:01:01, 04:02:01:04-04:02:15, 04:02:17-04:02:18, 04:02:20, 04:06-04:09, 04:11-04:21, 04:23, 04:25N-04:32, 04:34-04:37, 04:39-04:48, 04:50-04:74, 04:76-04:84	DQ4, Null, -								8	9							
*04:01:03, 04:02:16	-								8	9	12						
*04:02:19, 04:04-04:05, 04:33, 04:49	-								8								
*04:03:01-04:03:02	-						6		8	9							
*04:03:03	-						6			9	12						
*04:10	-		2						8	9				13			
*04:22, 04:24, 04:38	-								8	w							
*04:75	-					5			8	9							
*05:01:01:01-05:01:15, 05:01:17-05:03:23, 05:05:01-05:43:02, 05:45-05:51, 05:53, 05:55-05:59, 05:61-05:71, 05:73-05:76, 05:78-05:81, 05:84-05:97, 05:99-05:104, 05:106-05:113, 05:115, 05:117-05:127, 05:129-05:131, 05:133-05:145, 05:147-05:174, 05:177-05:206N, 05:208N-05:217, 05:219-05:243, 05:245-05:257, 05:259-05:261, 05:263-05:271, 06:325 ⁴	DQ5, Null, -				1						10						
*05:01:16, 05:03:25, 05:60, 05:82-05:83, 05:98, 05:116, 05:175, 05:218, 05:258	-										10						
*05:03:24, 05:04, 05:44, 05:52, 05:54, 05:72, 05:77, 05:105, 05:114, 05:128N, 05:132Q, 05:146, 05:207, 05:244, 05:262	DQ5, Null, -				1												
*05:176	-		1	2													
*06:01:01:01-06:01:06, 06:01:08-06:01:27, 06:43, 06:54N-06:58, 06:98, 06:100-06:105, 06:108, 06:120, 06:132, 06:140, 06:153:01-06:153:02, 06:157, 06:167, 06:177, 06:181, 06:194, 06:205, 06:214, 06:229, 06:239, 06:243, 06:245-06:247, 06:251, 06:257-06:258, 06:260, 06:263, 06:268, 06:274, 06:277, 06:285, 06:302, 06:305, 06:307, 06:309-06:310, 06:312, 06:321, 06:323, 06:330N, 06:342, 06:359, 06:382	DQ6, Null, -				2									15			
Well No.			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

DQ-DR 2021 8L5 DQ

*06:01:07, 06:03:11, 06:04:01:01-06:05:01, 06:07:01-06:07:02, 06:09:01:01-06:09:12, 06:12, 06:17, 06:21-06:22:03, 06:25, 06:28, 06:34-06:36, 06:38-06:39, 06:42, 06:45, 06:52-06:53:02, 06:69:01-06:69:02, 06:85-06:86, 06:88:01:01-06:89, 06:92:01-06:94, 06:118:01-06:119, 06:121, 06:129, 06:135, 06:142, 06:144N, 06:149, 06:155, 06:158N, 06:160, 06:164, 06:171, 06:180, 06:186, 06:189, 06:193N, 06:199, 06:202-06:204, 06:207, 06:210, 06:212, 06:217, 06:231, 06:241, 06:252N, 06:254, 06:261, 06:265-06:267, 06:275, 06:280-06:283, 06:287-06:288, 06:291-06:292, 06:299, 06:301, 06:303N, 06:313, 06:318, 06:320, 06:332, 06:339-06:340, 06:343, 06:348-06:349, 06:351, 06:353, 06:358, 06:361, 06:369, 06:375, 06:381	DQ6, DQ1, Null	2															
*06:02:01:01-06:02:12, 06:02:14-06:02:42, 06:02:44-06:02:47, 06:14:01-06:14:03, 06:16, 06:19:01-06:20, 06:23-06:24, 06:33, 06:37, 06:47-06:50, 06:51:02, 06:68, 06:70-06:78, 06:80-06:84:01:02, 06:95, 06:96:02-06:97, 06:107, 06:109, 06:111-06:117, 06:122, 06:124-06:127, 06:136-06:138, 06:147, 06:150-06:152, 06:159, 06:161, 06:163, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:188, 06:192, 06:197-06:198, 06:200-06:201, 06:208, 06:211, 06:213, 06:215-06:216N, 06:219, 06:224-06:228, 06:232, 06:236-06:237, 06:240, 06:242, 06:249, 06:255-06:256, 06:262, 06:264, 06:270:02-06:271, 06:273, 06:284, 06:286, 06:289-06:290, 06:293-06:298, 06:300, 06:304N, 06:306N, 06:308N, 06:311, 06:314-06:315, 06:317N, 06:324, 06:326, 06:333, 06:335, 06:338, 06:341N, 06:344, 06:347, 06:354-06:357, 06:363-06:364, 06:366, 06:370, 06:372, 06:374, 06:376, 06:379N-06:380	DQ6, Null, -	2											13	14			
*06:02:13, 06:15:01-06:15:02, 06:46, 06:79:01-06:79:02, 06:146:01-06:146:02, 06:235	DQ6, -	2											13				
*06:02:43, 06:03:10, 06:03:33, 06:51:01, 06:96:01, 06:322	-	2			6									14			
*06:03:01:01-06:03:09, 06:03:12-06:03:32, 06:03:34-06:03:40, 06:08:01-06:08:03, 06:10-06:11:04, 06:13:01-06:13:03, 06:18:01-06:18:02, 06:26N-06:27:02, 06:30-06:32:02, 06:40-06:41:01:03, 06:44, 06:59-06:65, 06:67, 06:87, 06:90-06:91, 06:99:01, 06:106, 06:110, 06:128, 06:130-06:131, 06:133-06:134, 06:141, 06:143, 06:145, 06:148, 06:154, 06:165, 06:170, 06:176, 06:184-06:185, 06:187, 06:190:01-06:191, 06:195-06:196, 06:206:01-06:206:02, 06:218, 06:221-06:223, 06:230, 06:233-06:234, 06:238, 06:244, 06:248, 06:250, 06:253, 06:259, 06:269-06:270:01, 06:272, 06:276, 06:278-06:279, 06:316, 06:319, 06:327-06:329, 06:331, 06:334, 06:336, 06:345N-06:346, 06:350, 06:352, 06:360, 06:362, 06:365, 06:367-06:368, 06:371, 06:373, 06:378	DQ6, DQ1, Null, -	2												14			
*06:05:02	DQ6	2															?
*06:06	DQ5	2															? ?
*06:29, 06:123, 06:139, 06:337	DQ6, -	2			5									14			
*06:66, 06:172	-	2			6												
*06:99:02-06:99:03	-	2												14	15		
*06:156, 06:162	-								10				13	14			
*06:168, 06:377	-	2			6												15
*06:169	-								10					14			
*06:209	-	2		4													15
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	