

Well No.	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
*08:01:06, 08:03:02:01-08:03:10, 08:03:12-08:03:16, 08:13, 08:18, 08:23, 08:27, 08:29-08:30:03, 08:33, 08:36:01, 08:38, 08:46-08:49, 08:51, 08:56, 08:58, 08:60N-08:62, 08:65-08:66, 08:68:02-08:69, 08:71, 08:74, 08:76, 08:81, 08:91, 08:95, 08:97, 08:102, 08:104, 08:110, 08:115, 08:119, 08:121, 08:128, 08:131	DR8, Null, -									25											36												
*08:03:11, 08:68:01	-									25										w	36												
*08:04:01:01-08:04:05, 08:04:07-08:04:08, 08:06:01:01-08:06:01:02, 08:22, 08:54, 08:59, 08:83, 08:94, 08:99, 08:122, 08:124, 08:126	DR8, -									25			31									36	37										
*08:04:06	-									25			31									37											
*08:07, 08:11:01:01-08:11:01:02, 08:52	DR8, -									25			31																				
*08:08	DR8									25			31		33																		
*08:09:01-08:09:02	DR8, -									25			31								35	36					42						
*08:10, 08:12, 08:14	DR8									25												36	37										
*08:15, 08:19, 08:25, 08:34	DR8, -									25																							
*08:17, 08:103	DR8, -									25			31									36	38										
*08:20	-												30	31	32						35	37					42						
*08:21	DR8									25			30	31							35	36					42						
*08:28	DR8									25				31								36	37	38									
*08:31, 11:05, 11:309 ⁴	DR8, DR11, -								28				31																				
*08:32, 08:53	-								25	29											35	36											
*08:35	-								25													35	36					42					
*08:36:02	-								25													35	36										
*08:37, 08:45:01-08:45:02	-								25													36	38										
*08:40	-				21				25													36											
*08:41	-								25	28			31																				
*08:57	-												31									36											
*08:63	-								25	28														38									
*08:67, 08:80, 08:98	-												31									36	37										
*08:70	-								25				31							w		36	37										
*08:75	-								25	28			31										37										
*08:87	-								25		29		31						w			36											
*08:92	-				21				25				31									36											
*08:101	-							24	25				31									36											
*08:113	-								25				31											38									
*08:129, 13:55 ⁴	-												30	31	32						35			39		42							
*09:01:02:01-09:01:05, 09:01:07-09:02:02, 09:04-09:22, 09:24-09:28, 09:30-09:52N, 09:54-09:60	DR9, Null, -								26								34																
*09:01:06, 09:03, 09:23, 09:29	DR9, -								26																								
*09:53	-								26	27																							
*10:01:01:01, 10:01:01:03-10:18, 10:20-10:46	DR10, -											27																					
*10:19	-											27																		45			
Well No.	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
*11:01:01:01, 11:01:01:03-11:01:06, 11:01:08-11:01:17, 11:01:20-11:01:27, 11:01:29-11:01:47, 11:01:49-11:01:50, 11:04:01:01-11:04:13, 11:04:15-11:04:16, 11:04:18-11:04:20, 11:04:22-11:04:24, 11:06:01-11:06:03, 11:09:01:01-11:10:02, 11:12:01-11:12:03, 11:15:01:01-11:15:02, 11:24:01-11:24:02, 11:27:01, 11:27:04-11:29:02, 11:32-11:33, 11:35, 11:38-11:39, 11:43:01:01-11:44, 11:46:01-11:47, 11:49:01-11:49:02, 11:51, 11:54:01-11:54:02, 11:58:01-11:58:02, 11:60-11:62:02, 11:66:01-11:66:02, 11:72, 11:74:01-11:75, 11:77-11:78, 11:81, 11:84:01, 11:90-11:92, 11:94-11:95, 11:97, 11:99-11:102:02, 11:106, 11:108-11:117:02, 11:120-11:121, 11:123, 11:126, 11:128-11:129, 11:133-11:134, 11:137, 11:140-11:141, 11:144, 11:146-11:147:02, 11:150, 11:152, 11:154-11:158, 11:160, 11:162-11:167, 11:169N, 11:172, 11:174-11:175, 11:178-11:181, 11:185-11:188, 11:190, 11:192, 11:194:01-11:199, 11:201, 11:203-11:208, 11:210-11:212, 11:214-11:228, 11:230, 11:232, 11:234-11:236, 11:238-11:240, 11:243-11:247, 11:250N, 11:252-11:255, 11:257-11:259, 11:262-11:263, 11:265, 11:269, 11:271-11:272Q, 11:274-11:275, 11:277-11:279, 11:281-11:282, 11:284-11:285:02, 11:287N, 11:289, 11:291-11:292, 11:294N-11:300, 11:302-11:308, 11:310-11:311, 11:313-11:318, 11:320-11:331, 11:333, 11:336-11:337	DR11, Null, -												28		30	31											41						
*11:01:07, 11:04:14, 11:37:01-11:37:02, 11:50, 11:56, 11:145, 11:170, 11:183, 11:189, 11:200, 11:249, 11:273, 11:332	DR11, -												28		30	31																	
*11:01:18-11:01:19, 11:01:48, 11:04:21, 11:08:01-11:08:03, 11:18-11:19:03, 11:42:01-11:42:02, 11:57, 11:130, 11:191, 11:248Q, 11:268, 11:280, 11:286	DR11, -												28		30													41					
*11:01:28	-												28		?	31											?						

Negative Control

Control

