

INTERPRETATION TABLE						
DQB1*02 PCR-SSP subtyping						
Allele-specific amplification patterns of the 0201 to 0205 alleles						
	Well					
	1	2	3	4	5	6
Length of spec.	120	85	145	140	145	180
PCR product						
Length of int.	515	430	515	430	430	430
pos. control¹						
5'-primer²	30	57	135	102	136	30
	5'-A Ag ^{3'}	5'-T gA ^{3'}	5'-T gA ^{3'}	5'-TC T ^{3'}	5'-AC g ^{3'}	5'-A Ag ^{3'}
3'-primer³	57	71	170	135	170	77
	5'-C gg ^{3'}	5'-Cg T ^{3'}	5'-gAC ^{3'}	5'-g gC ^{3'}	5'-gAC ^{3'}	5'-CCg ^{3'}
Well No.	1	2	3	4	5	6
DQB1 allele⁴						
*020101-020102	1		3			
*0202	1			4		
*0203		2		4		
*0204	1		3		5	
*0205	1		3			6
*030101-030204, 030302-0321, 0401-0403, 050101-0505, 060101-0634			3			
DQB1 allele⁴						
Well No.	1	2	3	4	5	6