

INTERPRETATION TABLE						
DQB1*02 PCR-SSP subtyping						
Allele-specific amplification patterns of the 0201 to 0205 alleles						
	Well					
	1	2	3	4	5	6
Length of spec.	120	85	145	140	145	180
PCR product						
Length of int.	515	430	515	430	430	430
pos. control¹						
5'-primer²	30(185)	57(266)	135(500)	102(400)	136(502)	30(185)
	5'-AAG 3'	5'-TgA 3'	5'-TgA 3'	5'-TCT 3'	5'-ACg 3'	5'-AAG 3'
3'-primer³	57(266)	71(309)	170(604)	135(500)	170(604)	77(326)
	5'-Cgg 3'	5'-CgT 3'	5'-gAC 3'	5'-ggC 3'	5'-gAC 3'	5'-CCg 3'
Well No.	1	2	3	4	5	6
DQB1 allele⁴						
*020101-020102	1		3			
*0202	1			4		
*0203		2		4		
*0204	1		3		5	
*0205	1		3			6
*030101-0323, 0325, 040101-0403, 050101- 0505, 060101-0635			3			
DQB1 allele⁴						
Well No.	1	2	3	4	5	6