

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DQB1 allele <sup>1,2,3</sup>																								
*02:01:01-02:01:17	1		3																					
*02:02:01-02:02:02	1			4																				
*02:03		2		4																				
*02:04	1		3		5																			
*02:05	1		3			6																		
*02:06	1			4			7																	
*02:07:01-02:07:02, 02:16 <sup>4</sup>	1		3						9															
*02:08	1		3							10														
*02:09, 02:24 <sup>5</sup>	1		3								11													
*02:10	1			4				8																
*02:11	1			4								12												
*02:12	1			4						10														
*02:13	1		3									12												
*02:14	1		3																	20				
*02:15, 02:29 <sup>5</sup>	1		3														18							
*02:17			3										13											
*02:18N, 02:34 <sup>5</sup>	1		3											14										
*02:19	1		3				7																	
*02:20N, 02:22 <sup>5</sup>	1		3													16								
*02:21, 02:35 <sup>5</sup>	1		3												15									
*02:23	1		3		5																21			
*02:25			3									12												
*02:26	1		3														17							
*02:27-02:28 <sup>5</sup>	1		3																			22		
*02:30	1		3					8																
*02:31	1		3																			21		
*02:32	1		3										13											
*02:33	1		3																	19				
*02:36	1		3																	19	20			
*02:37	1		3		5												17							
*02:38	1		3										13	14										
*02:39	1		3												15	16								
*02:40			3																			21		
*02:41	1		3																				23	
*03:01:01:01-03:23, 03:25-03:48, 03:50-03:78, 03:80-03:96, 03:98-03:150, 04:01:01-04:23, 05:01:01:01-05:13, 05:15-05:78, 06:01:01-06:37, 06:39-06:85, 06:87-06:101, 06:105-06:159			3																					
*03:24, 03:79							7																	
*03:49			3					8																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control