

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
DQB1 allele ^{1,2,3}																					
*02:01:01-02:01:14	1		3																		
*02:02:01-02:02:02	1			4																	
*02:03		2		4																	
*02:04	1		3		5																
*02:05	1		3			6															
*02:06	1			4			7														
*02:07, 02:16 ⁴	1		3						9												
*02:08	1		3							10											
*02:09, 02:24 ⁴	1		3								11										
*02:10	1			4				8													
*02:11	1			4								12									
*02:12	1			4						10											
*02:13	1		3									12									
*02:14	1		3																	20	
*02:15, 02:29 ⁴	1		3															18			
*02:17			3										13								
*02:18N, 02:34 ⁴	1		3											14							
*02:19	1		3				7														
*02:20N, 02:22 ⁴	1		3													16					
*02:21, 02:35 ⁴	1		3												15						
*02:23	1		3		5																21
*02:25			3									12									
*02:26	1		3														17				
*02:27-02:28 ⁴	1		3																		
*02:30	1		3					8													
*02:31	1		3																		21
*02:32	1		3										13								
*02:33	1		3																	19	
*03:01:01:01-03:23, 03:25-03:48, 03:50-03:78, 03:80-03:96, 03:98-03:113, 04:01:01-04:18, 05:01:01:01-05:13, 05:15-05:55, 06:01:01-06:37, 06:39-06:85, 06:87-06:101, 06:105-06:129			3																		
*03:24, 03:79							7														
*03:49			3					8													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

