

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DQB1 allele ^{1,2,3}																								
*02:01:01-02:01:23, 02:42-02:47, 02:49-02:52, 02:55-02:57, 02:60-02:61	1		3																					
*02:02:01-02:02:02	1			4																				
*02:03		2		4																				
*02:04	1		3		5																			
*02:05	1		3			6																		
*02:06	1			4			7																	
*02:07:01-02:07:02, 02:16 ⁴	1		3						9															
*02:08	1		3							10														
*02:09, 02:24 ⁴	1		3								11													
*02:10	1			4				8																
*02:11	1			4								12												
*02:12	1			4						10														
*02:13	1		3									12												
*02:14:01-02:14:02	1		3																		20			
*02:15, 02:29 ⁴	1		3															18						
*02:17			3										13											
*02:18N, 02:34 ⁴	1		3											14										
*02:19	1		3				7																	
*02:20N, 02:22 ⁴	1		3													16								
*02:21, 02:35 ⁴	1		3												15									
*02:23	1		3		5																	21		
*02:25			3									12												
*02:26	1		3															17						
*02:27-02:28 ⁴	1		3																			22		
*02:30	1		3					8																
*02:31	1		3																			21		
*02:32	1		3										13											
*02:33	1		3																19					
*02:36	1		3																19	20				
*02:37	1		3		5													17						
*02:38	1		3										13	14										
*02:39	1		3												15	16								
*02:40			3																			21		
*02:41, 02:53Q ⁴	1		3																				23	
*02:48	1						7											17						
*02:54	1		3											14		16								
*02:58N	1		3																18			21		
*02:59	1		3																	19			22	
*03:01:01:01-03:23:02, 03:25-03:48, 03:50-03:78, 03:80-03:96, 03:98-03:163, 03:166-03:191, 04:01:01-04:29, 05:01:01:01-05:02:10, 05:02:12-05:13, 05:15-05:83, 05:85-05:96, 06:01:01-06:37, 06:39-06:43, 06:45-06:46, 06:48-06:85, 06:87-06:101, 06:105-06:182			3																					
*03:24, 03:79							7											17						
*03:49			3					8																
*06:44, 06:47			3																	19			22	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control