

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
DQB1 allele ^{1,2,3}																																	
*02:01:01-02:01:24, 02:43-02:45, 02:47, 02:49, 02:52, 02:55-02:57, 02:60-02:61, 02:63, 02:66, 02:68-02:71, 02:73-02:76, 02:78	1		3																														
*02:02:01-01-02:02:03, 02:10, 02:65 ⁵	1			4																													
*02:03		2		4																													
*02:04	1		3		5																												
*02:05	1		3			6																											
*02:06	1			4			7																										
*02:07:01-02:07:02, 02:16 ⁴	1		3						9																								
*02:08	1		3							10																							
*02:09, 02:24 ⁴	1		3								11																						
*02:11	1			4								12																					
*02:12	1			4						10																							
*02:13	1		3									12																					
*02:14:01-02:14:02	1		3																	20													
*02:15, 02:29 ⁴	1		3															18															
*02:17			3										13																				
*02:18N, 02:34 ⁴	1		3											14																			
*02:19	1		3				7																										
*02:20N, 02:22 ⁴	1		3												16																		
*02:21, 02:35 ⁴	1		3												15																		
*02:23	1		3		5																21												
*02:25			3									12																					
*02:26	1			4													17																
*02:27-02:28 ⁴	1		3																		22												
*02:30	1		3									12								20													
*02:31	1		3																		21												
*02:32	1		3										13																				
*02:33	1		3																19														
*02:36	1		3																19	20													
*02:37	1		3		5												17																
*02:38	1		3										13	14																			
*02:39	1		3												15	16																	
*02:40			3																		21												
*02:41, 02:53Q ⁴	1		3																			23											
*02:42	1		3																				24										
*02:46	1		3																					25									
*02:48	1						7										17																
*02:50	1			4																						26							
*02:51	1		3																					24	25								
*02:54	1		3											14	16																		
*02:58N	1		3																											30			
*02:59	1		3																		22										31		
*02:62	1			4																								27					
*02:64	1			4																									28				
*02:67N	1		3																											29			
*02:72			3																									27					
*02:77		2	3																														
*03:01:01-01-03:20, 03:22-03:23:02, 03:25-01-03:78, 03:80-03:96, 03:98-03:163, 03:166-03:240, 03:243, 04:01:01-04:38, 05:01:01-01-05:02:10, 05:02:12-05:13, 05:15-05:34, 05:36-05:83, 05:85-05:131, 06:01:01-06:37, 06:39-06:43, 06:45-06:46, 06:48-06:85, 06:87-06:101, 06:105-06:215			3																														
*03:21, 05:35			3																										28				
*03:24, 03:79							7										17																
*06:44, 06:47			3																		22										31		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control