

INTERPRETATION TABLE				
DQB1*04 SSP subtyping				
Amplification patterns of the DQB1*04 alleles				
	Tube			
	1	2	3	
Length of spec.	195	205	205	Length of spec.
PCR product				PCR product
Length of int.	515	515	430	Length of int.
pos. control¹				pos. control¹
5'-primer(s)²	26	23	23	5'-primer(s)²
	5'-g gg ^{3'}	5'-g CT ^{3'}	5'-g Cg ^{3'}	
	26			
	5'-g gg ^{3'}			
3'-primer³	77	77	77	3'-primer³
	5'-AC g ^{3'}	5'-AC g ^{3'}	5'-AC g ^{3'}	
Tube No.	1	2	3	Tube No.
DQB1 allele				DQB1 allele
*0401	+	+		*0401
*0402	+		+	*0402
DQB1 allele				DQB1 allele
Tube No.	1	2	3	Tube No.