

INTERPRETATION TABLE				
DQB1*04 SSP subtyping				
Amplification patterns of the DQB1*04 alleles				
	Well			
	1	2	3	4
Length of spec.	210	205	205	195
PCR product				
Length of int.	515	515	430	430
pos. control¹				
5'-primer²	21	23	23	26
	5'-ACC3'	5'-g CT3'	5'-g Cg3'	5'-T CT3'
3'-primer³	77	77	77	77
	5'-AC g3'	5'-AC g3'	5'-AC g3'	5'-AC g3'
Well No.	1	2	3	4
DQB1 allele				
*0401	1	2		
*0402	1		3	
*0403	1		3	4
*0306				4
DQB1 allele				
Well No.	1	2	3	4