

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
DQB1 allele ^{1,2,3}															
*04:01:01-04:01:04	1	2													
*04:02:01-04:02:05	1		3												
*04:03:01	1		3	4											
*04:03:02	1			4											
*04:04, 04:12 ⁴	1		3		5										
*04:05-04:06 ⁴	1	2			5										
*04:07	1	2					7								
*04:08	1	2						8							
*04:09	1		3						9						
*04:10	1		3							10					
*04:11	1		3								11				
*04:13	1		3										13		
*04:14	1	2							9						
*04:15	1	2									11				
*04:16	1	2				6									
*04:17	1	2												14	
*04:18	1		3									12			
*03:06				4	5			?							
*03:07-03:08, 03:10:02-03:15, 03:17:01-03:18, 03:20, 03:23, 03:26, 03:37, 03:40, 03:48, 03:52-03:71, 03:74-03:78, 03:81-03:82, 03:101- 03:112								?							
*03:19								8							
*03:25				4	5										
*06:03:03										w					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

Negative Control