

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DQB1 allele^{1,2,3}																
*04:01:01-04:01:04	1	2														
*04:02:01-04:02:07, 04:19, 04:21, 04:26-04:27	1		3													
*04:02:08			3				7									
*04:03:01	1		3	4												
*04:03:02	1			4												
*04:04, 04:12 ⁴	1		3		5											
*04:05-04:06 ⁴	1	2			5											
*04:07	1	2										12				
*04:08	1	2						8								
*04:09	1		3						9							
*04:10	1		3							10						
*04:11	1		3								11					
*04:13	1		3										13			
*04:14	1	2							9							
*04:15	1	2									11					
*04:16	1	2				6										
*04:17	1	2												14		
*04:18	1		3									12				
*04:20	1				6											
*04:22, 04:24	1		w													
*04:23	1		3								11					15
*04:25N	1		3													15
*03:06				4	5			?								
*03:07-03:08, 03:10:02-03:15, 03:17:01-03:18, 03:20, 03:23, 03:26, 03:37, 03:40, 03:48, 03:52-03:71, 03:74-03:78, 03:81-03:82, 03:101-03:112, 03:118N-03:131, 03:133-03:163, 03:165-03:167								?								
*03:19								8								
*03:22, 03:96											11					15
*03:25				4	5											
*03:30, 03:72, 03:100, 06:02:01-06:02:23, 06:05:02, 06:10, 06:13:01-06:16, 06:19:01-06:20, 06:22:01-06:24, 06:29, 06:33, 06:37, 06:46-06:51:02, 06:68-06:84, 06:95-06:97, 06:106-06:107, 06:109, 06:111-06:117, 06:119, 06:122-06:123, 06:125-06:127, 06:130, 06:136, 06:138-06:139, 06:146-06:147, 06:150-06:152, 06:156, 06:159, 06:161-06:163, 06:166								7								
*03:132	1							?								
*05:01:19, 06:03:03										w						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control